CellSave Preservative 20 Tubes



 Version
 Revision Date:
 SDS Number:
 Date of last issue: 27.10.2016

 1.11
 09.01.2018
 100000011150
 Date of first issue: 18.09.2015

SECTION 1. PRODUCT AND COMPANY IDENTIFICATION

Product name : CellSave Preservative 20 Tubes

Substance name : CellSave Preservative 20 Tubes

Manufacturer or supplier's details

Company : Menarini Silicon Biosystems, Inc.

Address : 3401 Masons Mill Rd #100

Huntingdon Valley, PA 19006,

USA

Telephone : 1 (800) 381-4929

Emergency telephone

number

US: (303)-389-1805

International: +1 (303)-389-1805

E-mail address

Responsible/issuing person

Us-info@siliconbiosystems.com

Recommended use of the chemical and restrictions on use

Recommended use : Assay reagent

SECTION 2. HAZARDS IDENTIFICATION

GHS Classification

Skin sensitisation : Category 1

GHS label elements

Hazard pictograms



Signal word : Warning

Hazard statements : H317 May cause an allergic skin reaction.

Precautionary statements : **Prevention:**

P261 Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray. P272 Contaminated work clothing should not be allowed out of

the workplace.

P280 Wear protective gloves.

Response:

P302 + P352 IF ON SKIN: Wash with plenty of soap and water. P333 + P313 If skin irritation or rash occurs: Get medical ad-

vice/ attention.

P363 Wash contaminated clothing before reuse.

Disposal:

P501 Dispose of contents/ container to an approved waste

CellSave Preservative 20 Tubes



Version 1.11

Revision Date: 09.01.2018

SDS Number: 100000011150

Date of last issue: 27.10.2016 Date of first issue: 18.09.2015

disposal plant.

Other hazards which do not result in classification

None known.

SECTION 3. COMPOSITION/INFORMATION ON INGREDIENTS

Substance / Mixture : Mixture

Chemical nature : Liquid

Hazardous components

Chemical name	CAS-No.	Concentration (% w/w)
IMIDUREA	39236-46-9	>= 30 - < 60

SECTION 4. FIRST AID MEASURES

If inhaled : If breathed in, move person into fresh air.

Consult a physician.

In case of skin contact : Take off contaminated clothing and shoes immediately.

Wash off immediately with plenty of water. If symptoms persist, call a physician. Wash contaminated clothing before re-use.

In case of eye contact : Rinse immediately with plenty of water, also under the eyelids,

for at least 5 minutes. Remove contact lenses.

If eye irritation persists, consult a specialist.

If swallowed, rinse mouth with water (only if the person is con-

scious).

Call a physician immediately.

Most important symptoms and effects, both acute and

delayed

: No information available.

Notes to physician : Treat symptomatically.

SECTION 5. FIREFIGHTING MEASURES

Suitable extinguishing media : Use extinguishing measures that are appropriate to local cir-

cumstances and the surrounding environment.

Specific hazards during fire-

fighting

: No information available.

Specific extinguishing meth-

ods

: No information available.

CellSave Preservative 20 Tubes



Version 1.11

Revision Date: 09.01.2018

SDS Number: 100000011150 Date of last issue: 27.10.2016 Date of first issue: 18.09.2015

Special protective equipment

for firefighters

: In the event of fire, wear self-contained breathing apparatus.

SECTION 6. ACCIDENTAL RELEASE MEASURES

tive equipment and emer-

gency procedures

Personal precautions, protec-: In the event of an accidental release the emergency response team must respond based on a risk assessment and use personal protective equipment as appropriate.

Environmental precautions

: Should not be released into the environment.

Methods and materials for containment and cleaning up : Large spills: Dam up. Soak up with inert absorbent material.

Keep in properly labelled containers.

Small spills: Gently cover the spill with an absorbent towel or

pad.

Large spills + Small spills: Keep in suitable, closed containers for disposal. Treat recovered material as described in the sec-

tion "Disposal considerations".

SECTION 7. HANDLING AND STORAGE

Advice on protection against

fire and explosion

: No data available

: To avoid thermal decomposition, do not overheat. Advice on safe handling

Avoid inhalation, ingestion and contact with skin and eyes.

Use personal protective equipment as required.

Hygiene measures : Handle in accordance with good industrial hygiene and safety

practice.

Conditions for safe storage

To maintain product quality, do not store in heat or direct sun-

light.

Store in original container.

Keep containers tightly closed in a dry, cool and well-

ventilated place.

Keep away from heat and sources of ignition.

Keep locked up. Keep refrigerated.

Recommended storage tem-

perature

: 4 - 30 °C

SECTION 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Components with workplace control parameters

Contains no substances with occupational exposure limit values.

Engineering measures : All personal protective equipment should be based on a risk

assessment. Consult a Environment Health Safety expert if

necessary.

CellSave Preservative 20 Tubes



Version 1.11

Revision Date: 09.01.2018

SDS Number: 100000011150

Date of last issue: 27.10.2016 Date of first issue: 18.09.2015

Personal protective equipment

Respiratory protection

: Engineering controls should always be the primary method of

controlling exposures.

If respiratory protective equipment is needed for certain activities, the type as well as the corresponding protection factor will depend upon the risk assessment and air concentrations, hazards, physical and warning properties of substances pre-

sent.

No personal respiratory protective equipment normally re-

quired.

Hand protection

Remarks : Disposable gloves

Eye protection : No special precautions required.

Skin and body protection : No special precautions required.

Protective measures : The type of protective equipment must be selected based on

the Environmental Health and Safety risk assessment. Consult a Environmental Health and Safety expert if necessary.

SECTION 9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance : liquid

Colour : clear

SECTION 10. STABILITY AND REACTIVITY

Reactivity : None reasonably foreseeable.

Chemical stability : Stable under recommended storage conditions.

Possibility of hazardous reac-

tions

: No dangerous reaction known under conditions of normal use.

Conditions to avoid : To avoid thermal decomposition, do not overheat.

Incompatible materials : None known.

Hazardous decomposition

products

: None known.

SECTION 11. TOXICOLOGICAL INFORMATION

Acute toxicity

Components:

CellSave Preservative 20 Tubes



 Version
 Revision Date:
 SDS Number:
 Date of last issue: 27.10.2016

 1.11
 09.01.2018
 100000011150
 Date of first issue: 18.09.2015

IMIDUREA:

Acute oral toxicity : LD50 (Rat): 11,300 mg/kg

Acute inhalation toxicity : LC50 (Rat): > 5.5 mg/l

Exposure time: 1 h

Acute dermal toxicity : LD50 (Rabbit): > 5,000 mg/kg

Skin corrosion/irritation

Components:

IMIDUREA:

Result: No skin irritation

Serious eye damage/eye irritation

Components:

IMIDUREA:

Result: No eye irritation

Respiratory or skin sensitisation

Components:

IMIDUREA:

Method: Maximisation Test

Result: May cause sensitisation by skin contact.

Method: Local Lymph Node Assay (LLNA) in mice Result: May cause sensitisation by skin contact.

Chronic toxicity

Germ cell mutagenicity

Components:

IMIDUREA:

Genotoxicity in vitro : Test Type: Ames test

Metabolic activation: with and without metabolic activation

Result: negative

GLP: yes

: Test Type: Chromosome aberration test in vitro

Metabolic activation: with and without metabolic activation

Result: negative

Genotoxicity in vivo : Test Type: Micronucleus test

Species: Mouse

Application Route: Oral

Result: negative

Germ cell mutagenicity -

Assessment

: No information available.

CellSave Preservative 20 Tubes



Version 1.11

Revision Date: 09.01.2018

SDS Number: 100000011150

Date of last issue: 27.10.2016 Date of first issue: 18.09.2015

Carcinogenicity

Components: IMIDUREA:

Carcinogenicity -Assessment : No information available.

Reproductive toxicity

Components:

IMIDUREA:

Teratogenicity - Assessment : No information available.

STOT - single exposure

No data available

STOT - repeated exposure

No data available

Repeated dose toxicity

Components:

IMIDUREA:

Species: Rat NOAEL: 200 mg/kg LOAEL: 500 mg/kg Application Route: Oral

Species: Rabbit NOAEL: 200 mg/kg Application Route: Dermal

Aspiration toxicity

No data available

SECTION 12. ECOLOGICAL INFORMATION

Ecotoxicity

Components:

IMIDUREA:

Toxicity to fish : Remarks: No data available

Persistence and degradability

Components:

IMIDUREA:

Biodegradability : Remarks: No data available

CellSave Preservative 20 Tubes



 Version
 Revision Date:
 SDS Number:
 Date of last issue: 27.10.2016

 1.11
 09.01.2018
 100000011150
 Date of first issue: 18.09.2015

Bioaccumulative potential

Components: IMIDUREA:

Bioaccumulation : Remarks: No data available

Mobility in soil

Components:

IMIDUREA:

Distribution among : Remarks: No data available

environmental compartments

Other adverse effects

Components: IMIDUREA:

Results of PBT and vPvB

assessment

Additional ecological

information

: No information available.

: No data available

SECTION 13. DISPOSAL CONSIDERATIONS

Disposal methods

Waste from residues : In accordance with National, Federal, State and Local regula-

tions.

Contaminated packaging : Empty containers should be taken to an approved waste han-

dling site for recycling or disposal.

SECTION 14. TRANSPORT INFORMATION

International Regulations

UNRTDG

Not regulated as a dangerous good

IATA-DGR

Not regulated as a dangerous good

IMDG-Code

Not regulated as a dangerous good

Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable for product as supplied.

National Regulations

ADG

Not regulated as a dangerous good

SAFETY DATA SHEET CellSave Preservative 20 Tubes



 Version
 Revision Date:
 SDS Number:
 Date of last issue: 27.10.2016

 1.11
 09.01.2018
 100000011150
 Date of first issue: 18.09.2015

SECTION 15. REGULATORY INFORMATION

Safety, health and environmental regulations/legislation specific for the substance or mixture

Restricted to professional users.

R-phrase(s) : R43 May cause sensitisation by skin contact.

S-phrase(s) : S24 Avoid contact with skin.

S37 Wear suitable gloves.

S60 This material and its container must be

disposed of as hazardous waste.

SECTION 16. OTHER INFORMATION

Full text of other abbreviations

AICS - Australian Inventory of Chemical Substances; ANTT - National Agency for Transport by Land of Brazil; ASTM - American Society for the Testing of Materials; bw - Body weight; CMR -Carcinogen, Mutagen or Reproductive Toxicant; CPR - Controlled Products Regulations; DIN -Standard of the German Institute for Standardisation; DSL - Domestic Substances List (Canada); ECx - Concentration associated with x% response; ELx - Loading rate associated with x% response; EmS - Emergency Schedule; ENCS - Existing and New Chemical Substances (Japan); ErCx - Concentration associated with x% growth rate response; ERG - Emergency Response Guide; GHS - Globally Harmonized System; GLP - Good Laboratory Practice; IARC - International Agency for Research on Cancer; IATA - International Air Transport Association; IBC - International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk; IC50 - Half maximal inhibitory concentration; ICAO - International Civil Aviation Organization; IECSC - Inventory of Existing Chemical Substances in China; IMDG - International Maritime Dangerous Goods; IMO - International Maritime Organization; ISHL - Industrial Safety and Health Law (Japan); ISO - International Organisation for Standardization; KECI - Korea Existing Chemicals Inventory; LC50 - Lethal Concentration to 50 % of a test population; LD50 - Lethal Dose to 50% of a test population (Median Lethal Dose); MARPOL - International Convention for the Prevention of Pollution from Ships; n.o.s. - Not Otherwise Specified; Nch - Chilean Norm; NO(A)EC -No Observed (Adverse) Effect Concentration; NO(A)EL - No Observed (Adverse) Effect Level; NOELR - No Observable Effect Loading Rate; NOM - Official Mexican Norm; NTP - National Toxicology Program; NZIoC - New Zealand Inventory of Chemicals; OECD - Organization for Economic Co-operation and Development; OPPTS - Office of Chemical Safety and Pollution Prevention; PBT - Persistent, Bioaccumulative and Toxic substance; PICCS - Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR - (Quantitative) Structure Activity Relationship; REACH - Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; SADT - Self-Accelerating Decomposition Temperature; SDS - Safety Data Sheet; TCSI - Taiwan Chemical Substance Inventory; TDG - Transportation of Dangerous Goods; TSCA - Toxic Substances Control Act (United States); UN - United Nations; UNRTDG - United Nations Recommendations on the Transport of Dangerous Goods; vPvB - Very Persistent and Very Bioaccumulative; WHMIS -Workplace Hazardous Materials Information System

Date format : dd.mm.yyyy Numbers 123,456.78

SAFETY DATA SHEET CellSave Preservative 20 Tubes



 Version
 Revision Date:
 SDS Number:
 Date of last issue: 27.10.2016

 1.11
 09.01.2018
 100000011150
 Date of first issue: 18.09.2015

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

AU / EN