according to Regulation (EC) No. 1907/2006



CellSave Preservative 20 Tubes

Version Revision Date: SDS Number: Date of last issue: 2017-08-02 1.44 2018-01-09 100000011150 Date of first issue: 2015-02-04

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1 Product identifier

Trade name : CellSave Preservative 20 Tubes Substance name : CellSave Preservative 20 Tubes

1.2 Relevant identified uses of the substance or mixture and uses advised against

Use of the Sub-

stance/Mixture

: Assay reagent

1.3 Details of the supplier of the safety data sheet

Company : Menarini Silicon Biosystems, Inc.

3401 Masons Mill Rd #100 Huntingdon Valley, PA 19006,

USA

Telephone : 1 (800) 381-4929

Telefax

E-mail address : Us-info@siliconbiosystems.com

Responsible/issuing person

1.4 Emergency telephone number

US: (303) 389-1805

International: +1 (303) 389-1805

SECTION 2: Hazards identification

2.1 Classification of the substance or mixture

Classification (REGULATION (EC) No 1272/2008)

Skin sensitisation, Category 1 H317: May cause an allergic skin reaction.

2.2 Label elements

Labelling (REGULATION (EC) No 1272/2008)

Hazard pictograms

Signal word : Warning

Hazard statements : H317 May cause an allergic skin reaction.

Precautionary statements : Prevention:

according to Regulation (EC) No. 1907/2006



| CellSave Preservative 20 Tubes | | | | | | |
|--------------------------------|------------------------------|--------------------------|--|--|--|--|
| Version 1.44 | Revision Date: 2018-01-09 | SDS Number: 100000011150 | Date of last issue: 2017-08-02 Date of first issue: 2015-02-04 | | | |
| | | P261 | Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray. | | | |
| | | P272 | Contaminated work clothing should not be allowed out of the workplace. | | | |
| | | P280 | Wear protective gloves. | | | |
| | | Response: | | | | |
| | | P333 + P313 | If skin irritation or rash occurs: Get medical advice/ attention. | | | |
| | | P362 + P364 | Take off contaminated clothing and wash it before reuse. | | | |
| | | Disposal: | | | | |
| | | P501 | Dispose of contents/ container to an approved waste disposal plant. | | | |

Hazardous components which must be listed on the label: IMIDUREA

2.3 Other hazards

This substance/mixture contains no components considered to be either persistent, bioaccumulative and toxic (PBT), or very persistent and very bioaccumulative (vPvB) at levels of 0.1% or higher.

SECTION 3: Composition/information on ingredients

3.2 Mixtures

Chemical nature : Liquid

Hazardous components

| Chemical name | CAS-No. EC-No. Registration number | Classification (REGULATION (EC) No 1272/2008) | Concentration (%) |
|---------------|---|--|-------------------|
| IMIDUREA | 39236-46-9 254-372-6 | Skin Sens. 1; H317 | >= 25 - < 50 |

For explanation of abbreviations see section 16.

SECTION 4: First aid measures

4.1 Description of first aid measures

If inhaled : If breathed in, move person into fresh air.

Consult a physician.

In case of skin contact : Take off contaminated clothing and shoes immediately.

Wash off immediately with plenty of water. If symptoms persist, call a physician. Wash contaminated clothing before re-use.

according to Regulation (EC) No. 1907/2006



CellSave Preservative 20 Tubes

Date of last issue: 2017-08-02 Version Revision Date: SDS Number: 1.44 2018-01-09 100000011150 Date of first issue: 2015-02-04

Rinse immediately with plenty of water, also under the eyelids, In case of eye contact

> for at least 5 minutes. Remove contact lenses.

If eye irritation persists, consult a specialist.

If swallowed If swallowed, rinse mouth with water (only if the person is con-

scious).

Call a physician immediately.

4.2 Most important symptoms and effects, both acute and delayed

Symptoms : No information available.

4.3 Indication of any immediate medical attention and special treatment needed

Treatment : Treat symptomatically.

SECTION 5: Firefighting measures

5.1 Extinguishing media

Suitable extinguishing media Use extinguishing measures that are appropriate to local cir-

cumstances and the surrounding environment.

5.2 Special hazards arising from the substance or mixture

Specific hazards during fire-

fighting

: No information available.

5.3 Advice for firefighters

for firefighters

Special protective equipment : In the event of fire, wear self-contained breathing apparatus.

Further information : No information available.

SECTION 6: Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures

Personal precautions : In the event of an accidental release the emergency response

team must respond based on a risk assessment and use per-

sonal protective equipment as appropriate.

6.2 Environmental precautions

: Should not be released into the environment. Environmental precautions

6.3 Methods and material for containment and cleaning up

Methods for cleaning up : Large spills: Dam up. Soak up with inert absorbent material.

Keep in properly labelled containers.

according to Regulation (EC) No. 1907/2006



CellSave Preservative 20 Tubes

Version Revision Date: SDS Number: Date of last issue: 2017-08-02 1.44 2018-01-09 100000011150 Date of first issue: 2015-02-04

Small spills: Gently cover the spill with an absorbent towel or

pad.

Large spills + Small spills: Keep in suitable, closed containers for disposal. Treat recovered material as described in the sec-

tion "Disposal considerations".

6.4 Reference to other sections

For disposal information, see section 13

SECTION 7: Handling and storage

7.1 Precautions for safe handling

Advice on safe handling : To avoid thermal decomposition, do not overheat.

Avoid inhalation, ingestion and contact with skin and eyes.

Use personal protective equipment as required.

Advice on protection against

fire and explosion

: No data available

Hygiene measures : Handle in accordance with good industrial hygiene and safety

practice.

7.2 Conditions for safe storage, including any incompatibilities

Requirements for storage areas and containers

: To maintain product quality, do not store in heat or direct sunlight. Store in original container. Keep containers tightly closed in a dry, cool and well-ventilated place. Keep away from heat and sources of ignition. Keep locked up. Keep refrigerated.

Recommended storage tem-

perature

: 4 - 30 °C

7.3 Specific end use(s)

Specific use(s) : Consult the technical guidelines for the use of this sub-

stance/mixture.

SECTION 8: Exposure controls/personal protection

8.1 Control parameters

Contains no substances with occupational exposure limit values.

8.2 Exposure controls

Engineering measures

All personal protective equipment should be based on a risk assessment. Consult a Environment Health Safety expert if necessary.

Personal protective equipment

Eye protection : No special precautions required.

according to Regulation (EC) No. 1907/2006



CellSave Preservative 20 Tubes

 Version
 Revision Date:
 SDS Number:
 Date of last issue: 2017-08-02

 1.44
 2018-01-09
 100000011150
 Date of first issue: 2015-02-04

Hand protection

Remarks : Disposable gloves

Skin and body protection : No special precautions required.

Respiratory protection : Engineering controls should always be the primary method of

controlling exposures.

If respiratory protective equipment is needed for certain activities, the type as well as the corresponding protection factor will depend upon the risk assessment and air concentrations, hazards, physical and warning properties of substances pre-

sent.

No personal respiratory protective equipment normally re-

quired.

Protective measures : The type of protective equipment must be selected based on

the Environmental Health and Safety risk assessment. Consult a Environmental Health and Safety expert if necessary.

SECTION 9: Physical and chemical properties

9.1 Information on basic physical and chemical properties

Appearance : liquid

Colour : clear

9.2 Other information

No data available

SECTION 10: Stability and reactivity

10.1 Reactivity

None reasonably foreseeable.

10.2 Chemical stability

Stable under recommended storage conditions.

10.3 Possibility of hazardous reactions

Hazardous reactions : No dangerous reaction known under conditions of normal use.

10.4 Conditions to avoid

Conditions to avoid : To avoid thermal decomposition, do not overheat.

10.5 Incompatible materials

Materials to avoid : None known.

according to Regulation (EC) No. 1907/2006



CellSave Preservative 20 Tubes

Version Revision Date: SDS Number: Date of last issue: 2017-08-02 1.44 2018-01-09 100000011150 Date of first issue: 2015-02-04

10.6 Hazardous decomposition products

None known.

SECTION 11: Toxicological information

11.1 Information on toxicological effects

Acute toxicity

Components:

IMIDUREA:

Acute oral toxicity : LD50 (Rat): 11 300 mg/kg

Acute inhalation toxicity : LC50 (Rat): > 5,5 mg/l

Exposure time: 1 h

Acute dermal toxicity : LD50 (Rabbit): > 5 000 mg/kg

Skin corrosion/irritation

Components:

IMIDUREA:

Result: No skin irritation

Serious eye damage/eye irritation

Components:

IMIDUREA:

Result: No eye irritation

Respiratory or skin sensitisation

Components:

IMIDUREA:

Method: Maximisation Test

Result: May cause sensitisation by skin contact.

Method: Local Lymph Node Assay (LLNA) in mice Result: May cause sensitisation by skin contact.

Germ cell mutagenicity

Components:

IMIDUREA:

Genotoxicity in vitro : Test Type: Ames test

Metabolic activation: with and without metabolic activation

Result: negative GLP: yes

: Test Type: Chromosome aberration test in vitro

Metabolic activation: with and without metabolic activation

according to Regulation (EC) No. 1907/2006



CellSave Preservative 20 Tubes

 Version
 Revision Date:
 SDS Number:
 Date of last issue: 2017-08-02

 1.44
 2018-01-09
 100000011150
 Date of first issue: 2015-02-04

Result: negative

Genotoxicity in vivo : Test Type: Micronucleus test

Species: Mouse Application Route: Oral Result: negative

Germ cell mutagenicity- As-

sessment

: No information available.

Carcinogenicity

Components:

IMIDUREA:

Carcinogenicity - Assess-

ment

: No information available.

Reproductive toxicity

Components:

IMIDUREA:

Teratogenicity - Assessment : No information available.

STOT - single exposure

No data available

STOT - repeated exposure

No data available

Repeated dose toxicity

Components:

IMIDUREA:

Species: Rat NOAEL: 200 mg/kg LOAEL: 500 mg/kg Application Route: Oral

Species: Rabbit NOAEL: 200 mg/kg Application Route: Dermal

Aspiration toxicity

No data available

SECTION 12: Ecological information

12.1 Toxicity

Components:

IMIDUREA:

according to Regulation (EC) No. 1907/2006



CellSave Preservative 20 Tubes

 Version
 Revision Date:
 SDS Number:
 Date of last issue: 2017-08-02

 1.44
 2018-01-09
 100000011150
 Date of first issue: 2015-02-04

Toxicity to fish : Remarks: No data available

12.2 Persistence and degradability

Components: IMIDUREA:

Biodegradability : Remarks: No data available

12.3 Bioaccumulative potential

Components: IMIDUREA:

Bioaccumulation : Remarks: No data available

12.4 Mobility in soil

Components: IMIDUREA:

Distribution among environmental compartments

: Remarks: No data available

12.5 Results of PBT and vPvB assessment

Product:

Assessment : This substance/mixture contains no components considered

to be either persistent, bioaccumulative and toxic (PBT), or very persistent and very bioaccumulative (vPvB) at levels of

0.1% or higher...

Components: IMIDUREA:

Assessment : No information available..

12.6 Other adverse effects

Components: IMIDUREA:

Additional ecological infor-

mation

: Remarks: No data available

SECTION 13: Disposal considerations

13.1 Waste treatment methods

Product : In accordance with National, Federal, State and Local regula-

tions.

Contaminated packaging : Empty containers should be taken to an approved waste han-

according to Regulation (EC) No. 1907/2006



CellSave Preservative 20 Tubes

Version Revision Date: SDS Number: Date of last issue: 2017-08-02 1.44 2018-01-09 100000011150 Date of first issue: 2015-02-04

dling site for recycling or disposal.

SECTION 14: Transport information

14.1 UN number

Not regulated as a dangerous good

14.2 UN proper shipping name

Not regulated as a dangerous good

14.3 Transport hazard class(es)

Not regulated as a dangerous good

14.4 Packing group

Not regulated as a dangerous good

14.5 Environmental hazards

Not regulated as a dangerous good

14.6 Special precautions for user

Not applicable

14.7 Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable for product as supplied.

SECTION 15: Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

Seveso III: Directive 2012/18/EU of the European Parliament and of the Council on the control of major-accident hazards involving dangerous substances.

Not applicable

Other regulations : Restricted to professional users.

15.2 Chemical safety assessment

A Chemical Safety Assessment is not applicable (mixture)

SECTION 16: Other information

Full text of H-Statements

H317 : May cause an allergic skin reaction.

Date and Number Formats

This document uses the following notation for printing dates and numbers:

Date: Dec 31th, 2012 as 2012-12-31

according to Regulation (EC) No. 1907/2006



CellSave Preservative 20 Tubes

 Version
 Revision Date:
 SDS Number:
 Date of last issue: 2017-08-02

 1.44
 2018-01-09
 100000011150
 Date of first issue: 2015-02-04

Numbers: 123456,78 as 123 456,78

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

BE / EN