Capture enhancement reagent



Version Revision Date: SDS Number: Date of last issue: 2016-08-05 1.11 2016-10-27 100000010878 Date of first issue: 2015-07-14

1. PRODUCT AND COMPANY IDENTIFICATION

Product name : Capture enhancement reagent

Substance name : Capture enhancement reagent

7037

Manufacturer or supplier's details

Company : Janssen Diagnostics, LLC

Address : 700 US Highway Route 202

South Raritan, NJ 08869

US

Telephone : (877) 837-4339

Emergency telephone : CHEMTREC IN: 000-800-100-7141

number CHEMTREC International: +1 703-527-3887

E-mail address : SDSJanssen@its.jnj.com

Responsible/issuing person

Recommended use of the chemical and restrictions on use

Recommended use : Assay reagent

2. HAZARDS IDENTIFICATION

Manufacture, Storage and Import of Hazardous Chemicals Rules 1989

Classification

GHS Classification

Not a hazardous substance or mixture.

GHS label elements

Not a hazardous substance or mixture.

Other hazards which do not result in classification

None known.

3. COMPOSITION/INFORMATION ON INGREDIENTS

Substance / Mixture : Mixture

Chemical nature : Liquid

Hazardous components

Capture enhancement reagent



 Version
 Revision Date:
 SDS Number:
 Date of last issue: 2016-08-05

 1.11
 2016-10-27
 100000010878
 Date of first issue: 2015-07-14

Chemical name	CAS-No. Concentration	
		w/w)
sodium azide	26628-22-8	>= 0,1 - < 1

4. FIRST AID MEASURES

If inhaled : If breathed in, move person into fresh air.

Consult a physician.

In case of skin contact : Take off contaminated clothing and shoes immediately.

Wash off with plenty of water.

If symptoms persist, call a physician.

In case of eye contact : Rinse immediately with plenty of water, also under the eyelids,

for at least 5 minutes. Remove contact lenses.

If eye irritation persists, consult a specialist.

If swallowed, rinse mouth with water (only if the person is con-

scious)

Call a physician immediately.

Most important symptoms and effects, both acute and

delayed

No information available.

Notes to physician : Treat symptomatically.

5. FIREFIGHTING MEASURES

Suitable extinguishing media : Use extinguishing measures that are appropriate to local cir-

cumstances and the surrounding environment.

Specific hazards during fire-

fighting

No information available.

Specific extinguishing meth-

ods

No information available.

Special protective equipment:

for firefighters

In the event of fire, wear self-contained breathing apparatus.

6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment and emer-

gency procedures

: In the event of an accidental release the emergency response team must respond based on a risk assessment and use per-

sonal protective equipment as appropriate.

Environmental precautions : Should not be released into the environment.

Methods and materials for containment and cleaning up

: Large spills: Dam up. Soak up with inert absorbent material.

Keep in properly labelled containers.

Small spills: Gently cover the spill with an absorbent towel or

Capture enhancement reagent



Version 1.11

Revision Date: 2016-10-27

SDS Number: 100000010878

Date of last issue: 2016-08-05 Date of first issue: 2015-07-14

pad.

Large spills + Small spills: Keep in suitable, closed containers for disposal. Treat recovered material as described in the sec-

tion "Disposal considerations".

7. HANDLING AND STORAGE

Advice on protection against

fire and explosion

No data available

Advice on safe handling :

To avoid thermal decomposition, do not overheat.

Avoid inhalation, ingestion and contact with skin and eyes.

Use personal protective equipment as required.

Conditions for safe storage

To maintain product quality, do not store in heat or direct sun-

light.

Store in original container.

Keep containers tightly closed in a dry, cool and well-

ventilated place.

Keep away from heat and sources of ignition.

Keep locked up. Keep refrigerated.

Recommended storage tem: :

perature

2 - 8 °C

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Components with workplace control parameters

Components	CAS-No.	Value type (Form of exposure)	Control parameters / Permissible concentration	Basis
sodium azide	26628-22-8	C (Vapour)	0,11 ppm (Hydrazoic acid)	ACGIH
		С	0,29 mg/m3 (Sodium azide)	ACGIH

Engineering measures

All personal protective equipment should be based on a risk assessment. Consult a Environment Health Safety expert if

necessary.

Personal protective equipment

Respiratory protection

Engineering controls should always be the primary method of

controlling exposures.

If respiratory protective equipment is needed for certain activities, the type as well as the corresponding protection factor will depend upon the risk assessment and air concentrations, hazards, physical and warning properties of substances pre-

sent.

No personal respiratory protective equipment normally re-

quired.

Capture enhancement reagent



Version Revision Date: SDS Number: Date of last issue: 2016-08-05 1.11 2016-10-27 100000010878 Date of first issue: 2015-07-14

Hand protection

Remarks : Disposable gloves

Eye protection : No special precautions required.

Skin and body protection : No special precautions required.

Protective measures : The type of protective equipment must be selected based on

the Environmental Health and Safety risk assessment. Consult a Environmental Health and Safety expert if necessary.

Hygiene measures : Handle in accordance with good industrial hygiene and safety

practice.

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance : liquid

Colour : clear

Odour : No data available

10. STABILITY AND REACTIVITY

Reactivity : None reasonably foreseeable.

Chemical stability : Stable under recommended storage conditions.

Possibility of hazardous reac-

tions

No dangerous reaction known under conditions of normal use.

Conditions to avoid : To avoid thermal decomposition, do not overheat.

Incompatible materials : Oxidizing agents

Hazardous decomposition

products

None known.

11. TOXICOLOGICAL INFORMATION

Acute toxicity

Product:

Acute oral toxicity : Acute toxicity estimate: > 5 000 mg/kg

Method: Calculation method

Components:

sodium azide:

Acute oral toxicity : LD50 (Rat): 27 mg/kg

Capture enhancement reagent



Version 1.11

Revision Date: 2016-10-27

SDS Number: 100000010878

Date of last issue: 2016-08-05 Date of first issue: 2015-07-14

Skin corrosion/irritation

No data available

Serious eye damage/eye irritation

No data available

Respiratory or skin sensitisation

No data available

Germ cell mutagenicity

No data available

Carcinogenicity

No data available

Reproductive toxicity

No data available

STOT - single exposure

No data available

STOT - repeated exposure

No data available

Repeated dose toxicity

No data available

Aspiration toxicity

No data available

12. ECOLOGICAL INFORMATION

Ecotoxicity

Components:

sodium azide:

Toxicity to fish : LC50 (Lepomis macrochirus (Bluegill sunfish)): 0,7 mg/l

Exposure time: 96 h

Toxicity to daphnia and other :

aquatic invertebrates

EC50 (Daphnia pulex (Water flea)): 4,2 mg/l

Exposure time: 96 h

Toxicity to algae : IC50: 272 mg/l

Toxicity to bacteria : EC50 (Photobacterium phosphoreum): 38,5 mg/l

Persistence and degradability

No data available

Bioaccumulative potential

No data available

Capture enhancement reagent



Version 1.11

Revision Date: 2016-10-27

SDS Number: 100000010878

Date of last issue: 2016-08-05 Date of first issue: 2015-07-14

Mobility in soil

No data available

Other adverse effects

No data available

13. DISPOSAL CONSIDERATIONS

Disposal methods

Waste from residues

In accordance with National, Federal, State and Local regula-

tions.

Contaminated packaging

Empty containers should be taken to an approved waste han-

dling site for recycling or disposal.

14. TRANSPORT INFORMATION

International Regulations

UNRTDG

Not regulated as a dangerous good

IATA-DGR

Not regulated as a dangerous good

IMDG-Code

Not regulated as a dangerous good

Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable for product as supplied.

15. REGULATORY INFORMATION

Safety, health and environmental regulations/legislation specific for the substance or mixture

Restricted to professional users.

16. OTHER INFORMATION

Full text of other abbreviations

AICS - Australian Inventory of Chemical Substances; ANTT - National Agency for Transport by Land of Brazil; ASTM - American Society for the Testing of Materials; bw - Body weight; CMR - Carcinogen, Mutagen or Reproductive Toxicant; CPR - Controlled Products Regulations; DIN - Standard of the German Institute for Standardisation; DSL - Domestic Substances List (Canada); ECx - Concentration associated with x% response; ELx - Loading rate associated with x% response; EmS - Emergency Schedule; ENCS - Existing and New Chemical Substances (Japan); ErCx - Concentration associated with x% growth rate response; ERG - Emergency Response Guide; GHS - Globally Harmonized System; GLP - Good Laboratory Practice; IARC - International Agency for Research on Cancer; IATA - International Air Transport Association; IBC - International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk;

Capture enhancement reagent



Version Revision Date: SDS Number: Date of last issue: 2016-08-05 1.11 2016-10-27 100000010878 Date of first issue: 2015-07-14

IC50 - Half maximal inhibitory concentration; ICAO - International Civil Aviation Organization; IECSC - Inventory of Existing Chemical Substances in China; IMDG - International Maritime Dangerous Goods; IMO - International Maritime Organization; ISHL - Industrial Safety and Health Law (Japan); ISO - International Organisation for Standardization; KECI - Korea Existing Chemicals Inventory; LC50 - Lethal Concentration to 50 % of a test population; LD50 - Lethal Dose to 50% of a test population (Median Lethal Dose); MARPOL - International Convention for the Prevention of Pollution from Ships; n.o.s. - Not Otherwise Specified; Nch - Chilean Norm; NO(A)EC -No Observed (Adverse) Effect Concentration; NO(A)EL - No Observed (Adverse) Effect Level; NOELR - No Observable Effect Loading Rate; NOM - Official Mexican Norm; NTP - National Toxicology Program; NZIoC - New Zealand Inventory of Chemicals; OECD - Organization for Economic Co-operation and Development; OPPTS - Office of Chemical Safety and Pollution Prevention; PBT - Persistent, Bioaccumulative and Toxic substance; PICCS - Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR - (Quantitative) Structure Activity Relationship; REACH - Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; SADT - Self-Accelerating Decomposition Temperature; SDS - Safety Data Sheet; TCSI - Taiwan Chemical Substance Inventory; TDG - Transportation of Dangerous Goods; TSCA - Toxic Substances Control Act (United States); UN - United Nations; UNRTDG - United Nations Recommendations on the Transport of Dangerous Goods; vPvB - Very Persistent and Very Bioaccumulative; WHMIS -Workplace Hazardous Materials Information System

Date format : dd.mm.yyyy

Numbers : 123 456,78

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

IN / EN



Version Revision Date: SDS Number: Date of last issue: -

1.0 2016-11-15 100000013033 Date of first issue: 2016-11-15

1. PRODUCT AND COMPANY IDENTIFICATION

Substance name : Staining Reagent

Manufacturer or supplier's details

Company : Janssen Diagnostics, LLC

Address : 700 US Highway Route 202

South Raritan, NJ 08869

US

Telephone : (877) 837-4339

Emergency telephone : CHEMTREC IN: 000-800-100-7141

number CHEMTREC International: +1 703-527-3887

E-mail address : SDSJanssen@its.jnj.com

Responsible/issuing person

Recommended use of the chemical and restrictions on use

Recommended use : Assay reagent

2. HAZARDS IDENTIFICATION

Manufacture, Storage and Import of Hazardous Chemicals Rules 1989

Classification

GHS Classification

Not a hazardous substance or mixture. Not a hazardous substance or mixture.

GHS label elements

Not a hazardous substance or mixture. Not a hazardous substance or mixture.

Other hazards which do not result in classification

None known.

3. COMPOSITION/INFORMATION ON INGREDIENTS

Substance / Mixture : Mixture

Chemical nature : Liquid

Hazardous components



Version Revision Date: SDS Number: Date of last issue: -

1.0 2016-11-15 100000013033 Date of first issue: 2016-11-15

Chemical name	CAS-No. Concen	
sodium azide	26628-22-8	w/w) >= 0,1 - < 1

4. FIRST AID MEASURES

If inhaled : If breathed in, move person into fresh air.

Consult a physician.

In case of skin contact : Take off contaminated clothing and shoes immediately.

Wash off with plenty of water.

If symptoms persist, call a physician.

In case of eye contact : Rinse immediately with plenty of water, also under the eyelids,

for at least 5 minutes. Remove contact lenses.

If eye irritation persists, consult a specialist.

If swallowed : If swallowed, rinse mouth with water (only if the person is con-

scious).

Call a physician immediately.

Most important symptoms and effects, both acute and

delayed

No information available.

Notes to physician : Treat symptomatically.

5. FIREFIGHTING MEASURES

Suitable extinguishing media : Use extinguishing measures that are appropriate to local cir-

cumstances and the surrounding environment.

Specific hazards during fire-

fighting

No information available.

Hazardous combustion prod- :

ucts

No hazardous combustion products are known

Specific extinguishing meth-

ods

No information available.

Special protective equipment :

for firefighters

In the event of fire, wear self-contained breathing apparatus.

6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment and emer-

gency procedures

In the event of an accidental release the emergency response team must respond based on a risk assessment and use per-

sonal protective equipment as appropriate.

Environmental precautions : Should not be released into the environment.



Version Revision Date: SDS Number: Date of last issue: -

1.0 2016-11-15 100000013033 Date of first issue: 2016-11-15

Methods and materials for containment and cleaning up

Large spills: Dam up. Soak up with inert absorbent material.

Keep in properly labelled containers.

Small spills: Gently cover the spill with an absorbent towel or

pad.

Large spills + Small spills: Keep in suitable, closed containers for disposal. Treat recovered material as described in the sec-

tion "Disposal considerations".

7. HANDLING AND STORAGE

Advice on protection against

fire and explosion

No data available

Advice on safe handling : To avoid thermal decomposition, do not overheat.

Avoid inhalation, ingestion and contact with skin and eyes.

Use personal protective equipment as required.

Conditions for safe storage : To maintain product quality, do not store in heat or direct sun-

light.

Store in original container.

Keep containers tightly closed in a dry, cool and well-

ventilated place.

Keep away from heat and sources of ignition.

Keep locked up. Keep refrigerated.

Recommended storage tem- :

perature

2 - 8 °C

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Components with workplace control parameters

Components	CAS-No.	Value type	Control parame-	Basis
		(Form of	ters / Permissible	
		exposure)	concentration	
sodium azide	26628-22-8	C (Vapour)	0,11 ppm	ACGIH
			(Hydrazoic acid)	
		С	0,29 mg/m3	ACGIH
			(Sodium azide)	

Engineering measures : All personal protective equipment should be based on a risk

assessment. Consult a Environment Health Safety expert if

necessary.

Personal protective equipment

Respiratory protection : Engineering controls should always be the primary method of

controlling exposures.

If respiratory protective equipment is needed for certain activities, the type as well as the corresponding protection factor will depend upon the risk assessment and air concentrations, hazards, physical and warning properties of substances pre-

sent.



Version Revision Date: SDS Number: Date of last issue: -

1.0 2016-11-15 100000013033 Date of first issue: 2016-11-15

No personal respiratory protective equipment normally re-

quired.

Eye protection : No special precautions required.

Skin and body protection : No special precautions required.

Protective measures : The type of protective equipment must be selected based on

the Environmental Health and Safety risk assessment. Consult a Environmental Health and Safety expert if necessary.

Hygiene measures : Handle in accordance with good industrial hygiene and safety

practice.

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance : liquid

Colour : purple

Odour : odourless

pH : 7,5

Solubility(ies)

Water solubility : soluble

10. STABILITY AND REACTIVITY

Reactivity : None reasonably foreseeable.

Chemical stability : Stable under recommended storage conditions.

Possibility of hazardous reac-

tions

No dangerous reaction known under conditions of normal use.

Conditions to avoid : To avoid thermal decomposition, do not overheat.

Incompatible materials : Strong acids and strong bases

Reducing agents Oxidizing agents

Hazardous decomposition

products

None known.

11. TOXICOLOGICAL INFORMATION

Acute toxicity

Product:

Acute oral toxicity : Acute toxicity estimate: > 5 000 mg/kg



Version Revision Date: SDS Number: Date of last issue: -

1.0 2016-11-15 100000013033 Date of first issue: 2016-11-15

Method: Calculation method

Acute toxicity estimate: > 5 000 mg/kg

Method: Calculation method

Components:

sodium azide:

Acute oral toxicity : LD50 (Rat): 27 mg/kg

Skin corrosion/irritation

No data available

Serious eye damage/eye irritation

No data available

Respiratory or skin sensitisation

No data available

Germ cell mutagenicity

No data available

Carcinogenicity

No data available

Reproductive toxicity

No data available

STOT - single exposure

No data available

STOT - repeated exposure

No data available

Repeated dose toxicity

No data available

Aspiration toxicity

No data available

12. ECOLOGICAL INFORMATION

Ecotoxicity

Components:

sodium azide:

Toxicity to fish : LC50 (Lepomis macrochirus (Bluegill sunfish)): 0,7 mg/l

Exposure time: 96 h

Toxicity to daphnia and other :

aquatic invertebrates

EC50 (Daphnia pulex (Water flea)): 4,2 mg/l

Exposure time: 96 h

Toxicity to algae : IC50: 272 mg/l



Version Revision Date: SDS Number: Date of last issue: -

1.0 2016-11-15 100000013033 Date of first issue: 2016-11-15

Toxicity to bacteria : EC50 (Photobacterium phosphoreum): 38,5 mg/l

Persistence and degradability

No data available

Bioaccumulative potential

No data available

Mobility in soil

No data available

Other adverse effects

No data available

13. DISPOSAL CONSIDERATIONS

Disposal methods

Waste from residues : In accordance with National, Federal, State and Local regula-

tions.

Contaminated packaging : Empty containers should be taken to an approved waste han-

dling site for recycling or disposal.

14. TRANSPORT INFORMATION

International Regulations

UNRTDG

Not regulated as a dangerous good

IATA-DGR

Not regulated as a dangerous good

IMDG-Code

Not regulated as a dangerous good

Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable for product as supplied.

15. REGULATORY INFORMATION

Safety, health and environmental regulations/legislation specific for the substance or mixture

Restricted to professional users.



Version Revision Date: SDS Number: Date of last issue: -

1.0 2016-11-15 100000013033 Date of first issue: 2016-11-15

16. OTHER INFORMATION

Full text of other abbreviations

AICS - Australian Inventory of Chemical Substances; ANTT - National Agency for Transport by Land of Brazil; ASTM - American Society for the Testing of Materials; bw - Body weight; CMR -Carcinogen, Mutagen or Reproductive Toxicant; CPR - Controlled Products Regulations; DIN -Standard of the German Institute for Standardisation: DSL - Domestic Substances List (Canada): ECx - Concentration associated with x% response; ELx - Loading rate associated with x% response; EmS - Emergency Schedule; ENCS - Existing and New Chemical Substances (Japan); ErCx - Concentration associated with x% growth rate response; ERG - Emergency Response Guide; GHS - Globally Harmonized System; GLP - Good Laboratory Practice; IARC - International Agency for Research on Cancer; IATA - International Air Transport Association; IBC - International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk; IC50 - Half maximal inhibitory concentration; ICAO - International Civil Aviation Organization; IECSC - Inventory of Existing Chemical Substances in China; IMDG - International Maritime Dangerous Goods; IMO - International Maritime Organization; ISHL - Industrial Safety and Health Law (Japan): ISO - International Organisation for Standardization: KECI - Korea Existing Chemicals Inventory; LC50 - Lethal Concentration to 50 % of a test population; LD50 - Lethal Dose to 50% of a test population (Median Lethal Dose); MARPOL - International Convention for the Prevention of Pollution from Ships; n.o.s. - Not Otherwise Specified; Nch - Chilean Norm; NO(A)EC -No Observed (Adverse) Effect Concentration; NO(A)EL - No Observed (Adverse) Effect Level; NOELR - No Observable Effect Loading Rate; NOM - Official Mexican Norm; NTP - National Toxicology Program; NZIoC - New Zealand Inventory of Chemicals; OECD - Organization for Economic Co-operation and Development; OPPTS - Office of Chemical Safety and Pollution Prevention; PBT - Persistent, Bioaccumulative and Toxic substance; PICCS - Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR - (Quantitative) Structure Activity Relationship; REACH - Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; SADT - Self-Accelerating Decomposition Temperature; SDS - Safety Data Sheet; TCSI - Taiwan Chemical Substance Inventory; TDG - Transportation of Dangerous Goods; TSCA - Toxic Substances Control Act (United States); UN - United Nations; UNRTDG - United Nations Recommendations on the Transport of Dangerous Goods; vPvB - Very Persistent and Very Bioaccumulative; WHMIS -Workplace Hazardous Materials Information System

Date format : dd.mm.yyyy

Numbers : 123 456,78

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

IN / EN

Nucleic acid dye



Version Revision Date: SDS Number: Date of last issue: 2016-08-05 1.9 2016-10-27 100000010877 Date of first issue: 2015-07-14

1. PRODUCT AND COMPANY IDENTIFICATION

Product name : Nucleic acid dye

Substance name : Nucleic acid dye

7041

Manufacturer or supplier's details

Company : Janssen Diagnostics, LLC

Address : 700 US Highway Route 202

South Raritan, NJ 08869

US

Telephone : (877) 837-4339

Emergency telephone : CHEMTREC IN: 000-800-100-7141

number CHEMTREC International: +1 703-527-3887

E-mail address : SDSJanssen@its.jnj.com

Responsible/issuing person

Recommended use of the chemical and restrictions on use

Recommended use : Assay reagent

2. HAZARDS IDENTIFICATION

Manufacture, Storage and Import of Hazardous Chemicals Rules 1989

Classification

GHS Classification

Not a hazardous substance or mixture.

GHS label elements

Not a hazardous substance or mixture.

Other hazards which do not result in classification

None known.

3. COMPOSITION/INFORMATION ON INGREDIENTS

Substance / Mixture : Mixture

Chemical nature : Liquid

Nucleic acid dye



Version 1.9

Revision Date: 2016-10-27

SDS Number: 100000010877 Date of last issue: 2016-08-05 Date of first issue: 2015-07-14

Hazardous components

No hazardous ingredients

4. FIRST AID MEASURES

If inhaled If breathed in, move person into fresh air.

Consult a physician.

Take off contaminated clothing and shoes immediately. In case of skin contact

> Wash off immediately with plenty of water. If symptoms persist, call a physician. Wash contaminated clothing before re-use.

In case of eye contact Rinse immediately with plenty of water, also under the eyelids,

> for at least 5 minutes. Remove contact lenses.

If eye irritation persists, consult a specialist.

If swallowed If swallowed, rinse mouth with water (only if the person is con-

Call a physician immediately.

Most important symptoms and effects, both acute and

delayed

No information available.

Notes to physician Treat symptomatically.

5. FIREFIGHTING MEASURES

Suitable extinguishing media Use extinguishing measures that are appropriate to local cir-

cumstances and the surrounding environment.

Specific hazards during fire-

fighting

No information available.

Specific extinguishing meth-

No information available.

Special protective equipment:

for firefighters

In the event of fire, wear self-contained breathing apparatus.

6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protec- : tive equipment and emer-

gency procedures

In the event of an accidental release the emergency response team must respond based on a risk assessment and use per-

sonal protective equipment as appropriate.

Should not be released into the environment. **Environmental precautions**

Methods and materials for containment and cleaning up Large spills: Dam up. Soak up with inert absorbent material.

Keep in properly labelled containers.

Small spills: Gently cover the spill with an absorbent towel or

Nucleic acid dye



Version 1.9

Revision Date: 2016-10-27

SDS Number: 100000010877

Date of last issue: 2016-08-05 Date of first issue: 2015-07-14

pad.

Large spills + Small spills: Keep in suitable, closed containers for disposal. Treat recovered material as described in the sec-

tion "Disposal considerations".

7. HANDLING AND STORAGE

Advice on protection against

fire and explosion

No data available

Advice on safe handling

To avoid thermal decomposition, do not overheat.

Avoid inhalation, ingestion and contact with skin and eyes.

Use personal protective equipment as required.

Conditions for safe storage

To maintain product quality, do not store in heat or direct sun-

light.

Store in original container.

Keep containers tightly closed in a dry, cool and well-

ventilated place.

Keep away from heat and sources of ignition.

Keep locked up. Keep refrigerated.

Recommended storage tem: :

perature

2 - 8 °C

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Components with workplace control parameters

Contains no substances with occupational exposure limit values.

Engineering measures

All personal protective equipment should be based on a risk assessment. Consult a Environment Health Safety expert if

necessary.

Personal protective equipment

Respiratory protection

Engineering controls should always be the primary method of

controlling exposures.

If respiratory protective equipment is needed for certain activities, the type as well as the corresponding protection factor will depend upon the risk assessment and air concentrations, hazards, physical and warning properties of substances pre-

sent.

No personal respiratory protective equipment normally re-

quired.

Hand protection

Remarks : Disposable gloves

Eye protection : No special precautions required.

Skin and body protection : No special precautions required.

Nucleic acid dye



Version Revision Date: SDS Number: Date of last issue: 2016-08-05 1.9 2016-10-27 100000010877 Date of first issue: 2015-07-14

Protective measures : The type of protective equipment must be selected based on

the Environmental Health and Safety risk assessment. Consult a Environmental Health and Safety expert if necessary.

Hygiene measures : Handle in accordance with good industrial hygiene and safety

practice.

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance : liquid

Colour : clear, light yellow

Odour : odourless

pH : 7,5

Solubility(ies)

Water solubility : soluble

10. STABILITY AND REACTIVITY

Reactivity : None reasonably foreseeable.

Chemical stability : Stable under recommended storage conditions.

Possibility of hazardous reac-

tions

No dangerous reaction known under conditions of normal use.

Conditions to avoid : To avoid thermal decomposition, do not overheat.

Incompatible materials : None known.

Hazardous decomposition

products

None known.

11. TOXICOLOGICAL INFORMATION

Acute toxicity

No data available

Skin corrosion/irritation

No data available

Serious eye damage/eye irritation

No data available

Respiratory or skin sensitisation

No data available

Germ cell mutagenicity

No data available

Nucleic acid dye



Version 1.9

Revision Date: 2016-10-27

SDS Number: 100000010877

Date of last issue: 2016-08-05 Date of first issue: 2015-07-14

Carcinogenicity

No data available

Reproductive toxicity

No data available

STOT - single exposure

No data available

STOT - repeated exposure

No data available

Repeated dose toxicity

No data available

Aspiration toxicity

No data available

12. ECOLOGICAL INFORMATION

Ecotoxicity

No data available

Persistence and degradability

No data available

Bioaccumulative potential

No data available

Mobility in soil

No data available

Other adverse effects

No data available

13. DISPOSAL CONSIDERATIONS

Disposal methods

Waste from residues : In accordance with National, Federal, State and Local regula-

tions.

Contaminated packaging : Empty containers should be taken to an approved waste han-

dling site for recycling or disposal.

14. TRANSPORT INFORMATION

International Regulations

UNRTDG

Not regulated as a dangerous good

IATA-DGR

Not regulated as a dangerous good

IMDG-Code

Nucleic acid dye



Version 1.9

Revision Date: 2016-10-27

SDS Number: 100000010877

Date of last issue: 2016-08-05 Date of first issue: 2015-07-14

Not regulated as a dangerous good

Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable for product as supplied.

15. REGULATORY INFORMATION

Safety, health and environmental regulations/legislation specific for the substance or mixture

Restricted to professional users.

16. OTHER INFORMATION

Full text of other abbreviations

AICS - Australian Inventory of Chemical Substances; ANTT - National Agency for Transport by Land of Brazil; ASTM - American Society for the Testing of Materials; bw - Body weight; CMR -Carcinogen, Mutagen or Reproductive Toxicant; CPR - Controlled Products Regulations; DIN -Standard of the German Institute for Standardisation; DSL - Domestic Substances List (Canada); ECx - Concentration associated with x% response; ELx - Loading rate associated with x% response; EmS - Emergency Schedule; ENCS - Existing and New Chemical Substances (Japan); ErCx - Concentration associated with x% growth rate response; ERG - Emergency Response Guide; GHS - Globally Harmonized System; GLP - Good Laboratory Practice; IARC - International Agency for Research on Cancer; IATA - International Air Transport Association; IBC - International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk; IC50 - Half maximal inhibitory concentration; ICAO - International Civil Aviation Organization; IECSC - Inventory of Existing Chemical Substances in China; IMDG - International Maritime Dangerous Goods; IMO - International Maritime Organization; ISHL - Industrial Safety and Health Law (Japan); ISO - International Organisation for Standardization; KECI - Korea Existing Chemicals Inventory; LC50 - Lethal Concentration to 50 % of a test population; LD50 - Lethal Dose to 50% of a test population (Median Lethal Dose); MARPOL - International Convention for the Prevention of Pollution from Ships; n.o.s. - Not Otherwise Specified; Nch - Chilean Norm; NO(A)EC -No Observed (Adverse) Effect Concentration; NO(A)EL - No Observed (Adverse) Effect Level; NOELR - No Observable Effect Loading Rate; NOM - Official Mexican Norm; NTP - National Toxicology Program; NZIoC - New Zealand Inventory of Chemicals; OECD - Organization for Economic Co-operation and Development; OPPTS - Office of Chemical Safety and Pollution Prevention; PBT - Persistent, Bioaccumulative and Toxic substance; PICCS - Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR - (Quantitative) Structure Activity Relationship; REACH - Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; SADT - Self-Accelerating Decomposition Temperature; SDS - Safety Data Sheet; TCSI - Taiwan Chemical Substance Inventory; TDG - Transportation of Dangerous Goods; TSCA - Toxic Substances Control Act (United States): UN - United Nations: UNRTDG - United Nations Recommendations on the Transport of Dangerous Goods; vPvB - Very Persistent and Very Bioaccumulative; WHMIS -Workplace Hazardous Materials Information System

Date format : dd.mm.yyyy

Numbers : 123 456,78

Nucleic acid dye



 Version
 Revision Date:
 SDS Number:
 Date of last issue: 2016-08-05

 1.9
 2016-10-27
 100000010877
 Date of first issue: 2015-07-14

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

IN / EN

Permeabilization reagent



Version Revision Date: SDS Number: Date of last issue: 2016-08-05 1.11 2016-10-27 100000010887 Date of first issue: 2015-07-14

1. PRODUCT AND COMPANY IDENTIFICATION

Product name : Permeabilization reagent

Substance name : Permeabilization reagent

7038

Manufacturer or supplier's details

Company : Janssen Diagnostics, LLC

Address : 700 US Highway Route 202

South Raritan, NJ 08869

US

Telephone : (877) 837-4339

Emergency telephone : CHEMTREC IN: 000-800-100-7141

number CHEMTREC International: +1 703-527-3887

E-mail address : SDSJanssen@its.jnj.com

Responsible/issuing person

Recommended use of the chemical and restrictions on use

Recommended use : Assay reagent

2. HAZARDS IDENTIFICATION

Manufacture, Storage and Import of Hazardous Chemicals Rules 1989

Classification

GHS Classification

Not a hazardous substance or mixture.

GHS label elements

Not a hazardous substance or mixture.

Other hazards which do not result in classification

None known.

3. COMPOSITION/INFORMATION ON INGREDIENTS

Substance / Mixture : Mixture

Chemical nature : Liquid

Hazardous components

Permeabilization reagent



Version **Revision Date:** SDS Number: Date of last issue: 2016-08-05 2016-10-27 100000010887 Date of first issue: 2015-07-14 1.11

Chemical name	CAS-No.	Concentration (% w/w)
sodium azide	26628-22-8	>= 0,1 - < 1

4. FIRST AID MEASURES

If inhaled If breathed in, move person into fresh air.

Consult a physician.

In case of skin contact Take off contaminated clothing and shoes immediately.

Wash off with plenty of water.

If symptoms persist, call a physician.

In case of eye contact Remove contact lenses.

If eye irritation persists, consult a specialist.

Rinse immediately with plenty of water, also under the eyelids,

for at least 5 minutes.

If swallowed If swallowed, rinse mouth with water (only if the person is con-

Call a physician immediately.

Most important symptoms and effects, both acute and

delayed

No information available.

Notes to physician Treat symptomatically.

5. FIREFIGHTING MEASURES

Suitable extinguishing media Use extinguishing measures that are appropriate to local cir-

cumstances and the surrounding environment.

Specific hazards during fire-

fighting

No information available.

Specific extinguishing meth-

No information available.

Special protective equipment:

for firefighters

In the event of fire, wear self-contained breathing apparatus.

6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protec- : tive equipment and emer-

gency procedures

In the event of an accidental release the emergency response team must respond based on a risk assessment and use per-

sonal protective equipment as appropriate.

Environmental precautions Should not be released into the environment.

Methods and materials for containment and cleaning up Large spills: Dam up. Soak up with inert absorbent material.

Keep in properly labelled containers.

Small spills: Gently cover the spill with an absorbent towel or

Permeabilization reagent



Version 1.11

Revision Date: 2016-10-27

SDS Number: 100000010887

Date of last issue: 2016-08-05 Date of first issue: 2015-07-14

pad.

Large spills + Small spills: Keep in suitable, closed containers for disposal. Treat recovered material as described in the sec-

tion "Disposal considerations".

7. HANDLING AND STORAGE

Advice on protection against

fire and explosion

No data available

Advice on safe handling : To

To avoid thermal decomposition, do not overheat.

Avoid inhalation, ingestion and contact with skin and eyes.

Use personal protective equipment as required.

Conditions for safe storage

To maintain product quality, do not store in heat or direct sun-

light.

Store in original container.

Keep containers tightly closed in a dry, cool and well-

ventilated place.

Keep away from heat and sources of ignition.

Keep locked up. Keep refrigerated.

Recommended storage tem: :

perature

2 - 8 °C

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Components with workplace control parameters

Components	CAS-No.	Value type (Form of exposure)	Control parameters / Permissible concentration	Basis
sodium azide	26628-22-8	C (Vapour)	0,11 ppm (Hydrazoic acid)	ACGIH
		С	0,29 mg/m3 (Sodium azide)	ACGIH

Engineering measures

All personal protective equipment should be based on a risk assessment. Consult a Environment Health Safety expert if

necessary.

Personal protective equipment

Respiratory protection

Engineering controls should always be the primary method of

controlling exposures.

If respiratory protective equipment is needed for certain activities, the type as well as the corresponding protection factor will depend upon the risk assessment and air concentrations, hazards, physical and warning properties of substances pre-

sent.

No personal respiratory protective equipment normally re-

quired.

Permeabilization reagent



Version Revision Date: SDS Number: Date of last issue: 2016-08-05 1.11 2016-10-27 100000010887 Date of first issue: 2015-07-14

Hand protection

Remarks : Disposable gloves

Eye protection : No special precautions required.

Skin and body protection : No special precautions required.

Protective measures : The type of protective equipment must be selected based on

the Environmental Health and Safety risk assessment. Consult a Environmental Health and Safety expert if necessary.

Hygiene measures : Handle in accordance with good industrial hygiene and safety

practice.

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance : liquid

Colour : clear

Odour : odourless

pH : 7,5

Solubility(ies)

Water solubility : soluble

10. STABILITY AND REACTIVITY

Reactivity : None reasonably foreseeable.

Chemical stability : Stable under recommended storage conditions.

Possibility of hazardous reac-

tions

No dangerous reaction known under conditions of normal use.

Conditions to avoid : To avoid thermal decomposition, do not overheat.

Incompatible materials : Strong acids and strong bases

Reducing agents
Strong oxidizing agents

Hazardous decomposition

products

None known.

11. TOXICOLOGICAL INFORMATION

Acute toxicity

Product:

Acute oral toxicity : Acute toxicity estimate: > 5 000 mg/kg

Method: Calculation method

Permeabilization reagent



Version 1.11

Revision Date: 2016-10-27

SDS Number: 100000010887

Date of last issue: 2016-08-05 Date of first issue: 2015-07-14

Components:

sodium azide:

Acute oral toxicity : LD50 (Rat): 27 mg/kg

Skin corrosion/irritation

No data available

Serious eye damage/eye irritation

No data available

Respiratory or skin sensitisation

No data available

Germ cell mutagenicity

No data available

Carcinogenicity

No data available

Reproductive toxicity

No data available

STOT - single exposure

No data available

STOT - repeated exposure

No data available

Repeated dose toxicity

No data available

Aspiration toxicity

No data available

12. ECOLOGICAL INFORMATION

Ecotoxicity

Components:

sodium azide:

Toxicity to fish : LC50 (Lepomis macrochirus (Bluegill sunfish)): 0,7 mg/l

Exposure time: 96 h

Toxicity to daphnia and other :

aquatic invertebrates

EC50 (Daphnia pulex (Water flea)): 4,2 mg/l

Exposure time: 96 h

Toxicity to algae : IC50: 272 mg/l

Toxicity to bacteria : EC50 (Photobacterium phosphoreum): 38,5 mg/l

Permeabilization reagent



 Version
 Revision Date:
 SDS Number:
 Date of last issue: 2016-08-05

 1.11
 2016-10-27
 100000010887
 Date of first issue: 2015-07-14

Persistence and degradability

No data available

Bioaccumulative potential

No data available

Mobility in soil

No data available

Other adverse effects

No data available

13. DISPOSAL CONSIDERATIONS

Disposal methods

Waste from residues : In accordance with National, Federal, State and Local regula-

tions.

Contaminated packaging : Empty containers should be taken to an approved waste han-

dling site for recycling or disposal.

14. TRANSPORT INFORMATION

International Regulations

UNRTDG

Not regulated as a dangerous good

IATA-DGR

Not regulated as a dangerous good

IMDG-Code

Not regulated as a dangerous good

Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable for product as supplied.

15. REGULATORY INFORMATION

Safety, health and environmental regulations/legislation specific for the substance or mix-

Restricted to professional users.

16. OTHER INFORMATION

Full text of other abbreviations

AICS - Australian Inventory of Chemical Substances; ANTT - National Agency for Transport by Land of Brazil; ASTM - American Society for the Testing of Materials; bw - Body weight; CMR - Carcinogen, Mutagen or Reproductive Toxicant; CPR - Controlled Products Regulations; DIN - Standard of the German Institute for Standardisation; DSL - Domestic Substances List (Canada);

Permeabilization reagent



Version Revision Date: SDS Number: Date of last issue: 2016-08-05 1.11 2016-10-27 100000010887 Date of first issue: 2015-07-14

ECx - Concentration associated with x% response; ELx - Loading rate associated with x% response; EmS - Emergency Schedule; ENCS - Existing and New Chemical Substances (Japan); ErCx - Concentration associated with x% growth rate response; ERG - Emergency Response Guide; GHS - Globally Harmonized System; GLP - Good Laboratory Practice; IARC - International Agency for Research on Cancer; IATA - International Air Transport Association; IBC - International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk; IC50 - Half maximal inhibitory concentration; ICAO - International Civil Aviation Organization; IECSC - Inventory of Existing Chemical Substances in China; IMDG - International Maritime Dangerous Goods; IMO - International Maritime Organization; ISHL - Industrial Safety and Health Law (Japan); ISO - International Organisation for Standardization; KECI - Korea Existing Chemicals Inventory; LC50 - Lethal Concentration to 50 % of a test population; LD50 - Lethal Dose to 50% of a test population (Median Lethal Dose); MARPOL - International Convention for the Prevention of Pollution from Ships; n.o.s. - Not Otherwise Specified; Nch - Chilean Norm; NO(A)EC -No Observed (Adverse) Effect Concentration; NO(A)EL - No Observed (Adverse) Effect Level; NOELR - No Observable Effect Loading Rate; NOM - Official Mexican Norm; NTP - National Toxicology Program; NZIoC - New Zealand Inventory of Chemicals; OECD - Organization for Economic Co-operation and Development; OPPTS - Office of Chemical Safety and Pollution Prevention; PBT - Persistent, Bioaccumulative and Toxic substance; PICCS - Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR - (Quantitative) Structure Activity Relationship; REACH - Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; SADT - Self-Accelerating Decomposition Temperature; SDS - Safety Data Sheet; TCSI - Taiwan Chemical Substance Inventory; TDG - Transportation of Dangerous Goods; TSCA - Toxic Substances Control Act (United States); UN - United Nations; UNRTDG - United Nations Recommendations on the Transport of Dangerous Goods; vPvB - Very Persistent and Very Bioaccumulative; WHMIS -Workplace Hazardous Materials Information System

Date format : dd.mm.yyyy

Numbers : 123 456.78

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

IN / EN

Anti-EpCAM ferrofluid



Version Revision Date: SDS Number: Date of last issue: 2016-10-27 1.11 2016-11-15 100000010880 Date of first issue: 2015-07-14

1. PRODUCT AND COMPANY IDENTIFICATION

Product name : Anti-EpCAM ferrofluid

Substance name : Anti-EpCAM ferrofluid

7036

Manufacturer or supplier's details

Company : Janssen Diagnostics, LLC

Address : 700 US Highway Route 202

South Raritan, NJ 08869

US

Telephone : (877) 837-4339

Emergency telephone : CHEMTREC IN: 000-800-100-7141

number CHEMTREC International: +1 703-527-3887

E-mail address : SDSJanssen@its.jnj.com

Responsible/issuing person

Recommended use of the chemical and restrictions on use

Recommended use : Assay reagent

2. HAZARDS IDENTIFICATION

Manufacture, Storage and Import of Hazardous Chemicals Rules 1989

Classification

GHS Classification

Not a hazardous substance or mixture.

GHS label elements

Not a hazardous substance or mixture.

Other hazards which do not result in classification

None known.

3. COMPOSITION/INFORMATION ON INGREDIENTS

Substance / Mixture : Mixture

Chemical nature : Liquid

Hazardous components

Anti-EpCAM ferrofluid



 Version
 Revision Date:
 SDS Number:
 Date of last issue: 2016-10-27

 1.11
 2016-11-15
 100000010880
 Date of first issue: 2015-07-14

Chemical name	CAS-No.	Concentration (%
		w/w)
Anti-EpCAM mouse mAb conjugated to Ferrofluid	Not Assigned	< 0,1

4. FIRST AID MEASURES

If inhaled : If breathed in, move person into fresh air.

Consult a physician.

In case of skin contact : Take off contaminated clothing and shoes immediately.

Wash off immediately with plenty of water. If symptoms persist, call a physician. Wash contaminated clothing before re-use.

In case of eye contact : Rinse immediately with plenty of water, also under the eyelids,

for at least 5 minutes. Remove contact lenses.

If eye irritation persists, consult a specialist.

If swallowed : If swallowed, rinse mouth with water (only if the person is con-

scious).

Call a physician immediately.

Most important symptoms and effects, both acute and

delayed

No information available.

Notes to physician : Treat symptomatically.

5. FIREFIGHTING MEASURES

Suitable extinguishing media : Use extinguishing measures that are appropriate to local cir-

cumstances and the surrounding environment.

Specific hazards during fire-

fighting

No information available.

Specific extinguishing meth-

ods

No information available.

Special protective equipment:

for firefighters

In the event of fire, wear self-contained breathing apparatus.

6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment and emer-

gency procedures

In the event of an accidental release the emergency response team must respond based on a risk assessment and use per-

sonal protective equipment as appropriate.

Environmental precautions : Should not be released into the environment.

Methods and materials for containment and cleaning up

Large spills: Dam up. Soak up with inert absorbent material.

Keep in properly labelled containers.

Anti-EpCAM ferrofluid



Version 1.11

Revision Date: 2016-11-15

SDS Number: 100000010880

Date of last issue: 2016-10-27 Date of first issue: 2015-07-14

Small spills: Gently cover the spill with an absorbent towel or

pad.

Large spills + Small spills: Keep in suitable, closed containers for disposal. Treat recovered material as described in the sec-

tion "Disposal considerations".

7. HANDLING AND STORAGE

Advice on protection against

fire and explosion

No data available

Advice on safe handling : To avoid thermal decomposition, do not overheat.

Avoid inhalation, ingestion and contact with skin and eyes.

Use personal protective equipment as required.

Conditions for safe storage

To maintain product quality, do not store in heat or direct sun-

light.

Store in original container.

Keep containers tightly closed in a dry, cool and well-

ventilated place.

Keep away from heat and sources of ignition.

Keep locked up. Keep refrigerated.

Recommended storage tem- :

perature

2 - 8 °C

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Components with workplace control parameters

Components	CAS-No.	Value type (Form of exposure)	Control parameters / Permissible concentration	Basis
Anti-EpCAM mouse mAb conjugated to Ferrofluid	Not Assigned	PBOEL-HHC	2	J&J OEL/PBOEL HHC
	Further information: J&J has a hazard banding notation: PBOEL HHC. This substance is classified by J&J as being PBOEL HHC 2. This means that the OEL is estimated to be from 20 to 100 μg/m3			

Engineering measures

All personal protective equipment should be based on a risk assessment. Consult a Environment Health Safety expert if necessary.

Personal protective equipment

Respiratory protection

Engineering controls should always be the primary method of

controlling exposures.

If respiratory protective equipment is needed for certain activities, the type as well as the corresponding protection factor will depend upon the risk assessment and air concentrations, hazards, physical and warning properties of substances pre-

Anti-EpCAM ferrofluid



Version Revision Date: SDS Number: Date of last issue: 2016-10-27 1.11 2016-11-15 100000010880 Date of first issue: 2015-07-14

sent.

No personal respiratory protective equipment normally re-

quired.

Hand protection

Remarks : Disposable gloves

Eye protection : No special precautions required.

Skin and body protection : No special precautions required.

Protective measures : The type of protective equipment must be selected based on

the Environmental Health and Safety risk assessment. Consult a Environmental Health and Safety expert if necessary.

Hygiene measures : Handle in accordance with good industrial hygiene and safety

practice.

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance : liquid

Colour : brown

pH : 7,5

Solubility(ies)

Water solubility : soluble

10. STABILITY AND REACTIVITY

Reactivity : None reasonably foreseeable.

Chemical stability : Stable under recommended storage conditions.

Possibility of hazardous reac-

tions

No dangerous reaction known under conditions of normal use.

Conditions to avoid : To avoid thermal decomposition, do not overheat.

Incompatible materials : None known.

Hazardous decomposition

products

None known.

11. TOXICOLOGICAL INFORMATION

Acute toxicity

No data available

Skin corrosion/irritation

No data available

Anti-EpCAM ferrofluid



Version 1.11

Revision Date: 2016-11-15

SDS Number: 100000010880

Date of last issue: 2016-10-27 Date of first issue: 2015-07-14

Serious eye damage/eye irritation

No data available

Respiratory or skin sensitisation

No data available

Germ cell mutagenicity

No data available

Carcinogenicity

No data available

Reproductive toxicity

No data available

STOT - single exposure

No data available

STOT - repeated exposure

No data available

Repeated dose toxicity

No data available

Aspiration toxicity

No data available

12. ECOLOGICAL INFORMATION

Ecotoxicity

No data available

Persistence and degradability

No data available

Bioaccumulative potential

No data available

Mobility in soil

No data available

Other adverse effects

No data available

13. DISPOSAL CONSIDERATIONS

Disposal methods

Waste from residues : In accordance with National, Federal, State and Local regula-

tions.

Contaminated packaging : Empty containers should be taken to an approved waste han-

dling site for recycling or disposal.

Anti-EpCAM ferrofluid



 Version
 Revision Date:
 SDS Number:
 Date of last issue: 2016-10-27

 1.11
 2016-11-15
 100000010880
 Date of first issue: 2015-07-14

14. TRANSPORT INFORMATION

International Regulations

UNRTDG

Not regulated as a dangerous good

IATA-DGR

Not regulated as a dangerous good

IMDG-Code

Not regulated as a dangerous good

Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable for product as supplied.

15. REGULATORY INFORMATION

Safety, health and environmental regulations/legislation specific for the substance or mixture

Restricted to professional users.

16. OTHER INFORMATION

Full text of other abbreviations

AICS - Australian Inventory of Chemical Substances; ANTT - National Agency for Transport by Land of Brazil; ASTM - American Society for the Testing of Materials; bw - Body weight; CMR -Carcinogen, Mutagen or Reproductive Toxicant; CPR - Controlled Products Regulations; DIN -Standard of the German Institute for Standardisation; DSL - Domestic Substances List (Canada); ECx - Concentration associated with x% response; ELx - Loading rate associated with x% response; EmS - Emergency Schedule; ENCS - Existing and New Chemical Substances (Japan); ErCx - Concentration associated with x% growth rate response; ERG - Emergency Response Guide; GHS - Globally Harmonized System; GLP - Good Laboratory Practice; IARC - International Agency for Research on Cancer; IATA - International Air Transport Association; IBC - International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk; IC50 - Half maximal inhibitory concentration; ICAO - International Civil Aviation Organization; IECSC - Inventory of Existing Chemical Substances in China; IMDG - International Maritime Dangerous Goods; IMO - International Maritime Organization; ISHL - Industrial Safety and Health Law (Japan); ISO - International Organisation for Standardization; KECI - Korea Existing Chemicals Inventory; LC50 - Lethal Concentration to 50 % of a test population; LD50 - Lethal Dose to 50% of a test population (Median Lethal Dose); MARPOL - International Convention for the Prevention of Pollution from Ships; n.o.s. - Not Otherwise Specified; Nch - Chilean Norm; NO(A)EC -No Observed (Adverse) Effect Concentration; NO(A)EL - No Observed (Adverse) Effect Level; NOELR - No Observable Effect Loading Rate: NOM - Official Mexican Norm: NTP - National Toxicology Program; NZIoC - New Zealand Inventory of Chemicals; OECD - Organization for Economic Co-operation and Development; OPPTS - Office of Chemical Safety and Pollution Prevention; PBT - Persistent, Bioaccumulative and Toxic substance; PICCS - Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR - (Quantitative) Structure Activity Relationship; REACH - Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; SADT - Self-Accelerating Decomposition Temperature; SDS - Safety Data Sheet; TCSI - Taiwan Chemical

Anti-EpCAM ferrofluid



 Version
 Revision Date:
 SDS Number:
 Date of last issue: 2016-10-27

 1.11
 2016-11-15
 100000010880
 Date of first issue: 2015-07-14

Substance Inventory; TDG - Transportation of Dangerous Goods; TSCA - Toxic Substances Control Act (United States); UN - United Nations; UNRTDG - United Nations Recommendations on the Transport of Dangerous Goods; vPvB - Very Persistent and Very Bioaccumulative; WHMIS - Workplace Hazardous Materials Information System

Date format : dd.mm.yyyy

Numbers : 123 456,78

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

IN / EN

Dilution buffer



Version Revision Date: SDS Number: Date of last issue: 2016-08-05 1.11 2016-10-27 100000010879 Date of first issue: 2015-07-14

1. PRODUCT AND COMPANY IDENTIFICATION

Product name : Dilution buffer

Substance name : Dilution buffer

7039

Manufacturer or supplier's details

Company : Janssen Diagnostics, LLC

Address : 700 US Highway Route 202

South Raritan, NJ 08869

US

Telephone : (877) 837-4339

Emergency telephone : CHEMTREC IN: 000-800-100-7141

number CHEMTREC International: +1 703-527-3887

E-mail address : SDSJanssen@its.jnj.com

Responsible/issuing person

Recommended use of the chemical and restrictions on use

Recommended use : Assay reagent

2. HAZARDS IDENTIFICATION

Manufacture, Storage and Import of Hazardous Chemicals Rules 1989

Classification

GHS Classification

Not a hazardous substance or mixture.

GHS label elements

Not a hazardous substance or mixture.

Other hazards which do not result in classification

None known.

3. COMPOSITION/INFORMATION ON INGREDIENTS

Substance / Mixture : Mixture

Chemical nature : Liquid

Hazardous components

Dilution buffer



 Version
 Revision Date:
 SDS Number:
 Date of last issue: 2016-08-05

 1.11
 2016-10-27
 100000010879
 Date of first issue: 2015-07-14

Chemical name	CAS-No.	Concentration (%
		w/w)
sodium azide	26628-22-8	>= 0,1 - < 1

4. FIRST AID MEASURES

If inhaled : If breathed in, move person into fresh air.

Consult a physician.

In case of skin contact : Take off contaminated clothing and shoes immediately.

Wash off with plenty of water.

If symptoms persist, call a physician.

In case of eye contact : Rinse immediately with plenty of water, also under the eyelids,

for at least 5 minutes. Remove contact lenses.

If eye irritation persists, consult a specialist.

If swallowed : If swallowed, rinse mouth with water (only if the person is con-

scious)

Call a physician immediately.

Most important symptoms and effects, both acute and

delayed

No information available.

Notes to physician : Treat symptomatically.

5. FIREFIGHTING MEASURES

Suitable extinguishing media : Use extinguishing measures that are appropriate to local cir-

cumstances and the surrounding environment.

Specific hazards during fire-

fighting

No information available.

Specific extinguishing meth-

ods

No information available.

Special protective equipment:

for firefighters

In the event of fire, wear self-contained breathing apparatus.

6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment and emer-

gency procedures

: In the event of an accidental release the emergency response team must respond based on a risk assessment and use per-

sonal protective equipment as appropriate.

Environmental precautions : Should not be released into the environment.

Methods and materials for containment and cleaning up

: Large spills: Dam up. Soak up with inert absorbent material.

Keep in properly labelled containers.

Small spills: Gently cover the spill with an absorbent towel or

Dilution buffer



Version 1.11

Revision Date: 2016-10-27

SDS Number: 100000010879

Date of last issue: 2016-08-05 Date of first issue: 2015-07-14

pad.

Large spills + Small spills: Keep in suitable, closed containers for disposal. Treat recovered material as described in the sec-

tion "Disposal considerations".

7. HANDLING AND STORAGE

Advice on protection against

fire and explosion

No data available

Advice on safe handling

To avoid thermal decomposition, do not overheat.

Avoid inhalation, ingestion and contact with skin and eyes.

Use personal protective equipment as required.

Conditions for safe storage

To maintain product quality, do not store in heat or direct sun-

light.

Store in original container.

Keep containers tightly closed in a dry, cool and well-

ventilated place.

Keep away from heat and sources of ignition.

Keep locked up. Keep refrigerated.

Recommended storage tem: :

perature

2 - 8 °C

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Components with workplace control parameters

Components	CAS-No.	Value type (Form of exposure)	Control parameters / Permissible concentration	Basis
sodium azide	26628-22-8	C (Vapour)	0,11 ppm (Hydrazoic acid)	ACGIH
		С	0,29 mg/m3 (Sodium azide)	ACGIH

Engineering measures

All personal protective equipment should be based on a risk assessment. Consult a Environment Health Safety expert if necessary.

Personal protective equipment

Respiratory protection

Engineering controls should always be the primary method of

controlling exposures.

If respiratory protective equipment is needed for certain activities, the type as well as the corresponding protection factor will depend upon the risk assessment and air concentrations, hazards, physical and warning properties of substances pre-

sent.

No personal respiratory protective equipment normally re-

quired.

Dilution buffer



Version Revision Date: SDS Number: Date of last issue: 2016-08-05 1.11 2016-10-27 100000010879 Date of first issue: 2015-07-14

Hand protection

Remarks : Disposable gloves

Eye protection : No special precautions required.

Skin and body protection : No special precautions required.

Protective measures : The type of protective equipment must be selected based on

the Environmental Health and Safety risk assessment. Consult a Environmental Health and Safety expert if necessary.

Hygiene measures : Handle in accordance with good industrial hygiene and safety

practice.

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance : liquid

Colour : clear

Odour : odourless

pH : 7,5

Solubility(ies)

Water solubility : soluble

10. STABILITY AND REACTIVITY

Reactivity : None reasonably foreseeable.

Chemical stability : Stable under recommended storage conditions.

Possibility of hazardous reac-

tions

No dangerous reaction known under conditions of normal use.

Conditions to avoid : To avoid thermal decomposition, do not overheat.

Incompatible materials : Oxidizing agents

Hazardous decomposition

products

None known.

11. TOXICOLOGICAL INFORMATION

Acute toxicity

Product:

Acute oral toxicity : Acute toxicity estimate: > 5 000 mg/kg

Method: Calculation method

Dilution buffer



Version Revision Date: SDS Number: Date of last issue: 2016-08-05 1.11 2016-10-27 100000010879 Date of first issue: 2015-07-14

Components:

sodium azide:

Acute oral toxicity LD50 (Rat): 27 mg/kg

Skin corrosion/irritation

No data available

Serious eye damage/eye irritation

No data available

Respiratory or skin sensitisation

No data available

Germ cell mutagenicity

No data available

Carcinogenicity

No data available

Reproductive toxicity

No data available

STOT - single exposure

No data available

STOT - repeated exposure

No data available

Repeated dose toxicity

No data available

Aspiration toxicity

No data available

12. ECOLOGICAL INFORMATION

Ecotoxicity

Components:

sodium azide:

Toxicity to fish LC50 (Lepomis macrochirus (Bluegill sunfish)): 0,7 mg/l

Exposure time: 96 h

aquatic invertebrates

Toxicity to daphnia and other : EC50 (Daphnia pulex (Water flea)): 4,2 mg/l

Exposure time: 96 h

Toxicity to algae : IC50: 272 mg/l

Toxicity to bacteria EC50 (Photobacterium phosphoreum): 38,5 mg/l

Dilution buffer



Version 1.11 Revision Date: 2016-10-27

SDS Number: 100000010879

Date of last issue: 2016-08-05 Date of first issue: 2015-07-14

Persistence and degradability

No data available

Bioaccumulative potential

No data available

Mobility in soil

No data available

Other adverse effects

No data available

13. DISPOSAL CONSIDERATIONS

Disposal methods

Waste from residues

In accordance with National, Federal, State and Local regula-

tions.

Contaminated packaging

Empty containers should be taken to an approved waste han-

dling site for recycling or disposal.

14. TRANSPORT INFORMATION

International Regulations

UNRTDG

Not regulated as a dangerous good

IATA-DGR

Not regulated as a dangerous good

IMDG-Code

Not regulated as a dangerous good

Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable for product as supplied.

15. REGULATORY INFORMATION

Safety, health and environmental regulations/legislation specific for the substance or mixture

Restricted to professional users.

16. OTHER INFORMATION

Full text of other abbreviations

AICS - Australian Inventory of Chemical Substances; ANTT - National Agency for Transport by Land of Brazil; ASTM - American Society for the Testing of Materials; bw - Body weight; CMR - Carcinogen, Mutagen or Reproductive Toxicant; CPR - Controlled Products Regulations; DIN - Standard of the German Institute for Standardisation; DSL - Domestic Substances List (Canada);

Dilution buffer



Version Revision Date: SDS Number: Date of last issue: 2016-08-05 1.11 2016-10-27 100000010879 Date of first issue: 2015-07-14

ECx - Concentration associated with x% response; ELx - Loading rate associated with x% response; EmS - Emergency Schedule; ENCS - Existing and New Chemical Substances (Japan); ErCx - Concentration associated with x% growth rate response; ERG - Emergency Response Guide; GHS - Globally Harmonized System; GLP - Good Laboratory Practice; IARC - International Agency for Research on Cancer; IATA - International Air Transport Association; IBC - International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk; IC50 - Half maximal inhibitory concentration; ICAO - International Civil Aviation Organization; IECSC - Inventory of Existing Chemical Substances in China; IMDG - International Maritime Dangerous Goods; IMO - International Maritime Organization; ISHL - Industrial Safety and Health Law (Japan); ISO - International Organisation for Standardization; KECI - Korea Existing Chemicals Inventory; LC50 - Lethal Concentration to 50 % of a test population; LD50 - Lethal Dose to 50% of a test population (Median Lethal Dose); MARPOL - International Convention for the Prevention of Pollution from Ships; n.o.s. - Not Otherwise Specified; Nch - Chilean Norm; NO(A)EC -No Observed (Adverse) Effect Concentration; NO(A)EL - No Observed (Adverse) Effect Level; NOELR - No Observable Effect Loading Rate; NOM - Official Mexican Norm; NTP - National Toxicology Program; NZIoC - New Zealand Inventory of Chemicals; OECD - Organization for Economic Co-operation and Development; OPPTS - Office of Chemical Safety and Pollution Prevention; PBT - Persistent, Bioaccumulative and Toxic substance; PICCS - Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR - (Quantitative) Structure Activity Relationship; REACH - Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; SADT - Self-Accelerating Decomposition Temperature; SDS - Safety Data Sheet; TCSI - Taiwan Chemical Substance Inventory; TDG - Transportation of Dangerous Goods; TSCA - Toxic Substances Control Act (United States); UN - United Nations; UNRTDG - United Nations Recommendations on the Transport of Dangerous Goods; vPvB - Very Persistent and Very Bioaccumulative; WHMIS -Workplace Hazardous Materials Information System

Date format : dd.mm.yyyy

Numbers : 123 456.78

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

IN / EN

Cell fixative



Version Revision Date: SDS Number: Date of last issue: 2016-10-27 1.13 2016-11-15 100000010702 Date of first issue: 2015-07-14

1. PRODUCT AND COMPANY IDENTIFICATION

Product name : Cell fixative

Substance name : Cell fixative

7042

Manufacturer or supplier's details

Company : Janssen Diagnostics, LLC

Address : 700 US Highway Route 202

South Raritan, NJ 08869

US

Telephone : (877) 837-4339

Emergency telephone : CHEMTREC IN: 000-800-100-7141

number CHEMTREC International: +1 703-527-3887

E-mail address : SDSJanssen@its.jnj.com

Responsible/issuing person

Recommended use of the chemical and restrictions on use

Recommended use : Assay reagent

2. HAZARDS IDENTIFICATION

Manufacture, Storage and Import of Hazardous Chemicals Rules 1989

Classification

GHS Classification

Skin sensitisation : Category 1

GHS label elements

Hazard pictograms :

Signal word : Warning

Hazard statements : H317 May cause an allergic skin reaction.

Precautionary statements : Prevention:

P261 Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.

Cell fixative



Version **Revision Date:** SDS Number: Date of last issue: 2016-10-27 2016-11-15 100000010702 Date of first issue: 2015-07-14 1.13

P272 Contaminated work clothing should not be allowed out of

the workplace.

P280 Wear protective gloves.

Response:

P302 + P352 IF ON SKIN: Wash with plenty of water.

P333 + P313 If skin irritation or rash occurs: Get medical ad-

vice/ attention.

P362 + P364 Take off contaminated clothing and wash it before

reuse.

Disposal:

P501 Dispose of contents/ container to an approved waste

disposal plant.

Other hazards which do not result in classification

None known.

3. COMPOSITION/INFORMATION ON INGREDIENTS

Substance / Mixture Mixture

Chemical nature Liquid

Hazardous components

Chemical name	CAS-No.	Concentration (%
		w/w)
IMIDUREA	39236-46-9	>= 1 - < 10
Sodium chloride (NaCl)	7647-14-5	>= 1 - < 10
sodium azide	26628-22-8	>= 0,1 - < 1

4. FIRST AID MEASURES

If inhaled If breathed in, move person into fresh air.

Consult a physician.

In case of skin contact Take off contaminated clothing and shoes immediately.

Wash off with plenty of water.

If symptoms persist, call a physician.

Rinse immediately with plenty of water, also under the eyelids, In case of eye contact

> for at least 15 minutes. Remove contact lenses.

If eye irritation persists, consult a specialist.

If swallowed If swallowed, rinse mouth with water (only if the person is con-

scious).

Call a physician immediately.

Most important symptoms and effects, both acute and

delayed

No information available.

Cell fixative



Version **Revision Date:** Date of last issue: 2016-10-27 SDS Number: 2016-11-15 100000010702 Date of first issue: 2015-07-14 1.13

Notes to physician Treat symptomatically.

5. FIREFIGHTING MEASURES

Suitable extinguishing media Use extinguishing measures that are appropriate to local cir-

cumstances and the surrounding environment.

Specific hazards during fire-

fighting

No information available.

Hazardous combustion prod: :

ucts

No hazardous combustion products are known

Specific extinguishing meth-

ods

No information available.

Special protective equipment:

for firefighters

In the event of fire, wear self-contained breathing apparatus.

6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protec- : tive equipment and emer-

gency procedures

In the event of an accidental release the emergency response team must respond based on a risk assessment and use per-

sonal protective equipment as appropriate.

Should not be released into the environment. Environmental precautions

Methods and materials for containment and cleaning up

Large spills: Dam up. Soak up with inert absorbent material.

Keep in properly labelled containers.

Small spills: Gently cover the spill with an absorbent towel or

pad.

Large spills + Small spills: Keep in suitable, closed containers for disposal. Treat recovered material as described in the sec-

tion "Disposal considerations".

7. HANDLING AND STORAGE

Advice on protection against

fire and explosion

No data available

To avoid thermal decomposition, do not overheat. Advice on safe handling

Avoid inhalation, ingestion and contact with skin and eyes.

Use personal protective equipment as required.

Conditions for safe storage

To maintain product quality, do not store in heat or direct sun-

Store in original container.

Keep containers tightly closed in a dry, cool and well-

ventilated place.

Keep away from heat and sources of ignition.

Keep locked up. Keep refrigerated.

Cell fixative



Version **Revision Date:** SDS Number: Date of last issue: 2016-10-27 2016-11-15 100000010702 Date of first issue: 2015-07-14 1.13

Recommended storage tem- : 2 - 8 °C

perature

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Components with workplace control parameters

Components	CAS-No.	Value type	Control parame-	Basis
		(Form of	ters / Permissible	
		exposure)	concentration	
sodium azide	26628-22-8	C (Vapour)	0,11 ppm	ACGIH
			(Hydrazoic acid)	
		С	0,29 mg/m3	ACGIH
			(Sodium azide)	

Engineering measures All personal protective equipment should be based on a risk

assessment. Consult a Environment Health Safety expert if

necessary.

Personal protective equipment

Respiratory protection Engineering controls should always be the primary method of

controlling exposures.

If respiratory protective equipment is needed for certain activities, the type as well as the corresponding protection factor will depend upon the risk assessment and air concentrations, hazards, physical and warning properties of substances pre-

No personal respiratory protective equipment normally re-

quired.

Hand protection

Remarks Disposable gloves

Eye protection No special precautions required.

Skin and body protection No special precautions required.

Protective measures The type of protective equipment must be selected based on

> the Environmental Health and Safety risk assessment. Consult a Environmental Health and Safety expert if necessary.

Hygiene measures Handle in accordance with good industrial hygiene and safety

practice.

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance liquid

Colour clear

Odour odourless

Cell fixative



 Version
 Revision Date:
 SDS Number:
 Date of last issue: 2016-10-27

 1.13
 2016-11-15
 100000010702
 Date of first issue: 2015-07-14

pH : 7,5

Solubility(ies)

Water solubility : soluble

10. STABILITY AND REACTIVITY

Reactivity : None reasonably foreseeable.

Chemical stability : Stable under recommended storage conditions.

Possibility of hazardous reac-

tions

No dangerous reaction known under conditions of normal use.

Conditions to avoid : To avoid thermal decomposition, do not overheat.

Incompatible materials : Strong acids and strong bases

Reducing agents Oxidizing agents

Hazardous decomposition

products

None known.

11. TOXICOLOGICAL INFORMATION

Acute toxicity

Product:

Acute oral toxicity : Acute toxicity estimate: > 5 000 mg/kg

Method: Calculation method

Components:

IMIDUREA:

Acute oral toxicity : LD50 (Rat): 11 300 mg/kg

Acute inhalation toxicity : LC50 (Rat): > 5,5 mg/l

Exposure time: 1 h

Acute dermal toxicity : LD50 (Rabbit): > 5 000 mg/kg

Sodium chloride (NaCI):

Acute oral toxicity : LD50 Oral (Rat): 3 000 mg/kg

Assessment: The component/mixture is low toxic after single

ingestion.

Acute inhalation toxicity : Remarks: No data available

Acute dermal toxicity : Remarks: No data available

sodium azide:

Cell fixative



Version Revision Date: SDS Number: Date of last issue: 2016-10-27 1.13 2016-11-15 100000010702 Date of first issue: 2015-07-14

Acute oral toxicity : LD50 (Rat): 27 mg/kg

Skin corrosion/irritation

Components:

IMIDUREA:

Result: No skin irritation

Sodium chloride (NaCI): Remarks: No data available

Serious eye damage/eye irritation

Components:

IMIDUREA:

Result: No eye irritation

Sodium chloride (NaCl): Remarks: No data available

Respiratory or skin sensitisation

Components:

IMIDUREA:

Method: Maximisation Test

Result: May cause sensitisation by skin contact.

Method: Local Lymph Node Assay (LLNA) in mice Result: May cause sensitisation by skin contact.

Sodium chloride (NaCI):

Remarks: No data available

Germ cell mutagenicity

Components:

IMIDUREA:

Genotoxicity in vitro : Test Type: Ames test

Metabolic activation: with and without metabolic activation

Result: negative

GLP: yes

: Test Type: Chromosome aberration test in vitro

Metabolic activation: with and without metabolic activation

Result: negative

Genotoxicity in vivo : Test Type: Micronucleus test

Cell fixative



Version 1.13

Revision Date: 2016-11-15

SDS Number: 100000010702

Date of last issue: 2016-10-27 Date of first issue: 2015-07-14

Species: Mouse Application Route: Oral Remarks: negative

Germ cell mutagenicity -

Assessment

No information available.

Sodium chloride (NaCI):

Germ cell mutagenicity -

Assessment

No information available.

Carcinogenicity

Components:

IMIDUREA:

Carcinogenicity - Assess-

ment

No information available.

Sodium chloride (NaCI):

Carcinogenicity - Assess-

ment

No information available.

Reproductive toxicity

Components:

IMIDUREA:

Teratogenicity - Assessment : No information available.

Sodium chloride (NaCI):

Reproductive toxicity - As-

sessment

: No information available.

Teratogenicity - Assessment : No information available.

STOT - single exposure

Components:

Sodium chloride (NaCI):

Remarks: No data available

STOT - repeated exposure

No data available

Cell fixative



 Version
 Revision Date:
 SDS Number:
 Date of last issue: 2016-10-27

 1.13
 2016-11-15
 100000010702
 Date of first issue: 2015-07-14

Repeated dose toxicity

Components:

IMIDUREA:

Species: Rat NOAEL: 200 mg/kg LOAEL: 500 mg/kg Application Route: Oral

Species: Rabbit NOAEL: 200 mg/kg Application Route: Dermal

Aspiration toxicity

No data available

12. ECOLOGICAL INFORMATION

Ecotoxicity

Components:

IMIDUREA:

Toxicity to fish : Remarks: No data available

Sodium chloride (NaCI):

Toxicity to fish : LC50 (Fish): 6 750 mg/l

Exposure time: 96 h

EC50 (Daphnia (water flea)): 2 024 mg/l

Exposure time: 48 h

sodium azide:

Toxicity to fish : LC50 (Lepomis macrochirus (Bluegill sunfish)): 0,7 mg/l

Exposure time: 96 h

Toxicity to daphnia and other :

aquatic invertebrates

EC50 (Daphnia pulex (Water flea)): 4,2 mg/l

Exposure time: 96 h

Toxicity to algae : IC50: 272 mg/l

Toxicity to bacteria : EC50 (Photobacterium phosphoreum): 38,5 mg/l

Persistence and degradability

Components:

IMIDUREA:

Biodegradability : Remarks: No data available

Cell fixative



Version 1.13

Revision Date: 2016-11-15

SDS Number: 100000010702

Date of last issue: 2016-10-27 Date of first issue: 2015-07-14

Sodium chloride (NaCl):

Biodegradability : Remarks: No data available

Bioaccumulative potential

Components:

IMIDUREA:

Bioaccumulation : Remarks: No data available

Sodium chloride (NaCI):

Bioaccumulation : Remarks: No data available

Mobility in soil

Components:

IMIDUREA:

Distribution among environ-

mental compartments

Remarks: No data available

Sodium chloride (NaCI):

Mobility : Remarks: No data available

Other adverse effects

Components:

IMIDUREA:

Results of PBT and vPvB

assessment

No information available.

Additional ecological infor-

mation

No data available

Sodium chloride (NaCl):

Environmental fate and

pathways

: No data available

Results of PBT and vPvB

assessment

No information available.

Additional ecological infor-

mation

No data available

Cell fixative



 Version
 Revision Date:
 SDS Number:
 Date of last issue: 2016-10-27

 1.13
 2016-11-15
 100000010702
 Date of first issue: 2015-07-14

13. DISPOSAL CONSIDERATIONS

Disposal methods

Waste from residues : In accordance with National, Federal, State and Local regula-

tions.

Contaminated packaging : Empty containers should be taken to an approved waste han-

dling site for recycling or disposal.

14. TRANSPORT INFORMATION

International Regulations

UNRTDG

Not regulated as a dangerous good

IATA-DGR

Not regulated as a dangerous good

IMDG-Code

Not regulated as a dangerous good

Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable for product as supplied.

15. REGULATORY INFORMATION

Safety, health and environmental regulations/legislation specific for the substance or mixture

Restricted to professional users.

16. OTHER INFORMATION

Full text of other abbreviations

AICS - Australian Inventory of Chemical Substances; ANTT - National Agency for Transport by Land of Brazil; ASTM - American Society for the Testing of Materials; bw - Body weight; CMR - Carcinogen, Mutagen or Reproductive Toxicant; CPR - Controlled Products Regulations; DIN - Standard of the German Institute for Standardisation; DSL - Domestic Substances List (Canada); ECx - Concentration associated with x% response; ELx - Loading rate associated with x% response; EmS - Emergency Schedule; ENCS - Existing and New Chemical Substances (Japan); ErCx - Concentration associated with x% growth rate response; ERG - Emergency Response Guide; GHS - Globally Harmonized System; GLP - Good Laboratory Practice; IARC - International Agency for Research on Cancer; IATA - International Air Transport Association; IBC - International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk; IC50 - Half maximal inhibitory concentration; ICAO - International Civil Aviation Organization; IECSC - Inventory of Existing Chemical Substances in China; IMDG - International Maritime Dangerous Goods; IMO - International Maritime Organization; ISHL - Industrial Safety and Health Law (Japan); ISO - International Organisation for Standardization; KECI - Korea Existing Chemicals Inventory; LC50 - Lethal Concentration to 50 % of a test population; LD50 - Lethal Dose to

Cell fixative



Version Revision Date: SDS Number: Date of last issue: 2016-10-27 1.13 2016-11-15 100000010702 Date of first issue: 2015-07-14

50% of a test population (Median Lethal Dose); MARPOL - International Convention for the Prevention of Pollution from Ships; n.o.s. - Not Otherwise Specified; Nch - Chilean Norm; NO(A)EC - No Observed (Adverse) Effect Concentration; NO(A)EL - No Observed (Adverse) Effect Level; NOELR - No Observable Effect Loading Rate; NOM - Official Mexican Norm; NTP - National Toxicology Program; NZIoC - New Zealand Inventory of Chemicals; OECD - Organization for Economic Co-operation and Development; OPPTS - Office of Chemical Safety and Pollution Prevention; PBT - Persistent, Bioaccumulative and Toxic substance; PICCS - Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR - (Quantitative) Structure Activity Relationship; REACH - Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; SADT - Self-Accelerating Decomposition Temperature; SDS - Safety Data Sheet; TCSI - Taiwan Chemical Substance Inventory; TDG - Transportation of Dangerous Goods; TSCA - Toxic Substances Control Act (United States); UN - United Nations; UNRTDG - United Nations Recommendations on the Transport of Dangerous Goods; vPvB - Very Persistent and Very Bioaccumulative; WHMIS - Workplace Hazardous Materials Information System

Date format : dd.mm.yyyy

Numbers : 123 456,78

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

IN / EN