Capture enhancement reagent



Version Revision Date: SDS Number: Date of last issue: 2016/10/27 1.13 2018/01/09 100000010878 Date of first issue: 2015/09/22

1. PRODUCT AND COMPANY IDENTIFICATION

Product name : Capture enhancement reagent

Substance name : Capture enhancement reagent

7037

Manufacturer or supplier's details

Company name of supplier : Menarini Silicon Biosystems, Inc.

Address : 3401 Masons Mill Rd #100

Huntingdon Valley, PA

19006, USA

Telephone : 1 (800) 381-4929

Responsible/issuing person

E-mail address

Us-info@siliconbiosystems.com

Emergency telephone : US : (303)-389-1805

number International: +1 (303)-389-1805

Recommended use of the chemical and restrictions on use

Recommended use : Assay reagent

2. HAZARDS IDENTIFICATION

GHS Classification

Not a dangerous substance or mixture according to the Globally Harmonised System (GHS).

GHS label elements

Not a dangerous substance or mixture according to the Globally Harmonised System (GHS).

Other hazards which do not result in classification

None known.

3. COMPOSITION/INFORMATION ON INGREDIENTS

Substance / Mixture : Mixture

Chemical nature : Liquid

Hazardous components

| Chemical name | CAS-No. | Concentration (% w/w) | ENCS/ISHL number |
|---------------|------------|-----------------------|------------------|
| sodium azide | 26628-22-8 | >= 0.1 - < 1 | 1-482 |

For explanation of abbreviations see section 16.

4. FIRST AID MEASURES

Capture enhancement reagent



 Version
 Revision Date:
 SDS Number:
 Date of last issue: 2016/10/27

 1.13
 2018/01/09
 100000010878
 Date of first issue: 2015/09/22

If inhaled : If breathed in, move person into fresh air.

Consult a physician.

In case of skin contact : Take off contaminated clothing and shoes immediately.

Wash off with plenty of water.

If symptoms persist, call a physician.

In case of eye contact : Rinse immediately with plenty of water, also under the eyelids,

for at least 5 minutes. Remove contact lenses.

If eye irritation persists, consult a specialist.

If swallowed, rinse mouth with water (only if the person is con-

scious).

Call a physician immediately.

Most important symptoms and effects, both acute and

delayed

: No information available.

Notes to physician : Treat symptomatically.

5. FIREFIGHTING MEASURES

Suitable extinguishing media : Use extinguishing measures that are appropriate to local cir-

cumstances and the surrounding environment.

Specific hazards during fire-

fighting

: No information available.

Specific extinguishing meth-

ods

: No information available.

Special protective equipment

for firefighters

: In the event of fire, wear self-contained breathing apparatus.

6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment and emer-

gency procedures

: In the event of an accidental release the emergency response team must respond based on a risk assessment and use personal protective equipment as appropriate.

: Should not be released into the environment.

Environmental precautions : Shoul

Methods and materials for containment and cleaning up

: Large spills: Dam up. Soak up with inert absorbent material.

Keep in properly labelled containers.

Small spills: Gently cover the spill with an absorbent towel or

nad

Large spills + Small spills: Keep in suitable, closed containers for disposal. Treat recovered material as described in the sec-

tion "Disposal considerations".

Capture enhancement reagent



Version 1.13

Revision Date: 2018/01/09

SDS Number: 100000010878

Date of last issue: 2016/10/27 Date of first issue: 2015/09/22

7. HANDLING AND STORAGE

Handling

Advice on protection against

fire and explosion

: No data available

Advice on safe handling : To avoid thermal decomposition, do not overheat.

Avoid inhalation, ingestion and contact with skin and eyes.

Use personal protective equipment as required.

Avoidance of contact : Oxidizing agents

Hygiene measures : Handle in accordance with good industrial hygiene and safety

practice.

Storage

Conditions for safe storage : To maintain product quality, do not store in heat or direct sun-

light.

Store in original container.

Keep containers tightly closed in a dry, cool and well-

ventilated place.

Keep away from heat and sources of ignition.

Keep locked up. Keep refrigerated.

Recommended storage tem-

perature

: 2-8°C

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Threshold limit value and permissible exposure limits for each component in the work environment

| Components | CAS-No. | Value type (Form of exposure) | Control parameters / Permissible concentration | Basis |
|--------------|------------|-------------------------------------|--|-------|
| sodium azide | 26628-22-8 | C (Vapour) | 0.11 ppm (Hydrazoic acid) | ACGIH |
| | | С | 0.29 mg/m3 (Sodium azide) | ACGIH |

Engineering measures

: All personal protective equipment should be based on a risk assessment. Consult a Environment Health Safety expert if

necessary.

Personal protective equipment

Respiratory protection : Engineering controls should always be the primary method of

controlling exposures.

If respiratory protective equipment is needed for certain activities, the type as well as the corresponding protection factor will depend upon the risk assessment and air concentrations, hazards, physical and warning properties of substances pre-

Capture enhancement reagent



Version Revision Date: SDS Number: Date of last issue: 2016/10/27 1.13 2018/01/09 100000010878 Date of first issue: 2015/09/22

sent.

No personal respiratory protective equipment normally re-

quired.

Hand protection

Remarks : Disposable gloves

Eye protection : No special precautions required.

Skin and body protection : No special precautions required.

Protective measures : The type of protective equipment must be selected based on

the Environmental Health and Safety risk assessment. Consult a Environmental Health and Safety expert if necessary.

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance : liquid

Colour : clear

Odour : No data available

10. STABILITY AND REACTIVITY

Reactivity : None reasonably foreseeable.

Chemical stability : Stable under recommended storage conditions.

Possibility of hazardous reac-

tions

: No dangerous reaction known under conditions of normal use.

Conditions to avoid : To avoid thermal decomposition, do not overheat.

Incompatible materials : Oxidizing agents

Hazardous decomposition

products

None known.

11. TOXICOLOGICAL INFORMATION

Acute toxicity

Product:

Acute oral toxicity : Acute toxicity estimate: > 2,000 mg/kg

Method: Calculation method

Acute dermal toxicity : Acute toxicity estimate: > 2,000 mg/kg

Method: Calculation method

Components:

Capture enhancement reagent



Version 1.13 Revision Date: 2018/01/09

SDS Number: 100000010878

Date of last issue: 2016/10/27 Date of first issue: 2015/09/22

sodium azide:

Acute oral toxicity : LD50 (Rat): 27 mg/kg

Skin corrosion/irritation

No data available

Serious eye damage/eye irritation

No data available

Respiratory or skin sensitisation

No data available

Germ cell mutagenicity

No data available

Carcinogenicity

No data available

Reproductive toxicity

No data available

STOT - single exposure

No data available

STOT - repeated exposure

No data available

Repeated dose toxicity

No data available

Aspiration toxicity

No data available

12. ECOLOGICAL INFORMATION

Ecotoxicity

Components:

sodium azide:

Toxicity to fish : LC50 (Lepomis macrochirus (Bluegill sunfish)): 0.7 mg/l

Exposure time: 96 h

Toxicity to daphnia and other

aquatic invertebrates

: EC50 (Daphnia pulex (Water flea)): 4.2 mg/l

Exposure time: 96 h

Toxicity to algae : IC50: 272 mg/l

Toxicity to bacteria : EC50 (Photobacterium phosphoreum): 38.5 mg/l

Persistence and degradability

No data available

Capture enhancement reagent



Version 1.13

Revision Date: 2018/01/09

SDS Number: 100000010878

Date of last issue: 2016/10/27 Date of first issue: 2015/09/22

Bioaccumulative potential

No data available

Mobility in soil

No data available

Hazardous to the ozone layer

Not applicable

Other adverse effects

No data available

13. DISPOSAL CONSIDERATIONS

Disposal methods

Waste from residues : In accordance with National, Federal, State and Local regula-

tions.

Contaminated packaging : Empty containers should be taken to an approved waste han-

dling site for recycling or disposal.

14. TRANSPORT INFORMATION

International Regulations

UNRTDG

Not regulated as a dangerous good

IATA-DGR

Not regulated as a dangerous good

IMDG-Code

Not regulated as a dangerous good

Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable for product as supplied.

National Regulations

Refer to section 15 for specific national regulation.

15. REGULATORY INFORMATION

Related Regulations

Fire Service Law

Not applicable to dangerous materials / designated flammables.

Chemical Substance Control Law

Not applicable for Specified Chemical Substance, Monitoring Chemical Substance and Priority Assessment Chemical Substance.

Capture enhancement reagent



Version 1.13

Revision Date: 2018/01/09

SDS Number: 100000010878

Date of last issue: 2016/10/27 Date of first issue: 2015/09/22

Industrial Safety and Health Law

Harmful Substances Prohibited from Manufacture

Not applicable

Harmful Substances Required Permission for Manufacture

Not applicable

Substances Prevented From Impairment of Health

Not applicable

Circular concerning Information on Chemicals having Mutagenicity - Annex 2: Information on Existing Chemicals having Mutagenicity

Not applicable

Circular concerning Information on Chemicals having Mutagenicity - Annex 1: Information on Notified Substances having Mutagenicity

Not applicable

Substances Subject to be Notified Names

Not applicable

Substances Subject to be Indicated Names

Not applicable

Ordinance on Prevention of Hazards Due to Specified Chemical Substances

Not applicable

Ordinance on Prevention of Lead Poisoning

Not applicable

Ordinance on Prevention of Tetraalkyl Lead Poisoning

Not applicable

Ordinance on Prevention of Organic Solvent Poisoning

Not applicable

Enforcement Order of the Industrial Safety and Health Law - Attached table 1 (Dangerous Substances)

Not applicable

Poisonous and Deleterious Substances Control Law

Not applicable

Act on Confirmation, etc. of Release Amounts of Specific Chemical Substances in the Environment and Promotion of Improvements to the Management Thereof

Not applicable

High Pressure Gas Safety Act

Not applicable

Explosive Control Law

Not applicable

Vessel Safety Law

Not regulated as a dangerous good

Aviation Law

Not regulated as a dangerous good

Marine Pollution and Sea Disaster Prevention etc Law

Bulk transportation : Not classified as noxious liquid substance

Pack transportation : Not classified as marine pollutant

Capture enhancement reagent



Version 1.13

Revision Date: 2018/01/09

SDS Number: 100000010878

Date of last issue: 2016/10/27 Date of first issue: 2015/09/22

Waste Disposal and Public Cleansing Law

Industrial waste

16. OTHER INFORMATION

Full text of other abbreviations

AICS - Australian Inventory of Chemical Substances; ANTT - National Agency for Transport by Land of Brazil; ASTM - American Society for the Testing of Materials; bw - Body weight; CMR -Carcinogen, Mutagen or Reproductive Toxicant; CPR - Controlled Products Regulations; DIN -Standard of the German Institute for Standardisation; DSL - Domestic Substances List (Canada); ECx - Concentration associated with x% response; ELx - Loading rate associated with x% response; EmS - Emergency Schedule; ENCS - Existing and New Chemical Substances (Japan); ErCx - Concentration associated with x% growth rate response; ERG - Emergency Response Guide; GHS - Globally Harmonized System; GLP - Good Laboratory Practice; IARC - International Agency for Research on Cancer; IATA - International Air Transport Association; IBC - International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk; IC50 - Half maximal inhibitory concentration; ICAO - International Civil Aviation Organization; IECSC - Inventory of Existing Chemical Substances in China; IMDG - International Maritime Dangerous Goods; IMO - International Maritime Organization; ISHL - Industrial Safety and Health Law (Japan): ISO - International Organisation for Standardization: KECI - Korea Existing Chemicals Inventory; LC50 - Lethal Concentration to 50 % of a test population; LD50 - Lethal Dose to 50% of a test population (Median Lethal Dose); MARPOL - International Convention for the Prevention of Pollution from Ships; n.o.s. - Not Otherwise Specified; Nch - Chilean Norm; NO(A)EC -No Observed (Adverse) Effect Concentration; NO(A)EL - No Observed (Adverse) Effect Level; NOELR - No Observable Effect Loading Rate; NOM - Official Mexican Norm; NTP - National Toxicology Program; NZIoC - New Zealand Inventory of Chemicals; OECD - Organization for Economic Co-operation and Development; OPPTS - Office of Chemical Safety and Pollution Prevention; PBT - Persistent, Bioaccumulative and Toxic substance; PICCS - Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR - (Quantitative) Structure Activity Relationship; REACH - Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; SADT - Self-Accelerating Decomposition Temperature: SDS - Safety Data Sheet: TCSI - Taiwan Chemical Substance Inventory; TDG - Transportation of Dangerous Goods; TSCA - Toxic Substances Control Act (United States); UN - United Nations; UNRTDG - United Nations Recommendations on the Transport of Dangerous Goods; vPvB - Very Persistent and Very Bioaccumulative; WHMIS -Workplace Hazardous Materials Information System

Date format : yyyy/mm/dd Numbers 123,456.78

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

JP / EN

Staining Reagent



Version Revision Date: SDS Number: Date of last issue: 2016/11/15
1.1 2018/01/09 100000013033 Date of first issue: 2016/11/15

1. PRODUCT AND COMPANY IDENTIFICATION

Substance name : Staining Reagent

Manufacturer or supplier's details

Company name of supplier : Menarini Silicon Biosystems, Inc.

Address : 3401 Masons Mill Rd #100

Huntingdon Valley, PA 19006,

Us-info@siliconbiosystems.com

USA

Telephone : 1 (800) 381-4929

Responsible/issuing person

E-mail address

Emergency telephone : US : (303)-389-1805

number International: +1 (303)-389-1805

Recommended use of the chemical and restrictions on use

Recommended use : Assay reagent

2. HAZARDS IDENTIFICATION

GHS Classification

Not a dangerous substance or mixture according to the Globally Harmonised System (GHS). Not a dangerous substance or mixture according to the Globally Harmonised System (GHS).

GHS label elements

Not a dangerous substance or mixture according to the Globally Harmonised System (GHS). Not a dangerous substance or mixture according to the Globally Harmonised System (GHS).

Other hazards which do not result in classification

None known.

3. COMPOSITION/INFORMATION ON INGREDIENTS

Substance / Mixture : Mixture

Chemical nature : Liquid

Hazardous components

| Chemical name | CAS-No. | Concentration (% | ENCS/ISHL |
|---------------|------------|------------------|-----------|
| | | w/w) | number |
| sodium azide | 26628-22-8 | >= 0.1 - < 1 | 1-482 |

For explanation of abbreviations see section 16.

4. FIRST AID MEASURES



Version Revision Date: SDS Number: Date of last issue: 2016/11/15 2018/01/09 100000013033 1.1 Date of first issue: 2016/11/15

If inhaled : If breathed in, move person into fresh air.

Consult a physician.

In case of skin contact : Take off contaminated clothing and shoes immediately.

Wash off with plenty of water.

If symptoms persist, call a physician.

In case of eye contact : Rinse immediately with plenty of water, also under the eyelids,

> for at least 5 minutes. Remove contact lenses.

If eye irritation persists, consult a specialist.

If swallowed : If swallowed, rinse mouth with water (only if the person is con-

scious).

Call a physician immediately.

Most important symptoms and effects, both acute and

delayed

: No information available.

: Treat symptomatically. Notes to physician

5. FIREFIGHTING MEASURES

Suitable extinguishing media Use extinguishing measures that are appropriate to local cir-

cumstances and the surrounding environment.

Specific hazards during fire-

fighting

: No information available.

Hazardous combustion prod-

ucts

: No hazardous combustion products are known

Specific extinguishing meth-

: No information available.

Special protective equipment

for firefighters

: In the event of fire, wear self-contained breathing apparatus.

6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment and emergency procedures

: In the event of an accidental release the emergency response team must respond based on a risk assessment and use personal protective equipment as appropriate.

Environmental precautions : Should not be released into the environment.

Methods and materials for containment and cleaning up : Large spills: Dam up. Soak up with inert absorbent material.

Keep in properly labelled containers.

Small spills: Gently cover the spill with an absorbent towel or

Large spills + Small spills: Keep in suitable, closed containers for disposal. Treat recovered material as described in the sec-



Version 1.1

Revision Date: 2018/01/09

SDS Number: 100000013033

Date of last issue: 2016/11/15 Date of first issue: 2016/11/15

tion "Disposal considerations".

7. HANDLING AND STORAGE

Handling

Advice on protection against

fire and explosion

: No data available

Advice on safe handling

: To avoid thermal decomposition, do not overheat.

Avoid inhalation, ingestion and contact with skin and eyes.

Use personal protective equipment as required.

Avoidance of contact

: Strong acids and strong bases

Reducing agents Oxidizing agents

Hygiene measures

: Handle in accordance with good industrial hygiene and safety

practice.

Storage

Conditions for safe storage

: To maintain product quality, do not store in heat or direct sun-

light.

Store in original container.

Keep containers tightly closed in a dry, cool and well-

ventilated place.

Keep away from heat and sources of ignition.

Keep locked up. Keep refrigerated.

Recommended storage tem-

perature

: 2-8°C

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Threshold limit value and permissible exposure limits for each component in the work environment

| Components | CAS-No. | Value type (Form of exposure) | Control parameters / Permissible concentration | Basis |
|--------------|------------|-------------------------------------|--|-------|
| sodium azide | 26628-22-8 | C (Vapour) | 0.11 ppm (Hydrazoic acid) | ACGIH |
| | | С | 0.29 mg/m3 (Sodium azide) | ACGIH |

Engineering measures

: All personal protective equipment should be based on a risk assessment. Consult a Environment Health Safety expert if

necessary.

Personal protective equipment

Respiratory protection : Engineering controls should always be the primary method of



Version Revision Date: SDS Number: Date of last issue: 2016/11/15
1.1 2018/01/09 100000013033 Date of first issue: 2016/11/15

controlling exposures.

If respiratory protective equipment is needed for certain activities, the type as well as the corresponding protection factor will depend upon the risk assessment and air concentrations, hazards, physical and warning properties of substances pre-

sent.

No personal respiratory protective equipment normally re-

quired.

Eye protection : No special precautions required.

Skin and body protection : No special precautions required.

Protective measures : The type of protective equipment must be selected based on

the Environmental Health and Safety risk assessment. Consult a Environmental Health and Safety expert if necessary.

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance : liquid

Colour : purple

Odour : odourless

pH : 7.5

Solubility(ies)

Water solubility : soluble

10. STABILITY AND REACTIVITY

Reactivity : None reasonably foreseeable.

Chemical stability : Stable under recommended storage conditions.

Possibility of hazardous reac-

tions

No dangerous reaction known under conditions of normal use.

Conditions to avoid : To avoid thermal decomposition, do not overheat.

Incompatible materials : Strong acids and strong bases

Reducing agents Oxidizing agents

Hazardous decomposition

products

: None known.



Version 1.1 Revision Date: 2018/01/09

SDS Number: 100000013033

Date of last issue: 2016/11/15 Date of first issue: 2016/11/15

11. TOXICOLOGICAL INFORMATION

Acute toxicity

Product:

Acute oral toxicity : Acute toxicity estimate: > 2,000 mg/kg

Method: Calculation method

Acute toxicity estimate: > 2,000 mg/kg

Method: Calculation method

Acute dermal toxicity : Acute toxicity estimate: > 2,000 mg/kg

Method: Calculation method

Acute toxicity estimate: > 2,000 mg/kg

Method: Calculation method

Components:

sodium azide:

Acute oral toxicity : LD50 (Rat): 27 mg/kg

Skin corrosion/irritation

No data available

Serious eye damage/eye irritation

No data available

Respiratory or skin sensitisation

No data available

Germ cell mutagenicity

No data available

Carcinogenicity

No data available

Reproductive toxicity

No data available

STOT - single exposure

No data available

STOT - repeated exposure

No data available

Repeated dose toxicity

No data available

Aspiration toxicity

No data available



Version Revision Date: SDS Number: Date of last issue: 2016/11/15
1.1 2018/01/09 100000013033 Date of first issue: 2016/11/15

12. ECOLOGICAL INFORMATION

Ecotoxicity

Components:

sodium azide:

Toxicity to fish : LC50 (Lepomis macrochirus (Bluegill sunfish)): 0.7 mg/l

Exposure time: 96 h

Toxicity to daphnia and other

aquatic invertebrates

: EC50 (Daphnia pulex (Water flea)): 4.2 mg/l

Exposure time: 96 h

Toxicity to algae : IC50: 272 mg/l

Toxicity to bacteria : EC50 (Photobacterium phosphoreum): 38.5 mg/l

Persistence and degradability

No data available

Bioaccumulative potential

No data available

Mobility in soil

No data available

Hazardous to the ozone layer

Not applicable

Other adverse effects

No data available

13. DISPOSAL CONSIDERATIONS

Disposal methods

Waste from residues : In accordance with National, Federal, State and Local regula-

tions.

Contaminated packaging : Empty containers should be taken to an approved waste han-

dling site for recycling or disposal.

14. TRANSPORT INFORMATION

International Regulations

UNRTDG

Not regulated as a dangerous good

IATA-DGR

Not regulated as a dangerous good

IMDG-Code



Version Revision Date: SDS Number: Date of last issue: 2016/11/15
1.1 2018/01/09 100000013033 Date of first issue: 2016/11/15

Not regulated as a dangerous good

Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable for product as supplied.

National Regulations

Refer to section 15 for specific national regulation.

15. REGULATORY INFORMATION

Related Regulations

Fire Service Law

Not applicable to dangerous materials / designated flammables.

Chemical Substance Control Law

Not applicable for Specified Chemical Substance, Monitoring Chemical Substance and Priority Assessment Chemical Substance.

Not applicable for Specified Chemical Substance, Monitoring Chemical Substance and Priority Assessment Chemical Substance.

Industrial Safety and Health Law

Harmful Substances Prohibited from Manufacture

Not applicable

Harmful Substances Prohibited from Manufacture

Not applicable

Harmful Substances Required Permission for Manufacture

Not applicable

Harmful Substances Required Permission for Manufacture

Not applicable

Substances Prevented From Impairment of Health

Not applicable

Substances Prevented From Impairment of Health

Not applicable

Circular concerning Information on Chemicals having Mutagenicity - Annex 2: Information on Existing Chemicals having Mutagenicity

Not applicable

Circular concerning Information on Chemicals having Mutagenicity - Annex 2: Information on Existing Chemicals having Mutagenicity

Not applicable

Circular concerning Information on Chemicals having Mutagenicity - Annex 1: Information on Notified Substances having Mutagenicity

Not applicable

Circular concerning Information on Chemicals having Mutagenicity - Annex 1: Information on Notified Substances having Mutagenicity

Not applicable



Version Revision Date: SDS Number: Date of last issue: 2016/11/15
1.1 2018/01/09 100000013033 Date of first issue: 2016/11/15

Substances Subject to be Notified Names

Not applicable Not applicable

Substances Subject to be Indicated Names

Not applicable

Substances Subject to be Indicated Names

Not applicable

Ordinance on Prevention of Hazards Due to Specified Chemical Substances

Not applicable

Ordinance on Prevention of Hazards Due to Specified Chemical Substances

Not applicable

Ordinance on Prevention of Lead Poisoning

Not applicable

Ordinance on Prevention of Lead Poisoning

Not applicable

Ordinance on Prevention of Tetraalkyl Lead Poisoning

Not applicable

Ordinance on Prevention of Tetraalkyl Lead Poisoning

Not applicable

Ordinance on Prevention of Organic Solvent Poisoning

Not applicable

Ordinance on Prevention of Organic Solvent Poisoning

Not applicable

Enforcement Order of the Industrial Safety and Health Law - Attached table 1 (Dangerous Substances)

Not applicable

Enforcement Order of the Industrial Safety and Health Law - Attached table 1 (Dangerous Substances)

Not applicable

Poisonous and Deleterious Substances Control Law

Not applicable

Not applicable

Act on Confirmation, etc. of Release Amounts of Specific Chemical Substances in the Environment and Promotion of Improvements to the Management Thereof

Not applicable

Act on Confirmation, etc. of Release Amounts of Specific Chemical Substances in the Environment and Promotion of Improvements to the Management Thereof

Not applicable

High Pressure Gas Safety Act

Not applicable

Explosive Control Law

Not applicable

Vessel Safety Law

Not regulated as a dangerous good

Aviation Law

Not regulated as a dangerous good



Version Revision Date: SDS Number: Date of last issue: 2016/11/15
1.1 2018/01/09 100000013033 Date of first issue: 2016/11/15

Explosive Control Law

Not applicable

Vessel Safety Law

Not regulated as a dangerous good

Aviation Law

Not regulated as a dangerous good

Marine Pollution and Sea Disaster Prevention etc Law

Bulk transportation : Not classified as noxious liquid substance

Pack transportation : Not classified as marine pollutant

Bulk transportation : Not classified as noxious liquid substance

Pack transportation : Not classified as marine pollutant

Waste Disposal and Public Cleansing Law

Industrial waste

Industrial waste

16. OTHER INFORMATION

Full text of other abbreviations

AICS - Australian Inventory of Chemical Substances; ANTT - National Agency for Transport by Land of Brazil; ASTM - American Society for the Testing of Materials; bw - Body weight; CMR -Carcinogen, Mutagen or Reproductive Toxicant; CPR - Controlled Products Regulations; DIN -Standard of the German Institute for Standardisation: DSL - Domestic Substances List (Canada): ECx - Concentration associated with x% response; ELx - Loading rate associated with x% response; EmS - Emergency Schedule; ENCS - Existing and New Chemical Substances (Japan); ErCx - Concentration associated with x% growth rate response; ERG - Emergency Response Guide; GHS - Globally Harmonized System; GLP - Good Laboratory Practice; IARC - International Agency for Research on Cancer; IATA - International Air Transport Association; IBC - International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk; IC50 - Half maximal inhibitory concentration; ICAO - International Civil Aviation Organization; IECSC - Inventory of Existing Chemical Substances in China; IMDG - International Maritime Dangerous Goods; IMO - International Maritime Organization; ISHL - Industrial Safety and Health Law (Japan): ISO - International Organisation for Standardization: KECI - Korea Existing Chemicals Inventory; LC50 - Lethal Concentration to 50 % of a test population; LD50 - Lethal Dose to 50% of a test population (Median Lethal Dose); MARPOL - International Convention for the Prevention of Pollution from Ships; n.o.s. - Not Otherwise Specified; Nch - Chilean Norm; NO(A)EC -No Observed (Adverse) Effect Concentration; NO(A)EL - No Observed (Adverse) Effect Level; NOELR - No Observable Effect Loading Rate; NOM - Official Mexican Norm; NTP - National Toxicology Program; NZIoC - New Zealand Inventory of Chemicals; OECD - Organization for Economic Co-operation and Development; OPPTS - Office of Chemical Safety and Pollution Prevention; PBT - Persistent, Bioaccumulative and Toxic substance; PICCS - Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR - (Quantitative) Structure Activity Relationship; REACH - Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; SADT - Self-Accelerating Decomposition Temperature; SDS - Safety Data Sheet; TCSI - Taiwan Chemical Substance Inventory; TDG - Transportation of Dangerous Goods; TSCA - Toxic Substances Con-



 Version
 Revision Date:
 SDS Number:
 Date of last issue: 2016/11/15

 1.1
 2018/01/09
 100000013033
 Date of first issue: 2016/11/15

trol Act (United States); UN - United Nations; UNRTDG - United Nations Recommendations on the Transport of Dangerous Goods; vPvB - Very Persistent and Very Bioaccumulative; WHMIS - Workplace Hazardous Materials Information System

Date format : yyyy/mm/dd Numbers 123,456.78

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

JP / EN

Nucleic acid dye



Version Revision Date: SDS Number: Date of last issue: 2016/10/27 1.9 2018/01/09 100000010877 Date of first issue: 2015/09/21

1. PRODUCT AND COMPANY IDENTIFICATION

Product name : Nucleic acid dye

Substance name : Nucleic acid dye

7041

Manufacturer or supplier's details

Company name of supplier : Menarini Silicon Biosystems, Inc.

Address : 3401 Masons Mill Rd #100

Huntingdon Valley, PA 19006,

USA

Telephone : 1 (800) 381-4929

E-mail address

Responsible/issuing person

Us-info@siliconbiosystems.com

Emergency telephone

number

: US: (303)-389-1805

International: +1 (303)-389-1805

Recommended use of the chemical and restrictions on use

Recommended use : Assay reagent

2. HAZARDS IDENTIFICATION

GHS Classification

Not a dangerous substance or mixture according to the Globally Harmonised System (GHS).

GHS label elements

Not a dangerous substance or mixture according to the Globally Harmonised System (GHS).

Other hazards which do not result in classification

None known.

3. COMPOSITION/INFORMATION ON INGREDIENTS

Substance / Mixture : Mixture

Chemical nature : Liquid

Hazardous components

This product does not contain chemical(s) listed by Japanese regulations and/or classified at or above the cut-off value(s) according to GHS.

4. FIRST AID MEASURES

If inhaled : If breathed in, move person into fresh air.

Nucleic acid dye



Version Revision Date: SDS Number: Date of last issue: 2016/10/27 1.9 2018/01/09 100000010877 Date of first issue: 2015/09/21

Consult a physician.

In case of skin contact : Take off contaminated clothing and shoes immediately.

Wash off immediately with plenty of water. If symptoms persist, call a physician. Wash contaminated clothing before re-use.

In case of eye contact : Rinse immediately with plenty of water, also under the eyelids,

for at least 5 minutes. Remove contact lenses.

If eye irritation persists, consult a specialist.

If swallowed, rinse mouth with water (only if the person is con-

scious).

Call a physician immediately.

Most important symptoms and effects, both acute and

delayed

: No information available.

Notes to physician : Treat symptomatically.

5. FIREFIGHTING MEASURES

Suitable extinguishing media : Use extinguishing measures that are appropriate to local cir-

cumstances and the surrounding environment.

Specific hazards during fire-

fighting

: No information available.

Specific extinguishing meth-

ods

: No information available.

Special protective equipment

for firefighters

: In the event of fire, wear self-contained breathing apparatus.

6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment and emer-

gency procedures

: In the event of an accidental release the emergency response team must respond based on a risk assessment and use per-

sonal protective equipment as appropriate.

Environmental precautions : Should not be released into the environment.

Methods and materials for containment and cleaning up

: Large spills: Dam up. Soak up with inert absorbent material.

Keep in properly labelled containers.

Small spills: Gently cover the spill with an absorbent towel or

nad

Large spills + Small spills: Keep in suitable, closed containers for disposal. Treat recovered material as described in the sec-

tion "Disposal considerations".

Nucleic acid dye



Version 1.9

Revision Date: 2018/01/09

SDS Number: 100000010877

Date of last issue: 2016/10/27 Date of first issue: 2015/09/21

7. HANDLING AND STORAGE

Handling

Advice on protection against

fire and explosion

: No data available

Advice on safe handling : To avoid thermal decomposition, do not overheat.

Avoid inhalation, ingestion and contact with skin and eyes.

Use personal protective equipment as required.

Avoidance of contact : None known.

Hygiene measures : Handle in accordance with good industrial hygiene and safety

practice.

Storage

Conditions for safe storage : To maintain product quality, do not store in heat or direct sun-

light.

Store in original container.

Keep containers tightly closed in a dry, cool and well-

ventilated place.

Keep away from heat and sources of ignition.

Keep locked up. Keep refrigerated.

Recommended storage tem-

perature

: 2-8°C

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Threshold limit value and permissible exposure limits for each component in the work environment

Contains no substances with occupational exposure limit values.

Engineering measures : All personal protective equipment should be based on a risk

assessment. Consult a Environment Health Safety expert if

necessary.

Personal protective equipment

Respiratory protection : Engineering controls should always be the primary method of

controlling exposures.

If respiratory protective equipment is needed for certain activities, the type as well as the corresponding protection factor will depend upon the risk assessment and air concentrations, hazards, physical and warning properties of substances pre-

sent.

No personal respiratory protective equipment normally re-

quired.

Hand protection

Nucleic acid dye



Version Revision Date: SDS Number: Date of last issue: 2016/10/27 1.9 2018/01/09 100000010877 Date of first issue: 2015/09/21

Remarks : Disposable gloves

Eye protection : No special precautions required.

Skin and body protection : No special precautions required.

Protective measures : The type of protective equipment must be selected based on

the Environmental Health and Safety risk assessment. Consult a Environmental Health and Safety expert if necessary.

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance : liquid

Colour : clear, light yellow

Odour : odourless

pH : 7.5

Solubility(ies)

Water solubility : soluble

10. STABILITY AND REACTIVITY

Reactivity : None reasonably foreseeable.

Chemical stability : Stable under recommended storage conditions.

Possibility of hazardous reac-

tions

: No dangerous reaction known under conditions of normal use.

Conditions to avoid : To avoid thermal decomposition, do not overheat.

Incompatible materials : None known.

Hazardous decomposition

products

: None known.

11. TOXICOLOGICAL INFORMATION

Acute toxicity

No data available

Skin corrosion/irritation

No data available

Serious eye damage/eye irritation

No data available

Respiratory or skin sensitisation

No data available

Nucleic acid dye



Version 1.9

Revision Date: 2018/01/09

SDS Number: 100000010877

Date of last issue: 2016/10/27 Date of first issue: 2015/09/21

Germ cell mutagenicity

No data available

Carcinogenicity

No data available

Reproductive toxicity

No data available

STOT - single exposure

No data available

STOT - repeated exposure

No data available

Repeated dose toxicity

No data available

Aspiration toxicity

No data available

12. ECOLOGICAL INFORMATION

Ecotoxicity

No data available

Persistence and degradability

No data available

Bioaccumulative potential

No data available

Mobility in soil

No data available

Hazardous to the ozone layer

Not applicable

Other adverse effects

No data available

13. DISPOSAL CONSIDERATIONS

Disposal methods

Waste from residues : In accordance with National, Federal, State and Local regula-

tions.

Contaminated packaging : Empty containers should be taken to an approved waste han-

dling site for recycling or disposal.

Nucleic acid dye



Version 1.9

Revision Date: 2018/01/09

SDS Number: 100000010877

Date of last issue: 2016/10/27 Date of first issue: 2015/09/21

14. TRANSPORT INFORMATION

International Regulations

UNRTDG

Not regulated as a dangerous good

IATA-DGR

Not regulated as a dangerous good

IMDG-Code

Not regulated as a dangerous good

Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable for product as supplied.

National Regulations

Refer to section 15 for specific national regulation.

15. REGULATORY INFORMATION

Related Regulations

Fire Service Law

Not applicable

Chemical Substance Control Law

Not applicable for Specified Chemical Substance, Monitoring Chemical Substance and Priority Assessment Chemical Substance.

Industrial Safety and Health Law

Harmful Substances Prohibited from Manufacture

Not applicable

Harmful Substances Required Permission for Manufacture

Not applicable

Substances Prevented From Impairment of Health

Not applicable

Circular concerning Information on Chemicals having Mutagenicity - Annex 2: Information on Existing Chemicals having Mutagenicity

Not applicable

Circular concerning Information on Chemicals having Mutagenicity - Annex 1: Information on Notified Substances having Mutagenicity

Not applicable

Substances Subject to be Notified Names

Not applicable

Substances Subject to be Indicated Names

Not applicable

Ordinance on Prevention of Hazards Due to Specified Chemical Substances

Nucleic acid dye



Version Revision Date: SDS Number: Date of last issue: 2016/10/27 1.9 2018/01/09 100000010877 Date of first issue: 2015/09/21

Not applicable

Ordinance on Prevention of Lead Poisoning

Not applicable

Ordinance on Prevention of Tetraalkyl Lead Poisoning

Not applicable

Ordinance on Prevention of Organic Solvent Poisoning

Not applicable

Poisonous and Deleterious Substances Control Law

Not applicable

Act on Confirmation, etc. of Release Amounts of Specific Chemical Substances in the Environment and Promotion of Improvements to the Management Thereof

Not applicable

High Pressure Gas Safety Act

Not applicable

16. OTHER INFORMATION

Full text of other abbreviations

AICS - Australian Inventory of Chemical Substances; ANTT - National Agency for Transport by Land of Brazil; ASTM - American Society for the Testing of Materials; bw - Body weight; CMR -Carcinogen, Mutagen or Reproductive Toxicant; CPR - Controlled Products Regulations; DIN -Standard of the German Institute for Standardisation; DSL - Domestic Substances List (Canada); ECx - Concentration associated with x% response; ELx - Loading rate associated with x% response; EmS - Emergency Schedule; ENCS - Existing and New Chemical Substances (Japan); ErCx - Concentration associated with x% growth rate response; ERG - Emergency Response Guide; GHS - Globally Harmonized System; GLP - Good Laboratory Practice; IARC - International Agency for Research on Cancer; IATA - International Air Transport Association; IBC - International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk; IC50 - Half maximal inhibitory concentration: ICAO - International Civil Aviation Organization: IECSC - Inventory of Existing Chemical Substances in China; IMDG - International Maritime Dangerous Goods; IMO - International Maritime Organization; ISHL - Industrial Safety and Health Law (Japan); ISO - International Organisation for Standardization; KECI - Korea Existing Chemicals Inventory; LC50 - Lethal Concentration to 50 % of a test population; LD50 - Lethal Dose to 50% of a test population (Median Lethal Dose); MARPOL - International Convention for the Prevention of Pollution from Ships; n.o.s. - Not Otherwise Specified; Nch - Chilean Norm; NO(A)EC -No Observed (Adverse) Effect Concentration; NO(A)EL - No Observed (Adverse) Effect Level; NOELR - No Observable Effect Loading Rate; NOM - Official Mexican Norm; NTP - National Toxicology Program; NZIoC - New Zealand Inventory of Chemicals; OECD - Organization for Economic Co-operation and Development; OPPTS - Office of Chemical Safety and Pollution Prevention; PBT - Persistent, Bioaccumulative and Toxic substance; PICCS - Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR - (Quantitative) Structure Activity Relationship; REACH - Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; SADT - Self-Accelerating Decomposition Temperature; SDS - Safety Data Sheet; TCSI - Taiwan Chemical Substance Inventory; TDG - Transportation of Dangerous Goods; TSCA - Toxic Substances Control Act (United States); UN - United Nations; UNRTDG - United Nations Recommendations on the Transport of Dangerous Goods; vPvB - Very Persistent and Very Bioaccumulative; WHMIS -Workplace Hazardous Materials Information System

Date format : yyyy/mm/dd Numbers 123,456.78

Nucleic acid dye



Version Revision Date: SDS Number: Date of last issue: 2016/10/27 1.9 2018/01/09 100000010877 Date of first issue: 2015/09/21

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

JP / EN

Permeabilization reagent



Version Revision Date: SDS Number: Date of last issue: 2016/10/27 1.12 2018/01/09 100000010887 Date of first issue: 2015/09/21

1. PRODUCT AND COMPANY IDENTIFICATION

Product name : Permeabilization reagent

Substance name : Permeabilization reagent

7038

Manufacturer or supplier's details

Company name of supplier : Menarini Silicon Biosystems, Inc.

Address : 3401 Masons Mill Rd #100

Huntingdon Valley, PA 19006,

USA

Telephone : 1 (800) 381-4929

E-mail address Responsible/issuing person Us-info@siliconbiosystems.com

Emergency telephone

number

: US: (303)-389-1805 International: +1 (303)-389-1805

Recommended use of the chemical and restrictions on use

Recommended use : Assay reagent

2. HAZARDS IDENTIFICATION

GHS Classification

Not a dangerous substance or mixture according to the Globally Harmonised System (GHS).

GHS label elements

Not a dangerous substance or mixture according to the Globally Harmonised System (GHS).

Other hazards which do not result in classification

None known.

3. COMPOSITION/INFORMATION ON INGREDIENTS

Substance / Mixture : Mixture

Chemical nature : Liquid

Hazardous components

| Chemical name | CAS-No. | Concentration (% | ENCS/ISHL |
|---------------|------------|------------------|-----------|
| | | w/w) | number |
| sodium azide | 26628-22-8 | >= 0.1 - < 1 | 1-482 |

For explanation of abbreviations see section 16.

4. FIRST AID MEASURES

Permeabilization reagent



Version Revision Date: SDS Number: Date of last issue: 2016/10/27 1.12 2018/01/09 100000010887 Date of first issue: 2015/09/21

If inhaled : If breathed in, move person into fresh air.

Consult a physician.

In case of skin contact : Take off contaminated clothing and shoes immediately.

Wash off with plenty of water.

If symptoms persist, call a physician.

In case of eye contact : Remove contact lenses.

If eye irritation persists, consult a specialist.

Rinse immediately with plenty of water, also under the eyelids,

for at least 5 minutes.

If swallowed, rinse mouth with water (only if the person is con-

scious).

Call a physician immediately.

Most important symptoms and effects, both acute and

delayed

: No information available.

Notes to physician : Treat symptomatically.

5. FIREFIGHTING MEASURES

Suitable extinguishing media : Use extinguishing measures that are appropriate to local cir-

cumstances and the surrounding environment.

Specific hazards during fire-

fighting

: No information available.

Specific extinguishing meth-

ods

: No information available.

Special protective equipment

for firefighters

: In the event of fire, wear self-contained breathing apparatus.

6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment and emer-

gency procedures

: In the event of an accidental release the emergency response team must respond based on a risk assessment and use per-

sonal protective equipment as appropriate.

Environmental precautions : Should not be released into the environment.

Methods and materials for containment and cleaning up

: Large spills: Dam up. Soak up with inert absorbent material.

Keep in properly labelled containers.

Small spills: Gently cover the spill with an absorbent towel or

pad.

Large spills + Small spills: Keep in suitable, closed containers for disposal. Treat recovered material as described in the sec-

tion "Disposal considerations".

Permeabilization reagent



Version 1.12

Revision Date: 2018/01/09

SDS Number: 100000010887

Date of last issue: 2016/10/27 Date of first issue: 2015/09/21

7. HANDLING AND STORAGE

Handling

Advice on protection against

fire and explosion

: No data available

Advice on safe handling : To avoid thermal decomposition, do not overheat.

Avoid inhalation, ingestion and contact with skin and eyes.

Use personal protective equipment as required.

Avoidance of contact : Strong acids and strong bases

Reducing agents Strong oxidizing agents

Hygiene measures : Handle in accordance with good industrial hygiene and safety

practice.

Storage

Conditions for safe storage : To maintain product quality, do not store in heat or direct sun-

light.

Store in original container.

Keep containers tightly closed in a dry, cool and well-

ventilated place.

Keep away from heat and sources of ignition.

Keep locked up. Keep refrigerated.

Recommended storage tem-

perature

: 2-8°C

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Threshold limit value and permissible exposure limits for each component in the work environment

| Components | CAS-No. | Value type | Control parame- | Basis |
|--------------|------------|------------|--------------------|-------|
| | | (Form of | ters / Permissible | |
| | | exposure) | concentration | |
| sodium azide | 26628-22-8 | C (Vapour) | 0.11 ppm | ACGIH |
| | | | (Hydrazoic acid) | |
| | | С | 0.29 mg/m3 | ACGIH |
| | | | (Sodium azide) | |

Engineering measures : All personal protective equipment should be based on a risk

assessment. Consult a Environment Health Safety expert if

necessary.

Personal protective equipment

Respiratory protection : Engineering controls should always be the primary method of

controlling exposures.

If respiratory protective equipment is needed for certain activities, the type as well as the corresponding protection factor

Permeabilization reagent



Version Revision Date: SDS Number: Date of last issue: 2016/10/27 1.12 2018/01/09 100000010887 Date of first issue: 2015/09/21

will depend upon the risk assessment and air concentrations, hazards, physical and warning properties of substances pre-

sent.

No personal respiratory protective equipment normally re-

quired.

Hand protection

Remarks : Disposable gloves

Eye protection : No special precautions required.

Skin and body protection : No special precautions required.

Protective measures : The type of protective equipment must be selected based on

the Environmental Health and Safety risk assessment. Consult a Environmental Health and Safety expert if necessary.

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance : liquid

Colour : clear

Odour : odourless

pH : 7.5

Solubility(ies)

Water solubility : soluble

10. STABILITY AND REACTIVITY

Reactivity : None reasonably foreseeable.

Chemical stability : Stable under recommended storage conditions.

Possibility of hazardous reac-

tions

: No dangerous reaction known under conditions of normal use.

Conditions to avoid : To avoid thermal decomposition, do not overheat.

Incompatible materials : Strong acids and strong bases

Reducing agents Strong oxidizing agents

Hazardous decomposition

products

: None known.

Permeabilization reagent



Version 1.12

Revision Date: 2018/01/09

SDS Number: 100000010887

Date of last issue: 2016/10/27 Date of first issue: 2015/09/21

11. TOXICOLOGICAL INFORMATION

Acute toxicity

Product:

Acute oral toxicity : Acute toxicity estimate: > 2,000 mg/kg

Method: Calculation method

Acute dermal toxicity : Acute toxicity estimate: > 2,000 mg/kg

Method: Calculation method

Components:

sodium azide:

Acute oral toxicity : LD50 (Rat): 27 mg/kg

Skin corrosion/irritation

No data available

Serious eye damage/eye irritation

No data available

Respiratory or skin sensitisation

No data available

Germ cell mutagenicity

No data available

Carcinogenicity

No data available

Reproductive toxicity

No data available

STOT - single exposure

No data available

STOT - repeated exposure

No data available

Repeated dose toxicity

No data available

Aspiration toxicity

No data available

12. ECOLOGICAL INFORMATION

Ecotoxicity

Components:

sodium azide:

Toxicity to fish : LC50 (Lepomis macrochirus (Bluegill sunfish)): 0.7 mg/l

Exposure time: 96 h

Permeabilization reagent



 Version
 Revision Date:
 SDS Number:
 Date of last issue: 2016/10/27

 1.12
 2018/01/09
 100000010887
 Date of first issue: 2015/09/21

Toxicity to daphnia and other

aquatic invertebrates

: EC50 (Daphnia pulex (Water flea)): 4.2 mg/l

Exposure time: 96 h

Toxicity to algae : IC50: 272 mg/l

Toxicity to bacteria : EC50 (Photobacterium phosphoreum): 38.5 mg/l

Persistence and degradability

No data available

Bioaccumulative potential

No data available

Mobility in soil

No data available

Hazardous to the ozone layer

Not applicable

Other adverse effects

No data available

13. DISPOSAL CONSIDERATIONS

Disposal methods

Waste from residues : In accordance with National, Federal, State and Local regula-

tions.

Contaminated packaging : Empty containers should be taken to an approved waste han-

dling site for recycling or disposal.

14. TRANSPORT INFORMATION

International Regulations

UNRTDG

Not regulated as a dangerous good

IATA-DGR

Not regulated as a dangerous good

IMDG-Code

Not regulated as a dangerous good

Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable for product as supplied.

National Regulations

Refer to section 15 for specific national regulation.

Permeabilization reagent



Version 1.12

Revision Date: 2018/01/09

SDS Number: 100000010887

Date of last issue: 2016/10/27 Date of first issue: 2015/09/21

15. REGULATORY INFORMATION

Related Regulations

Fire Service Law

Not applicable to dangerous materials / designated flammables.

Chemical Substance Control Law

Not applicable for Specified Chemical Substance, Monitoring Chemical Substance and Priority Assessment Chemical Substance.

Industrial Safety and Health Law

Harmful Substances Prohibited from Manufacture

Not applicable

Harmful Substances Required Permission for Manufacture

Not applicable

Substances Prevented From Impairment of Health

Not applicable

Circular concerning Information on Chemicals having Mutagenicity - Annex 2: Information on Existing Chemicals having Mutagenicity

Not applicable

Circular concerning Information on Chemicals having Mutagenicity - Annex 1: Information on Notified Substances having Mutagenicity

Not applicable

Substances Subject to be Notified Names

Not applicable

Substances Subject to be Indicated Names

Not applicable

Ordinance on Prevention of Hazards Due to Specified Chemical Substances

Not applicable

Ordinance on Prevention of Lead Poisoning

Not applicable

Ordinance on Prevention of Tetraalkyl Lead Poisoning

Not applicable

Ordinance on Prevention of Organic Solvent Poisoning

Not applicable

Enforcement Order of the Industrial Safety and Health Law - Attached table 1 (Dangerous Substances)

Not applicable

Poisonous and Deleterious Substances Control Law

Not applicable

Act on Confirmation, etc. of Release Amounts of Specific Chemical Substances in the Environment and Promotion of Improvements to the Management Thereof

Not applicable

Explosive Control Law

Permeabilization reagent



Version Revision Date: SDS Number: Date of last issue: 2016/10/27 1.12 2018/01/09 100000010887 Date of first issue: 2015/09/21

Not applicable

Vessel Safety Law

Not regulated as a dangerous good

Aviation Law

Not regulated as a dangerous good

Marine Pollution and Sea Disaster Prevention etc Law

Bulk transportation : Not classified as noxious liquid substance

Pack transportation : Not classified as marine pollutant

Waste Disposal and Public Cleansing Law

Industrial waste

16. OTHER INFORMATION

Full text of other abbreviations

AICS - Australian Inventory of Chemical Substances: ANTT - National Agency for Transport by Land of Brazil; ASTM - American Society for the Testing of Materials; bw - Body weight; CMR -Carcinogen, Mutagen or Reproductive Toxicant; CPR - Controlled Products Regulations; DIN -Standard of the German Institute for Standardisation; DSL - Domestic Substances List (Canada); ECx - Concentration associated with x% response; ELx - Loading rate associated with x% response; EmS - Emergency Schedule; ENCS - Existing and New Chemical Substances (Japan); ErCx - Concentration associated with x% growth rate response; ERG - Emergency Response Guide; GHS - Globally Harmonized System; GLP - Good Laboratory Practice; IARC - International Agency for Research on Cancer; IATA - International Air Transport Association; IBC - International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk; IC50 - Half maximal inhibitory concentration; ICAO - International Civil Aviation Organization; IECSC - Inventory of Existing Chemical Substances in China; IMDG - International Maritime Dangerous Goods; IMO - International Maritime Organization; ISHL - Industrial Safety and Health Law (Japan); ISO - International Organisation for Standardization; KECI - Korea Existing Chemicals Inventory; LC50 - Lethal Concentration to 50 % of a test population; LD50 - Lethal Dose to 50% of a test population (Median Lethal Dose); MARPOL - International Convention for the Prevention of Pollution from Ships; n.o.s. - Not Otherwise Specified; Nch - Chilean Norm; NO(A)EC -No Observed (Adverse) Effect Concentration; NO(A)EL - No Observed (Adverse) Effect Level; NOELR - No Observable Effect Loading Rate; NOM - Official Mexican Norm; NTP - National Toxicology Program; NZIoC - New Zealand Inventory of Chemicals; OECD - Organization for Economic Co-operation and Development; OPPTS - Office of Chemical Safety and Pollution Prevention: PBT - Persistent, Bioaccumulative and Toxic substance: PICCS - Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR - (Quantitative) Structure Activity Relationship; REACH - Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; SADT - Self-Accelerating Decomposition Temperature; SDS - Safety Data Sheet; TCSI - Taiwan Chemical Substance Inventory; TDG - Transportation of Dangerous Goods; TSCA - Toxic Substances Control Act (United States); UN - United Nations; UNRTDG - United Nations Recommendations on the Transport of Dangerous Goods; vPvB - Very Persistent and Very Bioaccumulative; WHMIS -Workplace Hazardous Materials Information System

Date format : yyyy/mm/dd Numbers 123,456.78

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not

Permeabilization reagent



 Version
 Revision Date:
 SDS Number:
 Date of last issue: 2016/10/27

 1.12
 2018/01/09
 100000010887
 Date of first issue: 2015/09/21

to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

JP / EN

Anti-EpCAM ferrofluid



Version Revision Date: SDS Number: Date of last issue: 2016/11/15
1.11 2018/01/09 100000010880 Date of first issue: 2015/09/21

1. PRODUCT AND COMPANY IDENTIFICATION

Product name : Anti-EpCAM ferrofluid

Substance name : Anti-EpCAM ferrofluid

7036

Manufacturer or supplier's details

Company name of supplier : Menarini Silicon Biosystems, Inc.

Address : 3401 Masons Mill Rd #100

Huntingdon Valley, PA

Us-info@siliconbiosystems.com

19006, USA

Telephone : 1 (800) 381-4929

Responsible/issuing person

E-mail address

Emergency telephone : US: (303)-389-1805

number International: +1 (303)-389-1805

Recommended use of the chemical and restrictions on use

Recommended use : Assay reagent

2. HAZARDS IDENTIFICATION

GHS Classification

Not a dangerous substance or mixture according to the Globally Harmonised System (GHS).

GHS label elements

Not a dangerous substance or mixture according to the Globally Harmonised System (GHS).

Other hazards which do not result in classification

None known.

3. COMPOSITION/INFORMATION ON INGREDIENTS

Substance / Mixture : Mixture

Chemical nature : Liquid

Hazardous components

| Chemical name | CAS-No. | Concentration (% w/w) | ENCS/ISHL number |
|---|--------------|-----------------------|------------------|
| Anti-EpCAM mouse mAb conjugated to Ferrofluid | Not Assigned | < 0.1 | |

For explanation of abbreviations see section 16.

Anti-EpCAM ferrofluid



Date of last issue: 2016/11/15 Version Revision Date: SDS Number: 100000010880 Date of first issue: 2015/09/21 1.11 2018/01/09

4. FIRST AID MEASURES

If inhaled If breathed in, move person into fresh air.

Consult a physician.

In case of skin contact Take off contaminated clothing and shoes immediately.

> Wash off immediately with plenty of water. If symptoms persist, call a physician. Wash contaminated clothing before re-use.

In case of eye contact : Rinse immediately with plenty of water, also under the eyelids,

> for at least 5 minutes. Remove contact lenses.

If eye irritation persists, consult a specialist.

If swallowed : If swallowed, rinse mouth with water (only if the person is con-

scious).

Call a physician immediately.

Most important symptoms and effects, both acute and

delayed

: No information available.

Notes to physician : Treat symptomatically.

5. FIREFIGHTING MEASURES

Suitable extinguishing media : Use extinguishing measures that are appropriate to local cir-

cumstances and the surrounding environment.

Specific hazards during fire-

fighting

: No information available.

Specific extinguishing meth-

ods

: No information available.

Special protective equipment

for firefighters

: In the event of fire, wear self-contained breathing apparatus.

6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment and emer-

gency procedures

: In the event of an accidental release the emergency response team must respond based on a risk assessment and use personal protective equipment as appropriate.

: Should not be released into the environment. **Environmental precautions**

Methods and materials for containment and cleaning up Large spills: Dam up. Soak up with inert absorbent material.

Keep in properly labelled containers.

Small spills: Gently cover the spill with an absorbent towel or

Large spills + Small spills: Keep in suitable, closed containers

Anti-EpCAM ferrofluid



Version 1.11

Revision Date: 2018/01/09

SDS Number: 100000010880

Date of last issue: 2016/11/15 Date of first issue: 2015/09/21

for disposal. Treat recovered material as described in the section "Disposal considerations".

7. HANDLING AND STORAGE

Handling

Advice on protection against

fire and explosion

: No data available

Advice on safe handling : To avoid thermal decomposition, do not overheat.

Avoid inhalation, ingestion and contact with skin and eyes.

Use personal protective equipment as required.

Avoidance of contact : None known.

Hygiene measures : Handle in accordance with good industrial hygiene and safety

practice.

Storage

Conditions for safe storage : To maintain product quality, do not store in heat or direct sun-

light.

Store in original container.

Keep containers tightly closed in a dry, cool and well-

ventilated place.

Keep away from heat and sources of ignition.

Keep locked up. Keep refrigerated.

Recommended storage tem-

perature

: 2-8°C

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Threshold limit value and permissible exposure limits for each component in the work environment

| Components | CAS-No. | Value type | Control parame- | Basis |
|---------------------------|--|------------|--------------------|-----------|
| | | (Form of | ters / Permissible | |
| | | exposure) | concentration | |
| Anti-EpCAM mouse mAb con- | Not Assigned | PBOEL-HHC | 2 | J&J |
| jugated to Ferrofluid | _ | | | OEL/PBOEL |
| | | | | HHC |
| | Further information: J&J has a hazard banding notation: PBOEL HHC. This substance is classified by J&J as being PBOEL HHC 2. This means that the OEL is estimated to be from 20 to 100 µg/m3 | | | |
| | | | | |
| | | | | |

Engineering measures

: All personal protective equipment should be based on a risk assessment. Consult a Environment Health Safety expert if

necessary.

Personal protective equipment

Anti-EpCAM ferrofluid



Version Revision Date: SDS Number: Date of last issue: 2016/11/15
1.11 2018/01/09 100000010880 Date of first issue: 2015/09/21

Respiratory protection : Engineering controls should always be the primary method of

controlling exposures.

If respiratory protective equipment is needed for certain activities, the type as well as the corresponding protection factor will depend upon the risk assessment and air concentrations, hazards, physical and warning properties of substances pre-

sent.

No personal respiratory protective equipment normally re-

quired.

Hand protection

Remarks : Disposable gloves

Eye protection : No special precautions required.

Skin and body protection : No special precautions required.

Protective measures : The type of protective equipment must be selected based on

the Environmental Health and Safety risk assessment. Consult a Environmental Health and Safety expert if necessary.

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance : liquid

Colour : brown

pH : 7.5

Solubility(ies)

Water solubility : soluble

10. STABILITY AND REACTIVITY

Reactivity : None reasonably foreseeable.

Chemical stability : Stable under recommended storage conditions.

Possibility of hazardous reac-

tions

: No dangerous reaction known under conditions of normal use.

Conditions to avoid : To avoid thermal decomposition, do not overheat.

Incompatible materials : None known.

Hazardous decomposition

products

: None known.

11. TOXICOLOGICAL INFORMATION

Acute toxicity

Anti-EpCAM ferrofluid



Version 1.11

Revision Date: 2018/01/09

SDS Number: 100000010880

Date of last issue: 2016/11/15 Date of first issue: 2015/09/21

No data available

Skin corrosion/irritation

No data available

Serious eye damage/eye irritation

No data available

Respiratory or skin sensitisation

No data available

Germ cell mutagenicity

No data available

Carcinogenicity

No data available

Reproductive toxicity

No data available

STOT - single exposure

No data available

STOT - repeated exposure

No data available

Repeated dose toxicity

No data available

Aspiration toxicity

No data available

12. ECOLOGICAL INFORMATION

Ecotoxicity

No data available

Persistence and degradability

No data available

Bioaccumulative potential

No data available

Mobility in soil

No data available

Hazardous to the ozone layer

Not applicable

Other adverse effects

No data available

13. DISPOSAL CONSIDERATIONS

Disposal methods

Anti-EpCAM ferrofluid



 Version
 Revision Date:
 SDS Number:
 Date of last issue: 2016/11/15

 1.11
 2018/01/09
 100000010880
 Date of first issue: 2015/09/21

Waste from residues : In accordance with National, Federal, State and Local regula-

tions.

Contaminated packaging : Empty containers should be taken to an approved waste han-

dling site for recycling or disposal.

14. TRANSPORT INFORMATION

International Regulations

UNRTDG

Not regulated as a dangerous good

IATA-DGR

Not regulated as a dangerous good

IMDG-Code

Not regulated as a dangerous good

Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable for product as supplied.

National Regulations

Refer to section 15 for specific national regulation.

15. REGULATORY INFORMATION

Related Regulations

Fire Service Law

Not applicable

Chemical Substance Control Law

Not applicable for Specified Chemical Substance, Monitoring Chemical Substance and Priority Assessment Chemical Substance.

Industrial Safety and Health Law

Harmful Substances Prohibited from Manufacture

Not applicable

Harmful Substances Required Permission for Manufacture

Not applicable

Substances Prevented From Impairment of Health

Not applicable

Circular concerning Information on Chemicals having Mutagenicity - Annex 2: Information on Existing Chemicals having Mutagenicity

Not applicable

Circular concerning Information on Chemicals having Mutagenicity - Annex 1: Information on Notified Substances having Mutagenicity

Not applicable

Anti-EpCAM ferrofluid



Version Revision Date: SDS Number: Date of last issue: 2016/11/15
1.11 2018/01/09 100000010880 Date of first issue: 2015/09/21

Substances Subject to be Notified Names

Not applicable

Substances Subject to be Indicated Names

Not applicable

Ordinance on Prevention of Hazards Due to Specified Chemical Substances

Not applicable

Ordinance on Prevention of Lead Poisoning

Not applicable

Ordinance on Prevention of Tetraalkyl Lead Poisoning

Not applicable

Ordinance on Prevention of Organic Solvent Poisoning

Not applicable

Enforcement Order of the Industrial Safety and Health Law - Attached table 1 (Dangerous Substances)

Not applicable

Poisonous and Deleterious Substances Control Law

Not applicable

Act on Confirmation, etc. of Release Amounts of Specific Chemical Substances in the Environment and Promotion of Improvements to the Management Thereof

Not applicable

Explosive Control Law

Not applicable

Vessel Safety Law

Not regulated as a dangerous good

Aviation Law

Not regulated as a dangerous good

Marine Pollution and Sea Disaster Prevention etc Law

Bulk transportation : Not classified as noxious liquid substance

Pack transportation : Not classified as marine pollutant

Waste Disposal and Public Cleansing Law

Industrial waste

16. OTHER INFORMATION

Full text of other abbreviations

AICS - Australian Inventory of Chemical Substances; ANTT - National Agency for Transport by Land of Brazil; ASTM - American Society for the Testing of Materials; bw - Body weight; CMR - Carcinogen, Mutagen or Reproductive Toxicant; CPR - Controlled Products Regulations; DIN - Standard of the German Institute for Standardisation; DSL - Domestic Substances List (Canada); ECx - Concentration associated with x% response; ELx - Loading rate associated with x% response; EmS - Emergency Schedule; ENCS - Existing and New Chemical Substances (Japan); ErCx - Concentration associated with x% growth rate response; ERG - Emergency Response Guide; GHS - Globally Harmonized System; GLP - Good Laboratory Practice; IARC - International Agency for Research on Cancer; IATA - International Air Transport Association; IBC - International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk;

SAFETY DATA SHEET Anti-EpCAM ferrofluid



Version Revision Date: SDS Number: Date of last issue: 2016/11/15
1.11 2018/01/09 100000010880 Date of first issue: 2015/09/21

IC50 - Half maximal inhibitory concentration; ICAO - International Civil Aviation Organization; IECSC - Inventory of Existing Chemical Substances in China; IMDG - International Maritime Dangerous Goods; IMO - International Maritime Organization; ISHL - Industrial Safety and Health Law (Japan); ISO - International Organisation for Standardization; KECI - Korea Existing Chemicals Inventory; LC50 - Lethal Concentration to 50 % of a test population; LD50 - Lethal Dose to 50% of a test population (Median Lethal Dose); MARPOL - International Convention for the Prevention of Pollution from Ships; n.o.s. - Not Otherwise Specified; Nch - Chilean Norm; NO(A)EC -No Observed (Adverse) Effect Concentration; NO(A)EL - No Observed (Adverse) Effect Level; NOELR - No Observable Effect Loading Rate; NOM - Official Mexican Norm; NTP - National Toxicology Program; NZIoC - New Zealand Inventory of Chemicals; OECD - Organization for Economic Co-operation and Development: OPPTS - Office of Chemical Safety and Pollution Prevention: PBT - Persistent, Bioaccumulative and Toxic substance: PICCS - Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR - (Quantitative) Structure Activity Relationship; REACH - Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; SADT - Self-Accelerating Decomposition Temperature; SDS - Safety Data Sheet; TCSI - Taiwan Chemical Substance Inventory; TDG - Transportation of Dangerous Goods; TSCA - Toxic Substances Control Act (United States); UN - United Nations; UNRTDG - United Nations Recommendations on the Transport of Dangerous Goods; vPvB - Very Persistent and Very Bioaccumulative; WHMIS -Workplace Hazardous Materials Information System

Date format : yyyy/mm/dd Numbers 123,456.78

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

JP / EN

Dilution buffer



Version Revision Date: SDS Number: Date of last issue: 2016/10/27 1.12 2018/01/09 100000010879 Date of first issue: 2015/09/22

1. PRODUCT AND COMPANY IDENTIFICATION

Product name : Dilution buffer

Substance name : Dilution buffer

7039

Manufacturer or supplier's details

Company name of supplier : Menarini Silicon Biosystems, Inc.

Address : 3401 Masons Mill Rd #100

Huntingdon Valley, PA 19006,

Us-info@siliconbiosystems.com

USA

Telephone : 1 (800) 381-4929

Responsible/issuing person

E-mail address

Emergency telephone : US : (303)-389-1805

number International: +1 (303)-389-1805

Recommended use of the chemical and restrictions on use

Recommended use : Assay reagent

2. HAZARDS IDENTIFICATION

GHS Classification

Not a dangerous substance or mixture according to the Globally Harmonised System (GHS).

GHS label elements

Not a dangerous substance or mixture according to the Globally Harmonised System (GHS).

Other hazards which do not result in classification

None known.

3. COMPOSITION/INFORMATION ON INGREDIENTS

Substance / Mixture : Mixture

Chemical nature : Liquid

Hazardous components

| Chemical name | CAS-No. | Concentration (% | ENCS/ISHL |
|---------------|------------|------------------|-----------|
| | | w/w) | number |
| sodium azide | 26628-22-8 | >= 0.1 - < 1 | 1-482 |

For explanation of abbreviations see section 16.

4. FIRST AID MEASURES

Dilution buffer



Version Revision Date: SDS Number: Date of last issue: 2016/10/27 1.12 2018/01/09 100000010879 Date of first issue: 2015/09/22

If inhaled : If breathed in, move person into fresh air.

Consult a physician.

In case of skin contact : Take off contaminated clothing and shoes immediately.

Wash off with plenty of water.

If symptoms persist, call a physician.

In case of eye contact : Rinse immediately with plenty of water, also under the eyelids,

for at least 5 minutes. Remove contact lenses.

If eye irritation persists, consult a specialist.

If swallowed, rinse mouth with water (only if the person is con-

scious).

Call a physician immediately.

Most important symptoms and effects, both acute and

delayed

: No information available.

Notes to physician : Treat symptomatically.

5. FIREFIGHTING MEASURES

Suitable extinguishing media : Use extinguishing measures that are appropriate to local cir-

cumstances and the surrounding environment.

Specific hazards during fire-

fighting

: No information available.

Specific extinguishing meth-

ods

: No information available.

Special protective equipment

for firefighters

: In the event of fire, wear self-contained breathing apparatus.

6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment and emer-

gency procedures

: In the event of an accidental release the emergency response team must respond based on a risk assessment and use personal protective equipment as appropriate.

: Should not be released into the environment.

Environmental precautions : Should not be released into the

Methods and materials for containment and cleaning up

: Large spills: Dam up. Soak up with inert absorbent material.

Keep in properly labelled containers.

Small spills: Gently cover the spill with an absorbent towel or

nad

Large spills + Small spills: Keep in suitable, closed containers for disposal. Treat recovered material as described in the sec-

tion "Disposal considerations".

Dilution buffer



Version 1.12

Revision Date: 2018/01/09

SDS Number: 100000010879

Date of last issue: 2016/10/27 Date of first issue: 2015/09/22

7. HANDLING AND STORAGE

Handling

Advice on protection against

fire and explosion

: No data available

Advice on safe handling : To avoid thermal decomposition, do not overheat.

Avoid inhalation, ingestion and contact with skin and eyes.

Use personal protective equipment as required.

Avoidance of contact : Oxidizing agents

Hygiene measures : Handle in accordance with good industrial hygiene and safety

practice.

Storage

Conditions for safe storage : To maintain product quality, do not store in heat or direct sun-

light.

Store in original container.

Keep containers tightly closed in a dry, cool and well-

ventilated place.

Keep away from heat and sources of ignition.

Keep locked up. Keep refrigerated.

Recommended storage tem-

perature

: 2-8°C

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Threshold limit value and permissible exposure limits for each component in the work environment

| Components | CAS-No. | Value type (Form of exposure) | Control parameters / Permissible concentration | Basis |
|--------------|------------|-------------------------------------|--|-------|
| sodium azide | 26628-22-8 | C (Vapour) | 0.11 ppm (Hydrazoic acid) | ACGIH |
| | | С | 0.29 mg/m3 (Sodium azide) | ACGIH |

Engineering measures

: All personal protective equipment should be based on a risk assessment. Consult a Environment Health Safety expert if necessary.

Personal protective equipment

Respiratory protection : Engineering controls should always be the primary method of

controlling exposures.

If respiratory protective equipment is needed for certain activities, the type as well as the corresponding protection factor will depend upon the risk assessment and air concentrations, hazards, physical and warning properties of substances pre-

Dilution buffer



Version Revision Date: SDS Number: Date of last issue: 2016/10/27 1.12 2018/01/09 100000010879 Date of first issue: 2015/09/22

sent.

No personal respiratory protective equipment normally re-

quired.

Hand protection

Remarks : Disposable gloves

Eye protection : No special precautions required.

Skin and body protection : No special precautions required.

Protective measures : The type of protective equipment must be selected based on

the Environmental Health and Safety risk assessment. Consult a Environmental Health and Safety expert if necessary.

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance : liquid

Colour : clear

Odour : odourless

pH : 7.5

Solubility(ies)

Water solubility : soluble

10. STABILITY AND REACTIVITY

Reactivity : None reasonably foreseeable.

Chemical stability : Stable under recommended storage conditions.

Possibility of hazardous reac-

tions

: No dangerous reaction known under conditions of normal use.

Conditions to avoid : To avoid thermal decomposition, do not overheat.

Incompatible materials : Oxidizing agents

Hazardous decomposition

products

: None known.

11. TOXICOLOGICAL INFORMATION

Acute toxicity

Product:

Acute oral toxicity : Acute toxicity estimate: > 2,000 mg/kg

Method: Calculation method

Dilution buffer



Version Revision Date: SDS Number: Date of last issue: 2016/10/27 1.12 2018/01/09 100000010879 Date of first issue: 2015/09/22

Acute dermal toxicity : Acute toxicity estimate: > 2,000 mg/kg

Method: Calculation method

Components:

sodium azide:

Acute oral toxicity : LD50 (Rat): 27 mg/kg

Skin corrosion/irritation

No data available

Serious eye damage/eye irritation

No data available

Respiratory or skin sensitisation

No data available

Germ cell mutagenicity

No data available

Carcinogenicity

No data available

Reproductive toxicity

No data available

STOT - single exposure

No data available

STOT - repeated exposure

No data available

Repeated dose toxicity

No data available

Aspiration toxicity

No data available

12. ECOLOGICAL INFORMATION

Ecotoxicity

Components:

sodium azide:

Toxicity to fish : LC50 (Lepomis macrochirus (Bluegill sunfish)): 0.7 mg/l

Exposure time: 96 h

Toxicity to daphnia and other

aquatic invertebrates

: EC50 (Daphnia pulex (Water flea)): 4.2 mg/l

Exposure time: 96 h

Toxicity to algae : IC50: 272 mg/l

Toxicity to bacteria : EC50 (Photobacterium phosphoreum): 38.5 mg/l

Dilution buffer



Version 1.12

Revision Date: 2018/01/09

SDS Number: 100000010879

Date of last issue: 2016/10/27 Date of first issue: 2015/09/22

Persistence and degradability

No data available

Bioaccumulative potential

No data available

Mobility in soil

No data available

Hazardous to the ozone layer

Not applicable

Other adverse effects

No data available

13. DISPOSAL CONSIDERATIONS

Disposal methods

Waste from residues : In accordance with National, Federal, State and Local regula-

tions.

Contaminated packaging : Empty containers should be taken to an approved waste han-

dling site for recycling or disposal.

14. TRANSPORT INFORMATION

International Regulations

UNRTDG

Not regulated as a dangerous good

IATA-DGR

Not regulated as a dangerous good

IMDG-Code

Not regulated as a dangerous good

Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable for product as supplied.

National Regulations

Refer to section 15 for specific national regulation.

15. REGULATORY INFORMATION

Related Regulations

Fire Service Law

Not applicable to dangerous materials / designated flammables.

Dilution buffer



Version 1.12

Revision Date: 2018/01/09

SDS Number: 100000010879

Date of last issue: 2016/10/27 Date of first issue: 2015/09/22

Chemical Substance Control Law

Not applicable for Specified Chemical Substance, Monitoring Chemical Substance and Priority Assessment Chemical Substance.

Industrial Safety and Health Law

Harmful Substances Prohibited from Manufacture

Not applicable

Harmful Substances Required Permission for Manufacture

Not applicable

Substances Prevented From Impairment of Health

Not applicable

Circular concerning Information on Chemicals having Mutagenicity - Annex 2: Information on Existing Chemicals having Mutagenicity

Not applicable

Circular concerning Information on Chemicals having Mutagenicity - Annex 1: Information on Notified Substances having Mutagenicity

Not applicable

Substances Subject to be Notified Names

Not applicable

Substances Subject to be Indicated Names

Not applicable

Ordinance on Prevention of Hazards Due to Specified Chemical Substances

Not applicable

Ordinance on Prevention of Lead Poisoning

Not applicable

Ordinance on Prevention of Tetraalkyl Lead Poisoning

Not applicable

Ordinance on Prevention of Organic Solvent Poisoning

Not applicable

Enforcement Order of the Industrial Safety and Health Law - Attached table 1 (Dangerous Substances)

Not applicable

Poisonous and Deleterious Substances Control Law

Not applicable

Act on Confirmation, etc. of Release Amounts of Specific Chemical Substances in the Environment and Promotion of Improvements to the Management Thereof

Not applicable

High Pressure Gas Safety Act

Not applicable

Explosive Control Law

Not applicable

Vessel Safety Law

Not regulated as a dangerous good

Aviation Law

Not regulated as a dangerous good

Dilution buffer



Version Revision Date: SDS Number: Date of last issue: 2016/10/27 1.12 2018/01/09 100000010879 Date of first issue: 2015/09/22

Marine Pollution and Sea Disaster Prevention etc Law

Bulk transportation : Not classified as noxious liquid substance

Pack transportation : Not classified as marine pollutant

Waste Disposal and Public Cleansing Law

Industrial waste

16. OTHER INFORMATION

Full text of other abbreviations

AICS - Australian Inventory of Chemical Substances: ANTT - National Agency for Transport by Land of Brazil; ASTM - American Society for the Testing of Materials; bw - Body weight; CMR -Carcinogen, Mutagen or Reproductive Toxicant; CPR - Controlled Products Regulations; DIN -Standard of the German Institute for Standardisation; DSL - Domestic Substances List (Canada); ECx - Concentration associated with x% response; ELx - Loading rate associated with x% response; EmS - Emergency Schedule; ENCS - Existing and New Chemical Substances (Japan); ErCx - Concentration associated with x% growth rate response; ERG - Emergency Response Guide; GHS - Globally Harmonized System; GLP - Good Laboratory Practice; IARC - International Agency for Research on Cancer; IATA - International Air Transport Association; IBC - International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk; IC50 - Half maximal inhibitory concentration; ICAO - International Civil Aviation Organization; IECSC - Inventory of Existing Chemical Substances in China; IMDG - International Maritime Dangerous Goods; IMO - International Maritime Organization; ISHL - Industrial Safety and Health Law (Japan); ISO - International Organisation for Standardization; KECI - Korea Existing Chemicals Inventory; LC50 - Lethal Concentration to 50 % of a test population; LD50 - Lethal Dose to 50% of a test population (Median Lethal Dose); MARPOL - International Convention for the Prevention of Pollution from Ships; n.o.s. - Not Otherwise Specified; Nch - Chilean Norm; NO(A)EC -No Observed (Adverse) Effect Concentration; NO(A)EL - No Observed (Adverse) Effect Level; NOELR - No Observable Effect Loading Rate; NOM - Official Mexican Norm; NTP - National Toxicology Program; NZIoC - New Zealand Inventory of Chemicals; OECD - Organization for Economic Co-operation and Development; OPPTS - Office of Chemical Safety and Pollution Prevention: PBT - Persistent, Bioaccumulative and Toxic substance: PICCS - Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR - (Quantitative) Structure Activity Relationship; REACH - Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; SADT - Self-Accelerating Decomposition Temperature; SDS - Safety Data Sheet; TCSI - Taiwan Chemical Substance Inventory; TDG - Transportation of Dangerous Goods; TSCA - Toxic Substances Control Act (United States); UN - United Nations; UNRTDG - United Nations Recommendations on the Transport of Dangerous Goods; vPvB - Very Persistent and Very Bioaccumulative; WHMIS -Workplace Hazardous Materials Information System

Date format : yyyy/mm/dd Numbers 123,456.78

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

JP / EN

Cell fixative



Version Revision Date: SDS Number: Date of last issue: 2016/11/15 1.14 2018/01/09 100000010702 Date of first issue: 2015/09/22

1. PRODUCT AND COMPANY IDENTIFICATION

Product name : Cell fixative

Substance name : Cell fixative

7042

Manufacturer or supplier's details

Company name of supplier : Menarini Silicon Biosystems, Inc.

Address : 3401 Masons Mill Rd #100

Huntingdon Valley, PA

19006, USA

Telephone : 1 (800) 381-4929

E-mail address

Responsible/issuing person

Us-info@siliconbiosystems.com

Emergency telephone

number

US: (303)-389-1805

International: +1 (303)-389-1805

Recommended use of the chemical and restrictions on use

Recommended use : Assay reagent

2. HAZARDS IDENTIFICATION

GHS Classification

Skin sensitisation : Category 1

GHS label elements

Hazard pictograms



Signal word : Warning

Hazard statements : H317 May cause an allergic skin reaction.

Precautionary statements : **Prevention:**

P261 Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.
P272 Contaminated work clothing should not be allowed out of

the workplace.

P280 Wear protective gloves.

Response:

P302 + P352 IF ON SKIN: Wash with plenty of water.

P333 + P313 If skin irritation or rash occurs: Get medical ad-

vice/ attention.

P362 + P364 Take off contaminated clothing and wash it before

reuse.

Cell fixative



Version Revision Date: SDS Number: Date of last issue: 2016/11/15
1.14 2018/01/09 100000010702 Date of first issue: 2015/09/22

Disposal:

P501 Dispose of contents/ container to an approved waste

disposal plant.

Other hazards which do not result in classification

None known.

3. COMPOSITION/INFORMATION ON INGREDIENTS

Substance / Mixture : Mixture

Chemical nature : Liquid

Hazardous components

| Chemical name | CAS-No. | Concentration (% w/w) | ENCS/ISHL number |
|---------------|------------|-----------------------|------------------|
| IMIDUREA | 39236-46-9 | >= 1 - < 10 | |
| sodium azide | 26628-22-8 | >= 0.1 - < 1 | 1-482 |

For explanation of abbreviations see section 16.

4. FIRST AID MEASURES

If inhaled : If breathed in, move person into fresh air.

Consult a physician.

In case of skin contact : Take off contaminated clothing and shoes immediately.

Wash off with plenty of water.

If symptoms persist, call a physician.

In case of eye contact : Rinse immediately with plenty of water, also under the eyelids,

for at least 15 minutes. Remove contact lenses.

If eye irritation persists, consult a specialist.

If swallowed : If swallowed, rinse mouth with water (only if the person is con-

scious).

Call a physician immediately.

Most important symptoms and effects, both acute and

delayed

: No information available.

Notes to physician : Treat symptomatically.

5. FIREFIGHTING MEASURES

Suitable extinguishing media : Use extinguishing measures that are appropriate to local cir-

cumstances and the surrounding environment.

Specific hazards during fire-

fighting

: No information available.

Cell fixative



Version 1.14

Revision Date: 2018/01/09

SDS Number: 100000010702

Date of last issue: 2016/11/15 Date of first issue: 2015/09/22

Hazardous combustion prod-

ucts

: No hazardous combustion products are known

Specific extinguishing meth-

ods

: No information available.

Special protective equipment

for firefighters

: In the event of fire, wear self-contained breathing apparatus.

6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment and emergency procedures

: In the event of an accidental release the emergency response team must respond based on a risk assessment and use personal protective equipment as appropriate.

Environmental precautions : Should not be released into the environment.

Methods and materials for containment and cleaning up

: Large spills: Dam up. Soak up with inert absorbent material.

Keep in properly labelled containers.

Small spills: Gently cover the spill with an absorbent towel or

pad.

Large spills + Small spills: Keep in suitable, closed containers for disposal. Treat recovered material as described in the sec-

tion "Disposal considerations".

7. HANDLING AND STORAGE

Handling

Advice on protection against

fire and explosion

: No data available

Advice on safe handling : To avoid thermal decomposition, do not overheat.

Avoid inhalation, ingestion and contact with skin and eyes.

Use personal protective equipment as required.

Avoidance of contact : Strong acids and strong bases

Reducing agents Oxidizing agents

Hygiene measures :

Handle in acco

Storage

Conditions for safe storage : To maintain product quality, do not store in heat or direct sun-

light.

Store in original container.

Keep containers tightly closed in a dry, cool and well-

ventilated place.

Keep away from heat and sources of ignition.

Keep locked up.

Cell fixative



Version Revision Date: SDS Number: Date of last issue: 2016/11/15
1.14 2018/01/09 100000010702 Date of first issue: 2015/09/22

Keep refrigerated.

Recommended storage tem-

perature

: 2-8°C

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Threshold limit value and permissible exposure limits for each component in the work environment

| Components | CAS-No. | Value type | Control parame- | Basis |
|--------------|------------|------------|--------------------|-------|
| | | (Form of | ters / Permissible | |
| | | exposure) | concentration | |
| sodium azide | 26628-22-8 | C (Vapour) | 0.11 ppm | ACGIH |
| | | | (Hydrazoic acid) | |
| | | С | 0.29 mg/m3 | ACGIH |
| | | | (Sodium azide) | |

Engineering measures

: All personal protective equipment should be based on a risk assessment. Consult a Environment Health Safety expert if necessary.

Personal protective equipment

Respiratory protection

Engineering controls should always be the primary method of

controlling exposures.

If respiratory protective equipment is needed for certain activities, the type as well as the corresponding protection factor will depend upon the risk assessment and air concentrations, hazards, physical and warning properties of substances pre-

sent.

No personal respiratory protective equipment normally re-

quired.

Hand protection

Remarks : Disposable gloves

Eye protection : No special precautions required.

Skin and body protection : No special precautions required.

Protective measures : The type of protective equipment must be selected based on

the Environmental Health and Safety risk assessment. Consult a Environmental Health and Safety expert if necessary.

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance : liquid

Colour : clear

Odour : odourless

Cell fixative



Version Revision Date: SDS Number: Date of last issue: 2016/11/15 1.14 2018/01/09 100000010702 Date of first issue: 2015/09/22

pH : 7.5

Solubility(ies)

Water solubility : soluble

10. STABILITY AND REACTIVITY

Reactivity : None reasonably foreseeable.

Chemical stability : Stable under recommended storage conditions.

Possibility of hazardous reac-

tions

: No dangerous reaction known under conditions of normal use.

Conditions to avoid : To avoid thermal decomposition, do not overheat.

Incompatible materials : Strong acids and strong bases

Reducing agents Oxidizing agents

Hazardous decomposition

products

: None known.

11. TOXICOLOGICAL INFORMATION

Acute toxicity

Product:

Acute oral toxicity : Acute toxicity estimate: > 2,000 mg/kg

Method: Calculation method

Acute dermal toxicity : Acute toxicity estimate: > 2,000 mg/kg

Method: Calculation method

Components:

IMIDUREA:

Acute oral toxicity : LD50 (Rat): 11,300 mg/kg

Acute inhalation toxicity : LC50 (Rat): > 5.5 mg/l

Exposure time: 1 h

Acute dermal toxicity : LD50 (Rabbit): > 5,000 mg/kg

sodium azide:

Acute oral toxicity : LD50 (Rat): 27 mg/kg

Skin corrosion/irritation

Components:

IMIDUREA:

Result: No skin irritation

Cell fixative



Version Revision Date: SDS Number: Date of last issue: 2016/11/15
1.14 2018/01/09 100000010702 Date of first issue: 2015/09/22

Serious eye damage/eye irritation

Components:

IMIDUREA:

Result: No eye irritation

Respiratory or skin sensitisation

Components:

IMIDUREA:

Method: Maximisation Test

Result: May cause sensitisation by skin contact.

Method: Local Lymph Node Assay (LLNA) in mice Result: May cause sensitisation by skin contact.

Germ cell mutagenicity

Components:

IMIDUREA:

Genotoxicity in vitro : Test Type: Ames test

Metabolic activation: with and without metabolic activation

Result: negative

GLP: yes

: Test Type: Chromosome aberration test in vitro

Metabolic activation: with and without metabolic activation

Result: negative

Genotoxicity in vivo : Test Type: Micronucleus test

Species: Mouse

Application Route: Oral

Result: negative

Germ cell mutagenicity -

Assessment

: No information available.

Carcinogenicity

Components:

IMIDUREA:

Carcinogenicity - Assessment : No information available.

Reproductive toxicity

Components:

IMIDUREA:

Teratogenicity - Assessment : No information available.

STOT - single exposure

No data available

Cell fixative



Version Revision Date: SDS Number: Date of last issue: 2016/11/15 1.14 2018/01/09 100000010702 Date of first issue: 2015/09/22

STOT - repeated exposure

No data available

Repeated dose toxicity

Components:

IMIDUREA:

Species: Rat

NOAEL: 200 mg/kg LOAEL: 500 mg/kg Application Route: Oral

Species: Rabbit NOAEL: 200 mg/kg Application Route: Dermal

Aspiration toxicity

No data available

12. ECOLOGICAL INFORMATION

Ecotoxicity

Components:

IMIDUREA:

Toxicity to fish : Remarks: No data available

sodium azide:

Toxicity to fish : LC50 (Lepomis macrochirus (Bluegill sunfish)): 0.7 mg/l

Exposure time: 96 h

Toxicity to daphnia and other

aquatic invertebrates

: EC50 (Daphnia pulex (Water flea)): 4.2 mg/l

Exposure time: 96 h

Toxicity to algae : IC50: 272 mg/l

Toxicity to bacteria : EC50 (Photobacterium phosphoreum): 38.5 mg/l

Persistence and degradability

Components:

IMIDUREA:

Biodegradability : Remarks: No data available

Bioaccumulative potential

Components:

IMIDUREA:

Bioaccumulation : Remarks: No data available

Cell fixative



Version Revision Date: SDS Number: Date of last issue: 2016/11/15 1.14 2018/01/09 100000010702 Date of first issue: 2015/09/22

Mobility in soil

Components:

IMIDUREA:

Distribution among : Remarks: No data available

environmental compartments

Hazardous to the ozone layer

Not applicable

Other adverse effects

Components:

IMIDUREA:

Results of PBT and vPvB

assessment

Additional ecological

information

: No information available.

: No data available

13. DISPOSAL CONSIDERATIONS

Disposal methods

Waste from residues : In accordance with National, Federal, State and Local regula-

tions.

Contaminated packaging : Empty containers should be taken to an approved waste han-

dling site for recycling or disposal.

14. TRANSPORT INFORMATION

International Regulations

UNRTDG

Not regulated as a dangerous good

IATA-DGR

Not regulated as a dangerous good

IMDG-Code

Not regulated as a dangerous good

Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable for product as supplied.

National Regulations

Refer to section 15 for specific national regulation.

15. REGULATORY INFORMATION

Related Regulations

Cell fixative



Version Revision Date: SDS Number: Date of last issue: 2016/11/15 1.14 2018/01/09 100000010702 Date of first issue: 2015/09/22

Fire Service Law

Not applicable to dangerous materials / designated flammables.

Chemical Substance Control Law

Not applicable for Specified Chemical Substance, Monitoring Chemical Substance and Priority Assessment Chemical Substance.

Industrial Safety and Health Law

Harmful Substances Prohibited from Manufacture

Not applicable

Harmful Substances Required Permission for Manufacture

Not applicable

Substances Prevented From Impairment of Health

Not applicable

Circular concerning Information on Chemicals having Mutagenicity - Annex 2: Information on Existing Chemicals having Mutagenicity

Not applicable

Circular concerning Information on Chemicals having Mutagenicity - Annex 1: Information on Notified Substances having Mutagenicity

Not applicable

Substances Subject to be Notified Names

Not applicable

Substances Subject to be Indicated Names

Not applicable

Ordinance on Prevention of Hazards Due to Specified Chemical Substances

Not applicable

Ordinance on Prevention of Lead Poisoning

Not applicable

Ordinance on Prevention of Tetraalkyl Lead Poisoning

Not applicable

Ordinance on Prevention of Organic Solvent Poisoning

Not applicable

Enforcement Order of the Industrial Safety and Health Law - Attached table 1 (Dangerous Substances)

Not applicable

Poisonous and Deleterious Substances Control Law

Not applicable

Act on Confirmation, etc. of Release Amounts of Specific Chemical Substances in the Environment and Promotion of Improvements to the Management Thereof

Not applicable

High Pressure Gas Safety Act

Not applicable

Explosive Control Law

Not applicable

Vessel Safety Law

Not regulated as a dangerous good

Cell fixative



Version Revision Date: SDS Number: Date of last issue: 2016/11/15
1.14 2018/01/09 100000010702 Date of first issue: 2015/09/22

Aviation Law

Not regulated as a dangerous good

Marine Pollution and Sea Disaster Prevention etc Law

Bulk transportation : Not classified as noxious liquid substance

Pack transportation : Not classified as marine pollutant

Waste Disposal and Public Cleansing Law

Industrial waste

16. OTHER INFORMATION

Full text of other abbreviations

AICS - Australian Inventory of Chemical Substances: ANTT - National Agency for Transport by Land of Brazil; ASTM - American Society for the Testing of Materials; bw - Body weight; CMR -Carcinogen, Mutagen or Reproductive Toxicant; CPR - Controlled Products Regulations; DIN -Standard of the German Institute for Standardisation; DSL - Domestic Substances List (Canada); ECx - Concentration associated with x% response; ELx - Loading rate associated with x% response; EmS - Emergency Schedule; ENCS - Existing and New Chemical Substances (Japan); ErCx - Concentration associated with x% growth rate response; ERG - Emergency Response Guide; GHS - Globally Harmonized System; GLP - Good Laboratory Practice; IARC - International Agency for Research on Cancer; IATA - International Air Transport Association; IBC - International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk; IC50 - Half maximal inhibitory concentration; ICAO - International Civil Aviation Organization; IECSC - Inventory of Existing Chemical Substances in China; IMDG - International Maritime Dangerous Goods; IMO - International Maritime Organization; ISHL - Industrial Safety and Health Law (Japan); ISO - International Organisation for Standardization; KECI - Korea Existing Chemicals Inventory; LC50 - Lethal Concentration to 50 % of a test population; LD50 - Lethal Dose to 50% of a test population (Median Lethal Dose); MARPOL - International Convention for the Prevention of Pollution from Ships; n.o.s. - Not Otherwise Specified; Nch - Chilean Norm; NO(A)EC -No Observed (Adverse) Effect Concentration; NO(A)EL - No Observed (Adverse) Effect Level; NOELR - No Observable Effect Loading Rate; NOM - Official Mexican Norm; NTP - National Toxicology Program; NZIoC - New Zealand Inventory of Chemicals; OECD - Organization for Economic Co-operation and Development; OPPTS - Office of Chemical Safety and Pollution Prevention; PBT - Persistent, Bioaccumulative and Toxic substance; PICCS - Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR - (Quantitative) Structure Activity Relationship; REACH - Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; SADT - Self-Accelerating Decomposition Temperature; SDS - Safety Data Sheet; TCSI - Taiwan Chemical Substance Inventory; TDG - Transportation of Dangerous Goods; TSCA - Toxic Substances Control Act (United States); UN - United Nations; UNRTDG - United Nations Recommendations on the Transport of Dangerous Goods; vPvB - Very Persistent and Very Bioaccumulative; WHMIS -Workplace Hazardous Materials Information System

Date format : yyyy/mm/dd Numbers 123,456.78

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

Cell fixative



 Version
 Revision Date:
 SDS Number:
 Date of last issue: 2016/11/15

 1.14
 2018/01/09
 100000010702
 Date of first issue: 2015/09/22

JP / EN