

## CellSave Preservative 100 Tubes

Version	Revision Date:	SDS Number:	Date of last issue: 2017-08-02
1.44	2018-01-09	100000011151	Date of first issue: 2015-02-26

### SECTION 1: Identification of the substance/mixture and of the company/undertaking

#### 1.1 Product identifier

Trade name : CellSave Preservative 100 Tubes  
Substance name : CellSave Preservative 100 Tubes

#### 1.2 Relevant identified uses of the substance or mixture and uses advised against

Use of the Sub- : Assay reagent  
stance/Mixture

#### 1.3 Details of the supplier of the safety data sheet

Company : Menarini Silicon Biosystems, Inc.  
3401 Masons Mill Rd #100  
Huntingdon Valley, PA 19006  
USA

Telephone : 1 (800) 381-4929

Telefax :

E-mail address : Us-info@siliconbiosystems.com  
Responsible/issuing person

#### 1.4 Emergency telephone number

US: (303) 389-1805  
International:: +1 (303) 389-1805

### SECTION 2: Hazards identification

#### 2.1 Classification of the substance or mixture

##### Classification (REGULATION (EC) No 1272/2008)

Skin sensitisation, Category 1 H317: May cause an allergic skin reaction.

#### 2.2 Label elements

##### Labelling (REGULATION (EC) No 1272/2008)

Hazard pictograms :



Signal word : Warning

Hazard statements : H317 May cause an allergic skin reaction.

Precautionary statements : **Prevention:**

## CellSave Preservative 100 Tubes

Version 1.44	Revision Date: 2018-01-09	SDS Number: 100000011151	Date of last issue: 2017-08-02 Date of first issue: 2015-02-26
-----------------	------------------------------	-----------------------------	---

P261	Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves.
<b>Response:</b> P333 + P313	If skin irritation or rash occurs: Get medical advice/ attention.
P362 + P364	Take off contaminated clothing and wash it before reuse.
<b>Disposal:</b> P501	Dispose of contents/ container to an approved waste disposal plant.

Hazardous components which must be listed on the label:  
IMIDUREA

### 2.3 Other hazards

This substance/mixture contains no components considered to be either persistent, bioaccumulative and toxic (PBT), or very persistent and very bioaccumulative (vPvB) at levels of 0.1% or higher.

## SECTION 3: Composition/information on ingredients

### 3.2 Mixtures

Chemical nature : Liquid

#### Hazardous components

Chemical name	CAS-No. EC-No. Registration number	Classification (REGULATION (EC) No 1272/2008)	Concentration (%)
IMIDUREA	39236-46-9 254-372-6	Skin Sens. 1; H317	>= 25 - < 50

For explanation of abbreviations see section 16.

## SECTION 4: First aid measures

### 4.1 Description of first aid measures

- If inhaled : If breathed in, move person into fresh air.  
Consult a physician.
- In case of skin contact : Take off contaminated clothing and shoes immediately.  
Wash off immediately with plenty of water.  
If symptoms persist, call a physician.  
Wash contaminated clothing before re-use.

## CellSave Preservative 100 Tubes

Version	Revision Date:	SDS Number:	Date of last issue: 2017-08-02
1.44	2018-01-09	100000011151	Date of first issue: 2015-02-26

- In case of eye contact : Rinse immediately with plenty of water, also under the eyelids, for at least 5 minutes.  
Remove contact lenses.  
If eye irritation persists, consult a specialist.
- If swallowed : If swallowed, rinse mouth with water (only if the person is conscious).  
Call a physician immediately.

### 4.2 Most important symptoms and effects, both acute and delayed

Symptoms : No information available.

### 4.3 Indication of any immediate medical attention and special treatment needed

Treatment : Treat symptomatically.

## SECTION 5: Firefighting measures

### 5.1 Extinguishing media

Suitable extinguishing media : Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

### 5.2 Special hazards arising from the substance or mixture

Specific hazards during fire-fighting : No information available.

### 5.3 Advice for firefighters

Special protective equipment for firefighters : In the event of fire, wear self-contained breathing apparatus.

Further information : No information available.

## SECTION 6: Accidental release measures

### 6.1 Personal precautions, protective equipment and emergency procedures

Personal precautions : In the event of an accidental release the emergency response team must respond based on a risk assessment and use personal protective equipment as appropriate.

### 6.2 Environmental precautions

Environmental precautions : Should not be released into the environment.

### 6.3 Methods and material for containment and cleaning up

Methods for cleaning up : Large spills: Dam up. Soak up with inert absorbent material. Keep in properly labelled containers.

## CellSave Preservative 100 Tubes

Version	Revision Date:	SDS Number:	Date of last issue: 2017-08-02
1.44	2018-01-09	100000011151	Date of first issue: 2015-02-26

Small spills: Gently cover the spill with an absorbent towel or pad.  
Large spills + Small spills: Keep in suitable, closed containers for disposal. Treat recovered material as described in the section "Disposal considerations".

### 6.4 Reference to other sections

For disposal information, see section 13

## SECTION 7: Handling and storage

### 7.1 Precautions for safe handling

Advice on safe handling : To avoid thermal decomposition, do not overheat. Avoid inhalation, ingestion and contact with skin and eyes. Use personal protective equipment as required.

Advice on protection against fire and explosion : No data available

Hygiene measures : Handle in accordance with good industrial hygiene and safety practice.

### 7.2 Conditions for safe storage, including any incompatibilities

Requirements for storage areas and containers : To maintain product quality, do not store in heat or direct sunlight. Store in original container. Keep containers tightly closed in a dry, cool and well-ventilated place. Keep away from heat and sources of ignition. Keep locked up. Keep refrigerated.

Recommended storage temperature : 4 - 30 °C

### 7.3 Specific end use(s)

Specific use(s) : Consult the technical guidelines for the use of this substance/mixture.

## SECTION 8: Exposure controls/personal protection

### 8.1 Control parameters

Contains no substances with occupational exposure limit values.

### 8.2 Exposure controls

#### Engineering measures

All personal protective equipment should be based on a risk assessment. Consult a Environment Health Safety expert if necessary.

#### Personal protective equipment

Eye protection : No special precautions required.

## CellSave Preservative 100 Tubes

Version	Revision Date:	SDS Number:	Date of last issue: 2017-08-02
1.44	2018-01-09	100000011151	Date of first issue: 2015-02-26

Hand protection  
Remarks : Disposable gloves

Skin and body protection : No special precautions required.

Respiratory protection : Engineering controls should always be the primary method of controlling exposures.  
If respiratory protective equipment is needed for certain activities, the type as well as the corresponding protection factor will depend upon the risk assessment and air concentrations, hazards, physical and warning properties of substances present.  
No personal respiratory protective equipment normally required.

Protective measures : The type of protective equipment must be selected based on the Environmental Health and Safety risk assessment. Consult a Environmental Health and Safety expert if necessary.

### SECTION 9: Physical and chemical properties

#### 9.1 Information on basic physical and chemical properties

Appearance : liquid

Colour : clear

#### 9.2 Other information

No data available

### SECTION 10: Stability and reactivity

#### 10.1 Reactivity

None reasonably foreseeable.

#### 10.2 Chemical stability

Stable under recommended storage conditions.

#### 10.3 Possibility of hazardous reactions

Hazardous reactions : No dangerous reaction known under conditions of normal use.

#### 10.4 Conditions to avoid

Conditions to avoid : To avoid thermal decomposition, do not overheat.

#### 10.5 Incompatible materials

Materials to avoid : None known.

## CellSave Preservative 100 Tubes

Version	Revision Date:	SDS Number:	Date of last issue: 2017-08-02
1.44	2018-01-09	100000011151	Date of first issue: 2015-02-26

### 10.6 Hazardous decomposition products

None known.

## SECTION 11: Toxicological information

### 11.1 Information on toxicological effects

#### Acute toxicity

##### Components:

##### IMIDUREA:

Acute oral toxicity : LD50 (Rat): 11 300 mg/kg

Acute inhalation toxicity : LC50 (Rat): > 5,5 mg/l  
Exposure time: 1 h

Acute dermal toxicity : LD50 (Rabbit): > 5 000 mg/kg

#### Skin corrosion/irritation

##### Components:

##### IMIDUREA:

Result: No skin irritation

#### Serious eye damage/eye irritation

##### Components:

##### IMIDUREA:

Result: No eye irritation

#### Respiratory or skin sensitisation

##### Components:

##### IMIDUREA:

Method: Maximisation Test  
Result: May cause sensitisation by skin contact.

Method: Local Lymph Node Assay (LLNA) in mice  
Result: May cause sensitisation by skin contact.

#### Germ cell mutagenicity

##### Components:

##### IMIDUREA:

Genotoxicity in vitro : Test Type: Ames test  
Metabolic activation: with and without metabolic activation  
Result: negative  
GLP: yes

: Test Type: Chromosome aberration test in vitro  
Metabolic activation: with and without metabolic activation

## CellSave Preservative 100 Tubes

Version 1.44	Revision Date: 2018-01-09	SDS Number: 100000011151	Date of last issue: 2017-08-02 Date of first issue: 2015-02-26
-----------------	------------------------------	-----------------------------	---

Result: negative

Genotoxicity in vivo : Test Type: Micronucleus test  
Species: Mouse  
Application Route: Oral  
Result: negative

Germ cell mutagenicity- Assessment : No information available.

### Carcinogenicity

#### Components:

##### IMIDUREA:

Carcinogenicity - Assessment : No information available.

### Reproductive toxicity

#### Components:

##### IMIDUREA:

Teratogenicity - Assessment : No information available.

### STOT - single exposure

No data available

### STOT - repeated exposure

No data available

### Repeated dose toxicity

#### Components:

##### IMIDUREA:

Species: Rat  
NOAEL: 200 mg/kg  
LOAEL: 500 mg/kg  
Application Route: Oral

Species: Rabbit  
NOAEL: 200 mg/kg  
Application Route: Dermal

### Aspiration toxicity

No data available

---

## SECTION 12: Ecological information

### 12.1 Toxicity

#### Components:

##### IMIDUREA:

## CellSave Preservative 100 Tubes

Version	Revision Date:	SDS Number:	Date of last issue: 2017-08-02
1.44	2018-01-09	100000011151	Date of first issue: 2015-02-26

Toxicity to fish : Remarks: No data available

### 12.2 Persistence and degradability

**Components:**

**IMIDUREA:**

Biodegradability : Remarks: No data available

### 12.3 Bioaccumulative potential

**Components:**

**IMIDUREA:**

Bioaccumulation : Remarks: No data available

### 12.4 Mobility in soil

**Components:**

**IMIDUREA:**

Distribution among environmental compartments : Remarks: No data available

### 12.5 Results of PBT and vPvB assessment

**Product:**

Assessment : This substance/mixture contains no components considered to be either persistent, bioaccumulative and toxic (PBT), or very persistent and very bioaccumulative (vPvB) at levels of 0.1% or higher..

**Components:**

**IMIDUREA:**

Assessment : No information available..

### 12.6 Other adverse effects

**Components:**

**IMIDUREA:**

Additional ecological information : Remarks: No data available

---

## SECTION 13: Disposal considerations

### 13.1 Waste treatment methods

Product : In accordance with National, Federal, State and Local regulations.

Contaminated packaging : Empty containers should be taken to an approved waste han-



## CellSave Preservative 100 Tubes

Version	Revision Date:	SDS Number:	Date of last issue: 2017-08-02
1.44	2018-01-09	100000011151	Date of first issue: 2015-02-26

ding site for recycling or disposal.

### SECTION 14: Transport information

#### 14.1 UN number

Not regulated as a dangerous good

#### 14.2 UN proper shipping name

Not regulated as a dangerous good

#### 14.3 Transport hazard class(es)

Not regulated as a dangerous good

#### 14.4 Packing group

Not regulated as a dangerous good

#### 14.5 Environmental hazards

Not regulated as a dangerous good

#### 14.6 Special precautions for user

Not applicable

#### 14.7 Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable for product as supplied.

### SECTION 15: Regulatory information

#### 15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

Seveso III: Directive 2012/18/EU of the European Parliament and of the Council on the control of major-accident hazards involving dangerous substances.

Not applicable

Other regulations : Restricted to professional users.

#### 15.2 Chemical safety assessment

A Chemical Safety Assessment is not applicable (mixture)

### SECTION 16: Other information

#### Full text of H-Statements

H317 : May cause an allergic skin reaction.

#### Date and Number Formats

This document uses the following notation for printing dates and numbers:

Date: Dec 31th, 2012 as 2012-12-31

## CellSave Preservative 100 Tubes

Version	Revision Date:	SDS Number:	Date of last issue: 2017-08-02
1.44	2018-01-09	100000011151	Date of first issue: 2015-02-26

**Numbers:** 123456,78 as 123 456,78

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

GB / EN